


K.T. SURGICO

The Finest Quality Surgical Instruments



K.T. SURGICO

Manufacturers & Exporters of Quality Surgical, Dental & Manicure Instruments

Certified Under: cGMP - QSR BY FDA (U.S.A), ISO - 9001 : 2000 & **CE** Mark

WWW.KTSURGICO.COM

	Introduction	
	A Hemostat Forceps	05
	B Scissors & Scalpel	24
	C Needle Holder & Retractor	38
	D Diagnostic Set, General Instruments & Orthopedic	55
	E Neurosurgery & Ent	76
	F Cardiovascular, Suction Tubes & Intestinal / Rectal	90
	G Urology Gynecology & Obstetric	102
	H Titanium	112
	Suggested Surgical Instruments Sets	121
	Index	127



**K.T. SURGICO***Manufacturers & Exporters of Quality Surgical Instruments**K.T. Road, Rang Pura Sialkot, Pakistan***THE MARKET PLACEMENT OF SURGICAL INSTRUMENTS GLOBALLY.****WELCOME TO K.T. SURGICO**

For the last 30 years, we have been manufacturing Surgical Instruments. We work extremely hard to satisfy our customer needs by producing the finest quality instruments.

At **K. T. SURGICO**, we are pleased to inform you that we make sure the every piece we ship out, exceeds customers specification and expectations.

We are ISO9001-2000 approved facility and also certified by FDA authority USA as well as **CE** mark.

To cater the market needs, we produce our instruments in O.R. grade and Econo-line all O.R. grade instruments cover with life time warranty and Econo-line instruments cover with one full year warranty.

We have hand selected about 1100+ commonly used items for this catalogue. However, if you do not find, what you need in this catalogue, we would be more than happy to assist you in producing counter samples from your master sample, from a sketch / photograph or even a universal part.

We thank you for taking time to go through this catalogue and look forward to assist you any of your surgical instruments needs.

Mian Khalid Mehmood
C.E.O



A

- **Hemostats**
- **Forceps**

**Category A:**


In this category we have included the hemostat forceps & thumb forceps. You will find Kelly, Crile, Mosquito, Pean, Rochester and other hemostats in different sizes both straight & curved for dissecting, grasping, holding & handling the surgeries. In the same section you will see adson, adson-brown, standard & micro thumb forceps T/C & without T/C for tissue holding & dressing purposes.





Bulldog Clamps / Vascular



Atrauma Debaque



Diethrich/Micro
Coronary



KTS1207001
Titan (4 cm)



KTS1207101
Titan (4 cm)





KTS1207301
Titan (4 cm)



KTS1211000
John-Hopkins
(24.5 cm) 9 3/4"



Pennington / Siezing



Halsted Mosquito



1/2

Halsted-Mosquito




KTS1222112
5" (12.5cm)
Curved, serrated
KTS1222114
5 1/2" (14cm)
Curved, serrated



KTS1223112
5" (12.5cm)
Curved, 1:2 Teeth
KTS1223114
5 1/2" (14cm)
Curved, 1:2 Teeth



KTS1222012
5" (12.5cm)
straight, serrated
KTS1222014
5 1/2" (14cm)
straight, serrated




KTS1223012
5" (12.5cm)
Straight, 1:2 Teeth
KTS1223014
5 1/2" (14cm)
Straight, 1:2 Teeth



1/2

Hartmann / Mosquito



KTS1221010
4" (10cm)
Straight, serrated



KTS1221109
3 1/2" (9.0cm)
Curved, Serrated
KTS1221110
4" (10cm)
Curved, Serrated



1/2

Micro-Mosquito /
Petit Point





KTS1222210
4" (10 cm)
Straight, serrated
KTS1222212
4 3/4" (12 cm)
Straight, serrated




KTS1222310
4" (10 cm)
Curved, serrated
KTS1222312
4 3/4" (12 cm)
Curved, serrated



Vascular



Vascular




KTS1226216
6 1/4" (16cm)
Straight, serrated

Crile-Rankin



KTS1226316
6 1/4" (16cm)
Curved, serrated



KTS1228014
5 1/2" (14cm)
Straight, serrated
KTS1228016
6 1/4" (16cm)
Straight, serrated

Coller



KTS1228114
5 1/2" (14cm)
Curved, serrated
KTS1228116
6 1/4" (16cm)
Curved, serrated



KTS1225014
5 1/2" (14cm)
Straight, serrated

Crile



KTS1225114
5 1/2" (14cm)
Curved, serrated



KTS1226014
5 1/2" (14cm)
Straight, 1:2 Teeth



KTS1226114
5 1/2" (14cm)
Curved, 1:2 Teeth




KTS1231214
5 1/2" (14cm)
Straight, Serrated
KTS1231216
6 1/4" (16cm)
Straight, Serrated

Pean - Delicate



KTS1231314
5 1/2" (14cm)
Curved, Serrated
KTS1231316
6 1/4" (16cm)
Curved, Serrated


Vascular



KTS1229518
7 1/2" (19cm)
Curved, serrated

1/2

Jacobson



KTS1227014
5 1/2" (14cm)
Straight, serrated


1/2

Providence / Hospital



KTS1227114
5 1/2" (14cm)
Curved, serrated

A



KTS1231016
6 1/4" (16cm)
straight,
serrated

1/2

Rochester / Pean



KTS1231018
7 1/4" (18cm)
straight,
serrated



KTS1231020
8" (20cm)
straight,
serrated




KTS1231030
12" (30cm)
straight,
serrated




KTS1231116
6 1/4" (16cm)
Curved,
serrated



KTS1231118
7 1/4" (18cm)
Curved,
serrated



KTS1231120
8" (20cm)
Curved,
serrated




KTS1231122
8 1/2" (22cm)
Curved,
serrated



KTS1231124
9 1/2" (24cm)
Curved,
serrated




KTS1231126
10 1/4" (26cm)
Curved,
serrated



KTS1231130
12" (30cm)
Curved,
serrated



Grasping



KTS1232016
6 1/4" (16cm)
straight,
1:2 Teeth



KTS1232018
7" (18cm)
straight,
1:2 Teeth



KTS1232020
8" (20cm)
straight,
1:2 Teeth



KTS1232022
8 3/4" (22cm)
straight,
1:2 Teeth



KTS1232024
9 1/2" (24cm)
straight,
1:2 Teeth



KTS1232026
10 1/4" (26cm)
straight,
1:2 Teeth



KTS1232030
12" (30cm)
straight,
serrated



KTS1232116
6 1/4" (16cm)
Curved,
1:2 Teeth



KTS1232118
7" (18cm)
Curved,
1:2 Teeth



KTS1232120
8" (20cm)
Curved,
1:2 Teeth



Rochester / Ochsner / Kocher 1/2



1/2 Kocher / Delicate



3/4 Spencer-Wells

KTS1232014
5 1/2" (14cm)
straight,
1:2 Teeth

KTS1232213
5" (13cm)
Straight, 1:2 Teeth
KTS1232214
5 1/2" (14cm)
Straight, 1:2 Teeth
KTS1232216
6 1/4" (16cm)
Straight, 1:2 Teeth


KTS1232314
5 1/2" (14cm)
curved, 1:2 Teeth

KTS1233115
6" (15.5cm)
curved, serrated

Vascular / Dissecting & Ligature / Tonsil Hemostatic Forceps



Baby-Adson




Adson




Meeker



Schmidt / Sawtell




Rochester / Carmalt



Sarot




Dissecting and Ligature, Right Angles / Thoracic / Hysterectomy




3/4
Bridge



1/2
Mixer



KTS1247123
9" (23cm)
Fully Curved,
Serrated




1/2
Gemini




1/2
Baby-Mixer



1/2
Mini-Gemini



3/4
Mixer




1/2
Kantrowitz



1/2
Lahey



Hysterectomy / Tonsil, Duct / Towel Clamp / Pedicle




Towel Clamps / Tubing Forceps



Roeder
KTS1411713
5" (13cm)
Curved,
Serrated




Non-Perforating



Lorna



KTS1412110
4" (10cm)
KTS1412111
4 1/4" (11cm)
KTS1412113
5" (13cm)




Ball & Socket



KTS1392310
4" (10cm)
KTS1392313
5" (13cm)
Non-perforating




Peers-Bertman



KTS1416514
5 1/4" (13.25cm)
Serrated



Tubing Clamp / Laryngeal



Dressing / Tissue



Standard / Thumb



KTS1010213
5" (13cm)
Dressing, Serrated




KTS1010214
5 1/2" (14.5cm)
Dressing, Serrated



KTS1010216
6 1/4" (16cm)
Dressing, Serrated



KTS1010225
10" (25cm)
Dressing, Serrated



KTS1012014
5 1/2" (14.5cm)
Tissue, 1:2 Teeth



KTS1012016
6 1/4" (16cm)
Tissue, 1:2 Teeth



KTS1012025
10" (25cm)
Tissue, 1:2 Teeth



Model-USA



KTS1010414
5 1/2" (14.5cm)
Dressing, Serrated
KTS1010415
6" (15.25cm)
Dressing, Serrated
KTS1010418
7" (18cm)
Dressing, Serrated
KTS1010425
10" (25cm)
Dressing, Serrated



KTS1012214
5 1/2" (14.5cm)
Tissue, 1:2 Teeth
KTS1012215
6" (15.25cm)
Tissue, 1:2 Teeth
KTS1012225
10" (25cm)
Tissue, 1:2 Teeth



Narrow Tip



KTS1011014
5 1/2" (14.5cm)
Straight, Serrated
KTS1011016
6 1/4" (16cm)
Straight, Serrated
KTS1011020
8" (20cm)
Straight, Serrated



KTS1013014
5 1/2" (14.5cm)
Straight, 1:2 Teeth
KTS1013016
6 1/4" (16cm)
Straight, 1:2 Teeth
KTS1013020
8" (20cm)
Straight, 1:2 Teeth

Adson / Dressing / Tissue / Bayonet



Adson / Micro-Adson



KTS1017012
4 3/4" (12cm)
dressing, serrated



KTS1017612
4 3/4" (12cm)
dressing,
Micro serrated



KTS1017412
4 3/4" (12cm)
dressing, serrated
tungsten carbide



KTS1018012
4 3/4" (12cm)
Tissue, 1:2 Teeth




KTS1018612
4 3/4" (12cm)
Tissue, Micro 1:2 Teeth
KTS1018212
4 3/4" (12cm)
Tissue, 2:3 Teeth



Lucae / ENT



KTS4420614
5 1/2" (14cm)
Dressing, Serrated
KTS4420616
6 1/4" (16cm)
Dressing, Serrated



Lightweight-Adson



KTS1016812
4 3/4" (12cm)
Dressing,
Micro serrated



KTS1016912
4 3/4" (12cm)
Tissue,
Micro-1:2 Teeth



KTS1017312
4 3/4" (12cm)
Tissue, 1:2 Teeth



Adson



KTS1000112
4 3/4" (12cm)
Tungsten carbide,
cross serration



KTS1000212
4 3/4" (12cm)
Tungsten carbide,
cross serration-
1:2 Teeth



KTS1000215
6" (15.25cm)
Tungsten carbide,
cross serration-
1:2 Teeth




KTS1014212
4 3/4" (12cm)
Dressing,
cross serration




KTS1014412
4 3/4" (12cm)
Tissue, cross
serration-1:2 Teeth



Adson-Brown



KTS1018812
4 3/4" (12cm)
Tissue, 7:7 Teeth



KTS1019912
4 3/4" (12cm)
Tissue, 9:9 Teeth






1/2
Cushing

KTS1022017
7" (18cm)
Dressing w/Pin,
Serrated



KTS1022502
7" (18cm)
Tungsten carbide
Cross Serrated



KTS1023017
7" (18cm)
Tissue w/Pin,
1:2 Teeth



1/2
Oehler

KTS1026114
5.5" (14cm)
Dissecting-Tungsten
Carbide, Cross
Serrated-1:2 Teeth



KTS1026116
6.25" (16cm)
Dissecting-Tungsten
Carbide, Cross
Serrated-1:2 Teeth



KTS1026118
7" (18cm)
Dissecting-Tungsten
Carbide, Cross
Serrated-1:2 Teeth




KTS1026120
8" (20cm)
Dissecting-Tungsten
Carbide, Cross
Serrated-1:2 Teeth



KTS1026125
10" (25cm)
Dissecting-Tungsten
Carbide, Cross
Serrated-1:2 Teeth



KTS1026130
12" (30cm)
Dissecting-Tungsten
Carbide, Cross
Serrated-1:2 Teeth



1/2
Bonney Dissecting

KTS1132218
7" (18cm)
Straight Tungsten
Carbide, X-Serrations
1:2 Teeth



KTS1045118
7" (18cm)
Straight, Serrated-
1:2 Teeth




KTS1045218
7" (18cm)
Straight, Serrated-
2:3 Teeth



1/2
Brown

KTS1018915
6" (15cm)
Straight, 9:9 Teeth



3/4
Ewald / Hudson

KTS1015212
4 3/4" (12cm)
Dressing, Serrated




KTS1015412
4 3/4" (12cm)
Tissue, 1:2 Teeth








1/2
Ferris-Smith

KTS1024618
7" (18cm)
Straight, 1:2 Teeth
KTS1024818
7" (18cm)
Straight, 2:3 Teeth



Tissue / Dressing / Grasping / Splinter / Intestinal





Splinter

KTS1080209
3 1/2" (9cm)
Straight-No Pin,
Serrated
KTS1080211
4 1/2" (11.5cm)
Straight-No Pin,
Serrated



Carmalt


KTS1081610
4 1/4" (10.75cm)
Straight, Serrated
w/Pin

KTS1081710
4 1/4" (10.75cm)
Angled, Serrated
w/Pin



Jansen-
Gruenwald

KTS4615416
6 1/4" (16cm)
Bayonet, Serrated
KTS4615418
7" (18cm)
Bayonet, Serrated
KTS4615720
8" (20cm)
Bayonet, Serrated



De Bakey /
Lightweight

KTS1084015
6 1/4" (16cm)
Tissue,
width 1.5mm



DeBakey

KTS1130613
5" (12.75cm)
1.5mm
KTS1130616
6 1/4" (16cm)
1.5mm
KTS1130620
8" (20cm)
1.5mm
KTS1130624
9 1/4" (24cm)
1.5mm

KTS1130816
6 1/4" (16cm)
2.0mm
KTS1130820
8" (20cm)
2.0mm
KTS1130824
9 1/4" (24cm)
2.0mm
KTS1130830
12" (30.5cm)
2.0mm

KTS1130416
6 1/4" (16cm)
2.7mm
KTS1130420
8" (20cm)
2.7mm
KTS1133220
8" (20cm) Heavy
2.7mm
KTS1130424
9 1/4" (24cm)
2.7mm
KTS1130430
12" (30.5cm)
2.7mm

KTS1131020
8" (20cm)
3.5mm
KTS1131024
9 1/4" (24cm)
3.5mm
KTS1131030
12" (30.5cm)
3.5mm



Castroviejo / 1:2

KTS4315001
4" (10cm)
Straight, 1:2 Teeth
.12mm

KTS4315003
4" (10cm)
Straight, 1:2 Teeth
.3mm

KTS4315005
4" (10cm)
Straight, 1:2 Teeth
.5mm



Opthamology / Gyn / Suture



10 x
1/4
KTS4294407
2 3/4" (7cm)
Dressing-with
Pin, 1:2 Teeth

Modified Bonn



10 x
1/4
KTS4294607
2 3/4" (7cm)
Dressing-without
Pin, 1:2 Teeth

Modified Bonn



Bishop-Harman

KTS4294205
3 1/2" (8.5cm)
Lightweight,
1:2 Teeth

Fine
KTS4294208
3 1/2" (8.5cm)
1:2 Teeth

1/4
KTS4294005
3 1/2" (8.5cm)
Lightweight,
serrated,
fine 0.5mm

1/4
KTS4294008
3 1/2" (8.5cm)
Serrated
0.8mm

A



Jeweler's Forceps
KTS1083090
11.5cm 4 1/2"



KTS1083080
10.5cm 4 1/8"



KTS1083070
10.5cm 4 1/8"



KTS1083060
11cm 4 3/8"



KTS1083050
11cm 4 3/8"



KTS1083040
11cm 4 3/8"



KTS1083015
12cm 4 3/4"




KTS1083012
12cm 4 3/4"



KTS1083010
12cm 4 3/4"




Articulating Paper Forceps / Ophthalmology / Gyn / Suture / Dressing Forceps



KTS1060301



Miller
15cm 6"



KTS1060302



Miller
15cm 6"



Heywood-Smith



KTS6840021
8 1/4" (21cm)
Spread,
serrated



A



Bozemann



Foerster-Ballenger



KTS1617018
7" (17.75cm)
Straight, serrated



KTS1617025
10" (25.5cm)
Straight, serrated



KTS1617020
8" (20.25cm)
Straight, serrated



KTS1617225
10" (25.5cm)
Straight, Smooth



KTS1617120
8" (20.25cm)
Curved, serrated



KTS1617125
10" (25.5cm)
Curved, serrated



KTS1617118
7" (17.75cm)
Curved, serrated



KTS1617325
10" (25.5cm)
Curved, Smooth



B

- **Scissors**
- **Scalpels**

**Category B:**

We manufacture our 200 different scissors with T/C & without T/C from Iris to Mayo, Fergusson, Metzenbaum, Lister and many more. We have tried to include as many as possible different sizes & shapes of scissors to meet all your surgical needs. A complete range of most commonly used scalpel handles is in the end of this section.




Standard / Operating / Gyn / Nasal




1/2


Standard / Operating




KTS0810214
5 1/2"(14cm)
Straight, Blunt/Blunt
KTS0810215
6"(15.25cm)
Straight, Blunt/Blunt
KTS0810216
6 1/2"(16.5cm)
Straight, Blunt/Blunt



KTS0811014
5 1/2"(14cm)
Straight, Sharp/Blunt
KTS0811414
5 1/2"(14 cm)
Straight, Tungsten
carbide, Sharp/Sharp



KTS0812014
5 1/2"(14cm)
Straight,
Sharp/Sharp



1/2

Mayo

KTS0716014
5 1/2"(14cm)
Straight-Supercut
KTS0716017
6 3/4"(17cm)
Straight-Supercut
KTS0716117
6 3/4"(17cm)
Curved-Supercut



1/2

Jorgenson




KTS0865022
8 3/4"(22cm)
Curved



1/2

Mayo-Noble / Beveled Blades



KTS0818017
6 3/4"(17cm)
Straight



KTS0818117
6 3/4"(17cm)
Curved



1/2

Fergusson




KTS0801914
5 1/2"(14cm)
Angled
KTS0801918
7"(17.75cm)
Angled



B


Dissecting / Tenotomy



KTS0816014
5 1/2"(14cm)
Straight
KTS0816015
6"(15.25cm)
Straight
KTS0816017
6 3/4"(17cm)
Straight
KTS0816019
7 3/4"(19.75cm)
Straight
KTS0816023
9"(23cm)
Straight
KTS0816028
11"(28cm)
Straight

KTS0816114
5 1/2"(14cm)
Curved
KTS0816117
6 3/4"(17cm)
Curved
KTS0816119
7 3/4"(19.75cm)
Curved
KTS0816128
11"(28cm)
Curved

Mayo / Develed Blades



KTS0816414
5 1/2"(14cm)
Straight-Tungsten Carbide
KTS0816417
6 3/4"(17cm)
Straight-Tungsten Carbide
KTS0816423
9"(23cm)
Straight-Tungsten Carbide

KTS0816514
5 1/2"(14cm)
Curved-Tungsten Carbide
KTS0816517
6 3/4"(17cm)
Curved-Tungsten Carbide
KTS0816523
9"(23cm)
Curved-Tungsten Carbide

Mayo Dissecting



KTS0817017
6 3/4"(17cm)
Straight
KTS0817415
6"(15cm)
Straight-Tungsten Carbide
KTS0817417
6 3/4"(17cm)
Straight-Tungsten Carbide

KTS0817117
6 3/4"(17cm)
Curved
KTS0817517
6 3/4"(17cm)
Curved-Tungsten Carbide

Mayo-Stille /
Rounded Blades



KTS0828111
4 1/2"(11cm)
Curved
Blunt/Blunt

KTS0828011
4 1/2"(11cm)
Straight
Blunt/Blunt

KTS0728011
4 1/2"(11cm)
Straight-Supercut
Blunt/Blunt

Metzenbaum-Baby



Dissecting / Tenotomy / Vascular






Fig.	Item	Description	Tip
A	KTS0828015	5 1/2" (14cm) Straight	Blunt/Blunt
A	KTS0828018	7" (18cm) Straight	Blunt/Blunt
B	KTS0828914	5 1/2" (14cm) Curved	Sharp/Blunt
B	KTS0828115	5 1/2" (14cm) Curved	Blunt/Blunt
B	KTS0828118	7" (18cm) Curved	Blunt/Blunt
B	KTS0802718	7" (18cm) Curved	Sharp/Sharp
B	KTS0828120	8" (20cm) Curved	Blunt/Blunt
B	KTS0828123	9" (23cm) Curved	Blunt/Blunt
B	KTS0802723	9" (23cm) Curved	Sharp/Sharp
A	KTS0828415	5 1/2" (14cm) Straight Tungsten Carbide	Blunt/Blunt
B	KTS0828515	5 1/2" (14cm) Curved Tungsten Carbide	Blunt/Blunt
B	KTS0828518	7" (18cm) Curved Tungsten Carbide	Blunt/Blunt
B	KTS0828523	9" (23cm) Curved Tungsten Carbide	Blunt/Blunt
B	KTS0828530	12" (30cm) Curved Tungsten Carbide	Blunt/Blunt
B	KTS0728115	5 1/2" (14cm) Curved Supercut	Blunt/Blunt
B	KTS0728118	7" (18cm) Curved Supercut	Blunt/Blunt
B	KTS0728123	9" (23cm) Curved Supercut	Blunt/Blunt






Fig.	Item	Description	Tip
A	KTS0828315	5 1/2" (14cm) Curved	Blunt/Blunt
A	KTS0828318	7" (18cm) Curved	Blunt/Blunt
A	KTS0828328	11" (28cm) Curved	Blunt/Blunt
B	KTS0828615	5 1/2" (14cm) Straight Tungsten Carbide	Blunt/Blunt
A	KTS0828715	5 1/2" (14cm) Curved Tungsten Carbide	Blunt/Blunt
A	KTS0828718	7" (18cm) Curved Tungsten Carbide	Blunt/Blunt
A	KTS0828723	9" (23cm) Curved Tungsten Carbide	Blunt/Blunt
A	KTS0829311	4 1/2" (11.5cm) Curved Tungsten Carbide	Blunt/Blunt



Vascular/Nasal, O.B./Tenotomy




Hegemann-Diethrich




KTS0803025
7.5" (19cm)
25°



KTS0803045
7.5" (19cm)
45°




KTS0803060
7.5" (19cm)
60°



KTS0803090
7.5" (19cm)
90°



KTS0803012
7.5" (19cm)
125°



Lincoln



KTS0800217
6 3/4" (17cm)
Straight
Blunt / Blunt



KTS0800117
6 3/4" (17cm)
Curve
Blunt / Blunt

B




KTS0840918
7" (18cm)
Angled

Knight



Modified-USA



KTS0855010
4" (10cm)



Operating
Fine

KTS0800411
4 3/4" (12cm)
Straight,
sharp/sharp
KTS0800511
4 3/4" (12cm)
Curved,
sharp/sharp



O.B. / Eye / Tenotomy / Vascular








Fig.	Item	Description	Tip
A	KTS0834009	3 1/2" (09cm) Straight	Sharp/Sharp
A	KTS0834010	4 1/4" (10.5cm) Straight	Sharp/Sharp
B	KTS0734110	4 1/4" (10.5cm) Curved - Supercut	Sharp/Sharp
A	KTS0834011	4 1/2" (11cm) Straight	Sharp/Sharp
B	KTS0834109	3 1/2" (9cm) Curved	Sharp/Sharp
B	KTS0834110	4 1/4" (10.5cm) Curved	Sharp/Sharp
B	KTS0834111	4 1/2" (11.5cm) Curved	Sharp/Sharp
B	KTS0734111	4 1/2" (11.5cm) Curved - Supercut	Sharp/Sharp
A	KTS0834211	4 1/2" (11.5cm) Straight	Sharp/Sharp
B	KTS0834312	4 1/2" (11.5cm) Curved	Sharp/Sharp
A	KTS0834411	4 1/2" (11.5cm) Straight Tungsten Carbide	Sharp/Sharp
B	KTS0834511	4 1/2" (11.5cm) Curved Tungsten Carbide	Sharp/Sharp
C	KTS0837602	4 1/2" (11.5cm) Straight	Sharp/Blunt
D	KTS0837701	4 1/2" (11.5cm) Curved	Blunt/Blunt
D	KTS0837702	4 1/2" (11.5cm) Curved	Sharp/Blunt

Tenotomy / Oral / Stitch / Wire



Tenotomy-Stevens





Fig.	Item	Description	Tip
A	KTS0836111	4 1/2" (11.5cm) Curved	Sharp/Sharp
B	KTS0836210	4 1/8" (10.5cm) Straight	Blunt/Blunt
B	KTS0836211	4 1/2" (11.5cm) Straight	Blunt/Blunt
C	KTS0836310	4 1/8" (10.5cm) Curved	Blunt/Blunt
C	KTS0836311	4 1/2" (11.5cm) Curved	Blunt/Blunt
C	KTS0736311	4 1/2" (11.5cm) Curved - Supercut	Blunt/Blunt
C	KTS0836314	5 3/4" (14.5cm) Curved	Blunt/Blunt



Tenotomy




KTS0806109
3 3/4" (9.5cm)
Curved-Ribbon
Handles blunt/blunt


B




ENT-Dean




KTS0842717
6 3/4" (17cm)
Serrated



Spencer



KTS0880209
3 1/2" (9cm)




Spencer / Delicate

KTS0880111
4 1/4" (11cm)



Harvey



KTS0886212
5" (12.5cm)
Saw Edge




Bandage / Wire / Micro Scissors



$\frac{1}{2}$
Gorney Septal Scissors
KTS0895245
19.7cm 7 3/4"
double action



$\frac{1}{2}$
Utility Shears
KTS0885421
21.5cm 8 1/2"
Serrated blade



$\frac{1}{4}$
Micro Scissor
KTS0920418
18cm 7"

KTS0870622
straight
22cm 9"

KTS0870722
curved
22cm 9"



$\frac{1}{2}$
Caplan

B





Wolff
KTS3039518
18cm 7"
KTS3039724
24cm 9 1/2"



Stille
KTS3010320
20cm 8"
KTS3010323
23cm 9"
KTS3010326
26cm 10 1/4"
KTS3010337
37cm 14 1/2"



Hennig
KTS3035028
28cm 11"



Wire / Bandage / Ophthalmic



Wire-Universal




KTS0883112
4 3/4" (12cm)
Saw Edge
KTS0883512
4 3/4" (12cm)
Tungsten Carbide
Saw Edge



1/2
KTS0895014 5 3/4" (14.5cm)
143 C/290 F Autoclavable
KTS0895015 6" (15cm)
143 C/290 F Autoclavable
KTS0895016 6 1/4" (16cm)
143 C/290 F Autoclavable
KTS0895018 7" (18cm)
143 C/290 F Autoclavable



1/3
Bergmann
KTS0891123 9" (23cm)



Lister
Banding Scissors





Fig.	Item	Description
A	KTS0890109	3 1/2" (9cm)
A	KTS0890111	4 1/2" (11cm)
A	KTS0890114	5 1/2" (14cm)
A	KTS0890115	6" (15.5cm)
A	KTS0890118	7" (18cm)
A	KTS0890518	7" (18cm) Tungsten Carbide
A	KTS0890120	8" (20cm)
A	KTS0890121	8 1/4" (21cm)



Westcott

KTS0910211 4 1/2" (11cm)
Curved,
Sharp/Sharp
KTS0910311 4 1/2" (11cm)
Curved-Micro,
Sharp/Sharp
KTS0910411 4 1/2" (11cm)
Curved,
Blunt/Blunt



Esmarch

KTS0891320
8" (20cm)
KTS0891322
8 3/4" (22cm)



B

Handles / Amputation



Corwin Round

KTS0610116
6" (15.25cm)
10, 11, 12, 15



#3 Knife Handle

KTS0610300
5" (12.75cm)
10, 11, 12, 15



#3 Knife Handle/Graduated

KTS0610301
5" (12.75cm)
10, 11, 12, 15



#3 Knife Handle/Long

KTS0610302
9 1/4" (12cm)
10, 11, 12, 15



#3 Knife Handle / Long Angled

KTS0610303
8" (20.25cm)
10, 11, 12, 15

B



#7 Knife Handle

KTS0610700
6 1/2" (16.5cm)
10, 11, 12, 15



#4 Knife Handle

KTS0610400
5" (12.75cm)
20, 22, 25



Collin

KTS0642028
11 3/4" (28cm)



Collin

KTS0610901
No. 9



Fig. 10



Fig. 11



Fig. 12



Fig. 15



Fig. 20



Fig. 22



Knives



Dieffenbach



Fig. 1
KTS0639801



Fig. 2
KTS0639802



Fig. 3
KTS0639803



Fig. 4
KTS0639804



Fig. 5
KTS0639805



Fig. 6
KTS0639806



Fig. 7
KTS0639807



Fig. 8
KTS0639808



Fig. 9
KTS0639809



Fig. 15
KTS0639815



Fig. 16
KTS0639816

B



Knives




Fig. 1
KTS0640001

Fig. 2
KTS0640002

Fig. 3
KTS0640003




Fig. 4
KTS0640004

Fig. 5
KTS0640005

Fig. 6
KTS0640006



Fig. 7
KTS0640007

Fig. 8
KTS0640008

Fig. 9
KTS0640009

Fig. 10
KTS0640010



Fig. 1
KTS0640601

Fig. 2
KTS0640602

Fig. 3
KTS0640603




Fig. 4
KTS0640604

Fig. 5
KTS0640605

Fig. 6
KTS0640606




Collin



Dieffenbach



Knives



17cm 6 3/4"

1/4

Fig. 1
KTS0640401



Fig. 3
KTS0640403



Fig. 5
KTS0640405



Fig. 7
KTS0640407



Fig. 9
KTS0640409

B



1/4

Fig. 10
KTS0640410



Fig. 12
KTS0640412



Fig. 14
KTS0640414



Fig. 16
KTS0640416



Fig. 18
KTS0640418



C

- **Needle Holders**
- **Retractors**




Category C:

A wide range of needle holders in different shapes; straight, curved, plain and with or without teeth, T.C. for all kinds of suturing needs are in this section. In addition we have also included some micro needle holders for delicate needs. In the second portion of this section we have hand selected , self-retaining, rake, vein & ortho retractors.



Gyne



Baumgartner




KTS2415014
5 1/2" (14cm)
Serrated-Solid




KTS2415412
5" (13cm)
Tungsten Carbide
Serrated-Solid



KTS2415414
5 1/2" (14cm)
Tungsten Carbide
Serrated-Solid




Webster




KTS2410013
5" (13cm)
Smooth



KTS2410413
5" (13cm)
Tungsten
Carbide
Smooth



Halsey



KTS2413013
5" (13cm)
Smooth




KTS2413413
5" (13cm)
Tungsten
carbide
Smooth



KTS2414013
5" (13cm)
Serrated-Solid



KTS2414413
5" (13cm)
Tungsten
carbide
Serrated-Solid



Heaney-Gyne




KTS2424121
8 1/4" (21cm)
Serrated-Solid



KTS2424521
8 1/4" (21cm)
Tungsten Carbide
Serrated-Solid



KTS2424526
Tungsten Carbide
Serrated-Solid
10 1/4" (26cm)



Baby Crile-Wood



KTS2416115
6" (15.25cm)
Serrated-Solid



KTS2416215
6" (15cm)
Tungsten Carbide
Serrated-Solid



Brown



KTS2434413
5" (13cm)
Serrated-Solid



KTS2434418
7" (19cm)
Serrated-Solid



Mayo-Hegar



1/2

Crite Wood



KTS2416015
6" (15.25cm)
Serrated-Grooved



KTS2416018
7" (18cm)
Serrated-Grooved



KTS2416430
12" (30cm)
Serrated-Grooved



KTS2416415
6" (15cm)
Tungsten Carbide
Serrated-Solid



KTS2416418
7" (18cm)
Tungsten Carbide
Serrated-Solid



KTS2416420
8" (20cm)
Tungsten Carbide
Serrated-Solid



KTS2416423
9" (23cm)
Tungsten Carbide
Serrated-Solid



KTS2416615
6" (15.25cm)
Tungsten Carbide
Smooth

C



1/2

Mayo Hegar



KTS2418014
5 1/2" (14cm)
Serrated-Grooved



KTS2418016
6 1/4" (16cm)
Serrated-Grooved



KTS2418018
7" (18cm)
Serrated-Grooved



KTS2418020
8" (20cm)
Serrated-Grooved



KTS2418024
9 1/2" (24cm)
Serrated-Grooved



KTS2418026
10 1/4" (26cm)
Serrated-Grooved



KTS2418030
12" (30.5cm)
Serrated-Grooved



KTS2418414
5 1/2" (14cm)
Tungsten Carbide
Serrated-Solid



KTS2418416
6 1/4" (16cm)
Tungsten Carbide
Serrated-Solid



KTS2418418
7" (18cm)
Tungsten Carbide
Serrated-Solid



KTS2418420
8" (20cm)
Tungsten Carbide
Serrated-Solid



KTS2418424
9 1/2" (24cm)
Tungsten Carbide
Serrated-Solid



KTS2418426
10 1/4" (26cm)
Tungsten Carbide
Serrated-Solid



KTS2418430
12" (30.5cm)
Tungsten Carbide
Serrated-Solid




Vascular



 KTS2442014 5 1/2" (14cm) Serrated-Solid	 KTS2442017 6 3/4" (17cm) Serrated-Solid	 KTS2442019 7 1/2" (19cm) Serrated-Solid	 KTS2421418 7" (18cm) Tungsten Carbide Serrated-Solid	 KTS2421420 8" (20cm) Tungsten Carbide Serrated-Solid	 KTS2421423 9" (23cm) Tungsten Carbide Serrated-Solid
 KTS2441414 5 1/2" (14cm) Tungsten Carbide Serrated-Solid	 KTS2441417 6 3/4" (17cm) Tungsten Carbide Serrated-Solid	 KTS2441419 7 1/2" (19cm) Tungsten Carbide Serrated-Solid	 KTS2421426 10 1/4" (26cm) Tungsten Carbide Serrated-Solid	 KTS2421431 12 1/4" (31cm) Tungsten Carbide Serrated-Solid	



 1/2
 Olsen-Hegar/With Scissor


 1/2
 Debaquey-Vascular/Extra Light

 KTS2426027 10 1/2" (26.5cm) Serrated-Solid	 KTS2422614 5 1/2" (14cm) Tungsten Carbide Serrated-Delicate	 KTS2419615 6" (15cm) Tungsten Carbide Serrated-Solid	 KTS2419818 7" (18cm) Tungsten Carbide Serrated-Solid
 KTS2426427 10 1/2" (26.5cm) Tungsten Carbide Serrated-Solid			


 1/3
 Wangenstein


 1/2
 Eufate Pasque


 1/2
 Carroll


 1/2
 New Orleans











Micro



C



Ryder French-Eye/Standard/Micro

Micro / Needles



KTS2466114
5 1/2" (14cm)
Curved
Smooth

1/2


Barraquer



KTS2466314
5 1/2" (14cm)
Curved
Smooth

1/2

Barraquer/with Catch



1/2

Castroviejo



KTS2482214
5 1/2" (14cm)
Straight
Smooth



KTS2482414
5 1/2" (14cm)
Straight
Serrated-Solid



KTS2482514
5 1/2" (14cm)
Curved
Smooth



KTS2467414
5 1/2" (14cm)
Straight
Tungsten Carbide
Smooth




KTS2467614
5 1/2" (14cm)
Straight
Tungsten Carbide
Serrated-Solid



KTS2467714
5 1/2" (14cm)
Curved
Tungsten Carbide
Serrated-Solid



KTS2467618
7" (18cm)
Straight
Tungsten Carbide
Serrated-Solid



KTS2429818
7" (18cm)
Smooth



KTS2429821
8 1/4" (21cm)
Smooth



KTS2469421
8 1/4" (21cm)
Tungsten Carbide
Smooth

KTS2469521
8 1/4" (21cm)
Tungsten Carbide
Serrated-Solid

1/2

Jacobson



1/2

Reverdin Needle
19.5 cm/7 3/4"
22.5 cm/8 3/4"



1/4

KTS2620119
KTS2620123



1/4

KTS2620219
KTS2620223



1/4

KTS2620319
KTS2620323




1/4


KTS2620419
KTS2620423








Delicate Hooks/Hand Retractors



C



Rake/Vein, Hand Retractors



Rake/Hand Retractors, Ortho



Volkmann
KTS1828001
8 1/2" (21.5cm)
1-Prong-Sharp



KTS1828002
8 1/2" (21.5cm)
2-Prong-Sharp



KTS1828003
8 1/2" (21.5cm)
3-Prong-Sharp



KTS1828004
8 1/2" (21.5cm)
4-Prong-Sharp



KTS1828006
8 1/2" (21.5cm)
6-Prong-Sharp



KTS1828101
8 1/2" (21.5cm)
1-Prong-Blunt



KTS1828102
8 1/2" (21.5cm)
2-Prong-Blunt



KTS1828103
8 1/2" (21.5cm)
3-Prong-Blunt



KTS1828104
8 1/2" (21.5cm)
4-Prong-Blunt



KTS1828106
8 1/2" (21.5cm)
6-Prong-Blunt

C



KTS1833304
8 3/4" (22.25cm)
4-Prong-Sharp
11 x 18mm



KTS1833306
8 3/4" (22.25cm)
6-Prong-Sharp
11 x 27mm



KTS1836608
9.5" (24cm)
8-Prong-Sharp
26 x 40mm



KTS1836440
10" (25.5cm)
4-Prong-Dull
40 x 40mm



KTS3281106
6 1/4" (16cm)
Dull 6mm
KTS3281108
6 1/4" (16cm)
Sharp 8mm

1/3
Kocher

1/3
Koerte

1/3
Israel

1/2
Mini Hohmann



Hand Retractors



1/3

Hohmann

KTS3280110
9 1/2" (24cm)
Sharp 10mm
KTS3280118
9 1/2" (24cm)
Sharp 18mm
KTS3280724
9 1/2" (24cm)
Dull 24mm



1/3

Richardson-Eastman

KTS1814000
Set of 2
See Below
KTS1814001
10 1/4" (26cm)
Small
28 x 20mm & 36 x 28mm
KTS1814002
11" (28cm)
Large
49 x 38mm & 63 x 49mm



1/2

Lahey

KTS1830419
7 1/2" (19cm)
28 x 6mm




1/2

Green



1/4

KTS1834822
8 3/4" (22cm)
25 x 12mm



1/4

KTS1852301G
11 1/2" (29.25cm)
35 x 45mm

1/4

KTS1852302G
13" (32.5cm)
60 x 45mm

1/4

KTS1852303G
13" (32.5cm)
60 x 45mm



Hand Retractors



1/3

Richardson/Kelly



1/4

KTS1815231
9 1/2" (24cm)
Grip Handle
25 x 31mm



KTS1815251
9 1/2" (24cm)
Grip Handle
19 x 51mm



KTS1815851
9 1/2" (24cm)
Grip Handle
38 x 51mm



KTS1815876
9 1/2" (24cm)
Grip Handle
64 x 76mm

C



1/3

Richardson/Kelly



1/4

KTS1815001
9 1/2" (24cm)
Loop Handle
28 x 20mm




KTS1815002
9 1/2" (24cm)
Loop Handle
36 x 28mm



KTS1815003
9 1/2" (24cm)
Loop Handle
44 x 38mm



KTS1815004
9 1/2" (24cm)
Loop Handle
52 x 22mm



KTS1815626
9 1/2" (24cm)
Loop Handle
65 x 50mm



1/4

Harrington/Sweetheart

KTS1819140
12 1/2" (32cm)
127 x 40mm
KTS1819162
12 1/2" (32cm)
127 x 62mm
KTS1819164
12 1/2" (32cm)
174 x 64mm




Hand Retractors



KTS1837825
10" (25cm)
50 x 85mm

Doyen




KTS1817612
8 1/4" (21cm)
12mm

S-Shaped



KTS1852075
9" (22.5cm)
75 x 25mm

Meyerding



KTS1813121
8 1/4" (21cm)
Double Ended
Set of 2
21 x 24, 40 x 15
25 x 24, 45 x 15



Parker-Langenbeck/Army-Navy



Deaver



KTS1818325
12" (30cm)
25mm




KTS1818750
12" (30cm)
50mm



Hand Retractors / Self Retaining Retractors



C



Self Retaining Retractors



KTS1867820
8" (20cm)
Sharp 4:4
20mm x 20mm
KTS1867920
8" (20cm)
Blunt 4:4
19mm x 19mm

Anderson-Adson




KTS1869031
12" (31cm)
Sharp 4:4
19mm x 25mm
KTS1800631
12" (31cm)
Blunt 4:4
13mm x 25mm

Beckmann



KTS1869132
13" (33cm)
Blunt 7:7
47mm x 44mm

Beckmann/Eaton



Weitlaner



A



B



C

Fig.	Item	Description	Blade (Depth x width)
A	KTS1865010	4 1/4" (10.5cm) Sharp, 2:3	15mm x 12mm
A	KTS1865110	4 1/4" (10.5cm) Blunt, 2:3	15mm x 12mm
B	KTS1865413	5" (13cm) Sharp, 3:4	18.5mm x 18.5mm
B	KTS1865513	5" (13cm) Blunt, 3:4	18.5mm x 18.5mm
B	KTS1867013	5" (13cm) Sharp, 3:4	18.5mm x 18.5mm
B	KTS1867113	5" (13cm) Blunt, 3:4	18.5mm x 20mm
B	KTS1867016	6 1/2" (16.5cm) Sharp, 3:4	18.5mm x 18.5mm
B	KTS1867116	6 1/2" (16.5cm) Blunt, 3:4	18.5mm x 18.5mm
B	KTS1867020	8" (20cm) Sharp, 3:4	20mm x 20mm
B	KTS1867120	8" (20cm) Blunt, 3:4	20mm x 20mm
B	KTS1867025	10" (25cm) Sharp, 3:4	30mm x 20mm
B	KTS1867125	10" (25cm) Blunt, 3:4	23.5mm x 19mm
C	KTS1864826	10 1/2" (26cm) Sharp, 5:6	27mm x 32mm
C	KTS1864926	10 1/2" (26cm) Blunt, 5:6	27mm x 32mm




Self Retaining Retractors / Abdominal



KTS1800731
12" (31cm)
Sharp 4:4
25mm x 31mm
KTS1801031
12" (31cm)
Blunt 4:4
25mm x 31mm
KTS1801045
12" (31cm)
Blunt 4:4
25mm x 45mm

Beckmann/Adson



KTS1800513
5" (31cm)
Blunt 3:4
13mm x 14.5mm
KTS1800416
6 1/2" (16.5cm)
Sharp 3:4
17mm x 19mm
KTS1800516
6 1/2" (16.5cm)
Blunt 3:4
17mm x 19mm
KTS1800520
8" (20cm)
Blunt 3:4
21mm x 210mm

Cone



KTS1880209
3 1/2" (9cm)
Sharp 6mm
KTS1880211
4 1/2" (11cm)
Sharp 10mm
KTS1880214
5 1/2" (14cm)
Sharp 10mm
KTS1880218
7 1/2" (19cm)
Sharp 10mm


Gelpi

C



KTS1882018
7" (18cm)
Left
55 x 25mm
KTS1882118
7" (18cm)
Right
55 x 25mm

Markham/Meyerding



KTS1891800
10" (25cm)
45 x 80mm
70 x 35mm
100 x 35mm

Balfour/US Type



KTS1876690
3 1/2" (90mm)
24 x 20mm
27 x 27mm


Balfour/Baby



Self Retaining Retractors/Abdominal




KTS1875402
Single Retractor Blade
110 x 70mm



KTS1875403
Set of 2 Retractor Blade
80 x 45mm


KTS1875400
Set includes frame
& 3 blades
KTS1875401
Frame Only
(not pictured)

O'Sullivan/O'Connor



KTS1876418
7" (18cm)
Spread
45 x 80mm
63 x 35mm

Balfour



KTS1876545
Center Blade
45 x 80mm



Self Retaining Retractors / Abdominal, Malleable






Fig.	Item	Description	Tip
A	KTS1816812	8" (20cm)	12mm
A	KTS1816817	8" (20cm)	17mm
A	KTS1816820	8" (20cm)	25mm
A	KTS1816825	10" (25cm)	25mm
A	KTS1817025	13" (33cm)	25mm
A	KTS1817032	13" (33cm)	32mm
A	KTS1817040	13" (33cm)	40mm
A	KTS1817050	13" (33cm)	50mm



D

- **DIAGNOSTIC SET**
- **GENERAL INSTRUMENTS**
- **ORTHOPEDIC**

**D****Category D:**

Orthopedic surgeries require from very light to very heavy duty instruments. And yes, we have carefully selected a wide range of instruments to suit all orthopedic surgery needs. A complete line of diagnostic sets, stethoscope, percussion hammers & general instruments are part of this category.




Varicose Vein Prob Set, Cotton Applicators



Doyen
KTS2011114
5 1/2" (14.5cm)



Nelaton
KTS2011716
6 1/4" (16cm)



Lockhart-Mummery
KTS2012101
8" (20cm)



Nabatoff
KTS2015000
varicose vein probe
set complete



KTS2015002
Plastic Probe tip




KTS2015003
3mm Ø
Olive




KTS2015006
6mm Ø
Olive




KTS2015009
9mm Ø
Olive



KTS2015012
12mm Ø
Olive



KTS2015015
15mm Ø
Olive




KTS2015080
Handle, T-form



KTS2015090
Traction Cable



Probes / Cuticle Nail Nippers / Hygiene



Bowman-Lacrimal /
Double Ended-Malleable

Fig.	Item	Description	Tip
A	KTS4261599	4 7/8" (12.4cm)	0000/000
A	KTS4261600	4 7/8" (12.4cm)	00/0
A	KTS4261612	4 7/8" (12.4cm)	1/2
A	KTS4261634	4 7/8" (12.4cm)	3/4
A	KTS4261656	4 7/8" (12.4cm)	5/6
A	KTS4261678	4 7/8" (12.4cm)	7/8



*Nail Nippers / Ring Cutters / Umbilical Cord Clamp Cutter /
Dermal / Needle / Needle Diagnostics*



Nail Nipper



Nail Nipper



Ring Cutter
KTS3317500
6 1/2" (17cm)



KTS3317301
Saw Blade Only
for KTS3317300



Ring Cutter
KTS3317300
6 3/4" (17cm)



Waller
KTS6622300
5" (12.5cm)
Tungsten Carbide



Fox



Verres



Stethoscopes / Percussion Hammers



1/2
Buck Hammer
KTS0230818
7" (18cm)



1/2
Dejerine
KTS0230420
8" (20 cm)



1/2
Taylor Hammer
KTS0228817
7" (18cm)
KTS0228820
8" (20cm)



1/2
Wartenberg
KTS0235200
19cm



Ø 47mm

Duplex
KTS0231528
super lightweight



Ø 37mm

Duplex-Baby
KTS0231529
super lightweight



Ø 37mm

Anaestophon
KTS0231530
super lightweight



D



Diagnostic Sets (color available)
KTS0025710



Diagnostic Sets
KTS0025712

D



Diagnostic Sets
KTS0025714



Diagnostic Sets
KTS0025716



Diagnostic Sets
KTS0025718





Diagnostic Sets
KTS0025720



Diagnostic Sets
KTS0025722



Guedel-Negus / Miller

KTS0422500
with three laryngoscope
blades No. 1 to 3
complete, in case.

KTS0422600
With four laryngoscope
blades No. 1 to 4
complete, in case.



KTS0422501
for children
working length
70mm



KTS0422502
for adolescents
working length
90mm



KTS0422503
for women
working length
110mm



KTS0422504
for men
working length
175mm





McIntosh

KTS0420300
with three laryngoscope
blades No. 1 to 3
complete, in case.

KTS0420400
With four laryngoscope
blades No. 1 to 4
complete, in case.



KTS0421000
for babies
working length
55mm



KTS0421001
for children
working length
75mm



KTS0421002
for adolescents
working length
90mm



KTS0421004
for men
working length
135mm



KTS0421003
for women
working length
110mm



Foregger

KTS0422300
with three laryngoscope
blades No. 2 to 4
complete, in case.



KTS0422001
for children
working length
70mm



KTS0422002
for adolescents
working length
90mm



KTS0422003
for women
working length
110mm



KTS0422004
for men
working length
175mm



*Shears / Cast Spreader / Bone Holding Clamp /
Bone Reduction Forceps / Wire Cutter Double Action*



Model USA




KTS3036623
9" (23cm)
3 Prong




Plaster-Bruns

KTS0893723
9 1/2" (24cm)
Saw Edge



Lewin



KTS3240217
7" (17.5cm)



Ortho

KTS3240913
5 1/4" (13cm)



Ortho



KTS3241113
5 1/4" (13.5cm)
Angled Serrated



Wire Cutter

KTS3354215
6" (15cm)
W / Spring



Diamant Pin Cutter

KTS3354016
6 1/4" (16cm)
Single Action 2 mm



D


Wire Cutter Double Action / Wire Pin Cutting / Bolt Cutters / Ortho Pliers



Wire Cutter




KTS3357218
7" (18cm)
Double Action
Tungsten Carbide
Soft = 2.00mm
Hard = 1.5mm



Wire Cutter



KTS3354517
6 3/4" (17cm)
Double Action
Tungsten Carbide
Soft = 2.00mm



KTS3354522
8 3/4" (22cm)
Double Action
Tungsten Carbide
2.5mm



Wire Cutter



Lewin / Bolt /
Pin Cutter

KTS3304056
22" (56cm)
6.3mm



Wire Cutter


KTS3200218
5.5" (14cm)
1.2mm




Needle Nose Pliers




KTS3356614
5 1/2" (14cm)
Cross Serrated



Needle Nose Pliers



KTS3352219
7 1/2" (19cm)
Serrated



Wire Twisting



KTS2480818
7 1/4" (18.5cm)
Tungsten Carbide




*Locking Grip Pliers / Osteotomes /
Retractor / Chiesel*



KTS3353318
6" (15cm)
With Spring
serrated

Locking Pliers



KTS3353217
6 1/2" (17cm)
Grooved, Serrated

Flat Nosed Parallel Pliers



Fig.	Item	Description	Width
A	KTS3265403	5" (13cm) Straight	3mm
A	KTS3265404	5" (13cm) Straight	4mm
A	KTS3265406	5" (13cm) Straight	6mm
A	KTS3265408	5" (13cm) Straight	8mm
A	KTS3265410	5" (13cm) Straight	10mm
A	KTS3265412	5" (13cm) Straight	12mm
A	KTS3265415	5" (13cm) Straight	15mm
A	KTS3265303	6 3/4" (17cm) Straight	3mm
A	KTS3265304	6 3/4" (17cm) Straight	4mm
A	KTS3265305	6 3/4" (17cm) Straight	5mm
A	KTS3265306	6 3/4" (17cm) Straight	6mm
A	KTS3265308	6 3/4" (17cm) Straight	8mm
A	KTS3265310	6 3/4" (17cm) Straight	10mm
A	KTS3265312	6 3/4" (17cm) Straight	12mm
A	KTS3265315	6 3/4" (17cm) Straight	15mm
A	KTS3265319	6 3/4" (17cm) Straight	19mm

Mini Lambotte




Fig.	Item	Description	Width
A	KTS3265706	9" (23cm) Straight	1/4" 6mm
A	KTS3265712	9" (23cm) Straight	1/2" 12mm
A	KTS3265719	9" (23cm) Straight	3/4" 19mm
A	KTS3265725	9" (23cm) Straight	1" 25mm
A	KTS3265732	9" (23cm) Straight	1 1/4" 32mm
A	KTS3265738	9" (23cm) Straight	1 1/2" 38mm
B	KTS3265806	9" (23cm) Curved	1/4" 6mm
B	KTS3265812	9" (23cm) Curved	1/2" 12mm
B	KTS3265819	9" (23cm) Curved	3/4" 19mm
B	KTS3265825	9" (23cm) Curved	1" 25mm
B	KTS3265832	9" (23cm) Curved	1 1/4" 32mm
B	KTS3265838	9" (23cm) Curved	1 1/2" 38mm

Lambotte



Chandler


KTS3276118
7 1/2" (19cm)
13mm
KTS3276120
8" (20cm)
20mm
KTS3276123
9" (23cm)
25mm
KTS3276125
10" (25cm)
31mm



Retractor / Bone Holding Clamps / Elevators Ortho / Mallets



Gerster-Lowmann



1/4

KTS3244120
8" (20cm)
1 x 2 prongs



Key
(Fig. A)

Fig.	Item	Description	Width
A	KTS3275903	7 1/2" (19cm)	3mm
A	KTS3275906	7 1/2" (19cm)	6mm
A	KTS3275909	7 1/2" (19cm)	9mm
A	KTS3275912	7 1/2" (19cm)	12mm
A	KTS3275915	7 1/2" (19cm)	15mm
A	KTS3275919	7 1/2" (19cm)	19mm
A	KTS3275925	7 1/2" (19cm)	25mm



Bennett
KTS3277171
10 1/4" (26cm)
65mm




Bergmann
KTS3274101
9 1/2" (24cm)
10 1/4OZ (300g)




Gerzog
KTS3268625
7" (18cm)
8 1/2OZ (240g)



Ombredanne
KTS3274224
9 1/2" (24cm)
19 1/2OZ (550g)




Cottle
KTS3274624
7 1/2" (19cm)
8 1/2OZ (240g)



Collin
KTS3269222
8" (23cm)
33mm Dia.
7 3/4OZ (220g)



Mallets / Curettes



Hajek

KTS3274822
8 3/4" (22.25cm)
27mm Dia.
5OZ (140g)

Ortho / With Plastic Caps

KTS3277214
7" (18cm)
25mm Dia.
5OZ (140g)


KTS3277215
Pair of replacement
caps for KTS3277214

Partsch

KTS3277620
7" (18cm)
22mm Dia.
7OZ (200g)

Ortho

KTS3269145
10 1/2" (26.5cm)
160 OZ (450g)




Ortho


KTS3269090
10 1/2" (26.5cm)
32 OZ (900g)
KTS3269045
10 1/2" (26.5cm)
16 OZ (450g)



Spratt /
Hollow Handle



Volkman



D

Curettes / Rasp



1/2

Bruns



KTS3247197
6 3/4" (17cm)
000



KTS3247198
6 3/4" (17cm)
00



KTS3247200
6 3/4" (17cm)
0




KTS3247201
6 3/4" (17cm)
1



KTS3247202
6 3/4" (17cm)
2



KTS3247203
6 3/4" (17cm)
3



1/4
KTS3247204
6 3/4" (17cm)
4



KTS3247205
6 3/4" (17cm)
5



KTS3247206
6 3/4" (17cm)
6



KTS3261397
9" (22.5cm)
000



KTS3261398
9" (22.5cm)
00



KTS3261300
9" (22.5cm)
0



KTS3261301
9" (22.5cm)
1



KTS3261302
9" (22.5cm)
2



KTS3261303
9" (22.5cm)
3



1/4
KTS3261304
9" (22.5cm)
4



KTS3261305
9" (22.5cm)
5



KTS3261306
9" (22.5cm)
6



1/3

Hibbs-Spratt /
Angled

D



1/3

Hibbs-Spratt /
Straight



KTS3261297
9" (22.5cm)
000



KTS3261298
9" (22.5cm)
00



KTS3261200
9" (22.5cm)
0



KTS3261201
9" (22.5cm)
1



KTS3261202
9" (22.5cm)
2



KTS3261203
9" (22.5cm)
3



KTS3261204
9" (22.5cm)
4



KTS3261205
9" (22.5cm)
5



KTS3261206
9" (22.5cm)
6



KTS3249600
5 1/4" (13.25cm)
Double Ended
00/0

1/2

Williger



KTS3249213
5" (13cm)
Double Ended
6.2mm & 6.7mm

1/2

Volkman



KTS3249216
6 1/4" (16cm)
Double Ended
7mm & 9mm



KTS3249220
8" (20cm)
Double Ended
9.1mm & 11.1mm



Rasp / Curettes / Rongeurs Single Action



KTS3297327
10 1/2" (27cm)
Double Ended
KTS3297530
12" (30cm)
Double Ended

Putti



KTS9900006
7 1/4" (18cm)
Double Ended

Miller



Friedman &
Mini-Friedman



KTS3210112
5" (13cm)
Curved
13mm x 2.5mm



KTS3210114
5 1/2" (14cm)
Curved
12mm x 3mm



Universal



KTS3215613
5 1/4" (13.25cm)
Angled
17mm x 3.5mm



1/2

Luer




KTS3213515
6" (15.25cm)
Curved-Small
9mm x 3.5mm




1/2

Lempert



KTS3215116
6 1/4" (16cm)
Curved, 10mm x 4.5mm
KTS3215418
7 1/4" (18cm)
Straight, 12mm x 3mm
KTS3215518
7 1/4" (18cm)
Curved, 12mm x 3mm



1/2

Adson

KTS3214220
8" (20cm)
Straight,
15mm x 7.5mm



1/2

Bane




KTS9911050
7" (18cm)
Slightly Curved,
16.5mm x 4mm



D

Rongeurs Double Action



Kleinert-Kutz



Ruskin-Mini




Ruskin



Beyer



Zaufal-Jansen



Semb



Stille




Rongeurs Double Action / Bone Cutters-Single Action



Stille-Ruskin



Leksell-Stille




Stille-Luer




Echlin / Stille-Luer





Liston



Mayfield



Littauer / Liston



Bone Cutters Double Action / Pliers



1/2

Böhler

KTS3221115
6" (15.25cm)
Straight,
5mm x 15mm



KTS3221315
6" (15.25cm)
Slight Curve,
5mm x 15mm

1/2

Walton-Liston



KTS3230015
6" (15.25cm)
Straight,
16mm



KTS3230215
6" (15.25cm)
Angled,
16mm

1/2

Ruskin-Liston



KTS3228618
7 1/4" (18.5cm)
21mm

D



1/3

Stille-Listo

KTS3229027
10 1/2" (27cm)
Straight, 36mm




1/3

Liston

KTS3230523
9 1/4" (23.5cm)
Angled, 28mm

1/2

Slip Joint



KTS3306019
7 3/4" (19.5cm)



1/3

Liston-Key
(Horsley)

KTS3230327
10 1/2" (27cm)
Angled, 28mm




*Bone Hook / Lead Hand / Tendon Passing Pliers
Retractors / Tendon Forceps*



Bone Hook
KTS3283030
9" (23cm)
3cm Dia.



Blount
Knee Retractor
KTS3282717
7" (17.5cm)




Needle Nosed



KTS3390214
5 1/2" (14cm)
Serrated-Groove



Tendon



KTS1299215
6" (15.25cm)
6mm



Lead Hand
KTS3245001
Adult Size



Carroll
KTS2367911
5" (13cm)
Curved Shaft
2mm x 6mm

D



Colipers, Dynamometers, Flexometer



Graduated Tape
Measurer

KTS0235015
60" (1.5cm)
inch & cm markings



Graduated
Metal Ruler

KTS0234015
6" (15cm)
cm & mm
markings



Jameson
KTS0023401



Townley
KTS0023403
11.5cm



Möltgen
KTS0023405



Carroll
KTS0023407
14.5cm



Caliper
KTS0023409



Vernier
KTS0023411
22.5cm



Calipers



Castroviejo
KTS0023413
8cm
Measuring range
15mm



Castroviejo
KTS0023415
8cm
Measuring range
20mm



Castroviejo
KTS0023417
8cm
Measuring range
20mm



Castroviejo-Epker
KTS0023419
9cm
Measuring range
40mm



Castroviejo-Epker
KTS0023421
18cm
Measuring range
40mm




Castroviejo-Epker
KTS0023423
18cm
Measuring range
40mm



E

- **NEUROSURGERY**
- **ENT**

**Category E:**

Ear, Nose and Throat are very delicate parts of human body. All ENT and neuro surgeries require delicate and complicated instruments. We manufacture dissectors, Ronguers, spreaders, Mouth gags and Tongue depressors to meet the surgery needs of these human parts.



Dissector/Laminectomy Rongeurs, Lamina Spreader







Fig.	Item	Description	Tip
A	KTS4051201	7" (18cm) 40° Up	1mm
B	KTS4051202	7" (18cm) 40° Up	2mm
C	KTS4051203	7" (18cm) 40° Up	3mm
D	KTS4051204	7" (18cm) 40° Up	4mm
E	KTS4051205	7" (18cm) 40° Up	5mm
F	KTS4051801	8" (20cm) 40° Up	1mm
G	KTS4051802	8" (18cm) 40° Up	2mm
H	KTS4051803	8" (20cm) 40° Up	3mm
I	KTS4051804	8" (20cm) 40° Up	4mm
J	KTS4051805	8" (20cm) 40° Up	5mm




Laminectomy Rongeurs



KTS4042202
7" (18cm)
Up
3 x 10mm





KTS4042203
7" (18cm)
Down
3 x 10mm



KTS4042001
5" (12.5cm)
Straight
3 x 10mm
KTS4042201
7" (18cm)
Straight
3 x 10mm

1/4



KTS4043201
7" (18cm)
Straight
3 x 10mm
KTS4043301
8" (20cm)
Straight
3 x 10mm
KTS4043401
10" (25cm)
Straight
3 x 10mm


1/4



KTS4043202
7" (18cm)
Up
3 x 10mm



Laminectomy Rongeurs



Spurling

1/2

KTS4044101

6" (15cm)

Straight

4 x 10mm

KTS4044201

7" (18cm)

Straight



4 x 10mm

KTS4044301

8" (20cm)


Straight

4 x 10mm



Schlesinger

1/2



KTS4045602

6" (15cm)

Straight

2 x 10mm

Serrated




KTS4045603

6" (15cm)

Straight

3 x 10mm

Serrated



KTS4045604

6" (15cm)


Straight

4 x 10mm

Serrated





Mouth Gags / Tongue Depressor



E



Tongue Dilator, Hook / Tonsil / Tenaculum / Seizing Forceps



ENT Forceps / Tonsil Snares




KTS5020219
7 1/2" (20cm)
Straight




KTS5020319
7 1/2" (20cm)
Curved



KTS5025228
11 (28cm)
With Ratchet




KTS5025903
Package of 12
0.3mm
KTS5025955
Package of 12
0.55mm



KTS4426108
3 1/2" (8.5cm)
Alligator
2mm x 8mm


KTS4426114
5 1/2" (14cm)
Alligator
2mm x 8mm



KTS5026000
Straight Tip




Micro Ear Scissors / Mirrors / Ear Specula




For adults
KTS4401001
Set 2-4
KTS4401000
Se1-4

Hartmann
KTS4401030 3.0 mm
KTS4401035 3.5 mm
KTS4401040 4.0 mm
KTS4401050 5.0 mm
KTS4401065 6.5 mm
KTS4401075 7.5 mm



Micro-Ear Scissors



KTS4431040
3 1/4" (8cm)
4mm



Laryngeal Mirror
with Handle

KTS5227208
00 8mm

KTS5227210
0 10mm

KTS5227212
1 12mm

KTS5227214
2 14mm

KTS5227216
3 16mm

KTS5227218
4 18mm

KTS5227220
5 20mm

KTS5227222
6 22mm


KTS5227224
7 24mm

KTS5227226
8 26mm


E




Grasper / Suction / Tuning Forks



KTS4499401
3 1/4" (8cm)
Straight-Serrated
4 x 0.6mm




KTS4499601
3 1/4" (8cm)
Straight-Serrated
6 x 0.8mm





KTS4499602
3 1/4" (8cm)
Straight-Serrated
8.5 x 1.8mm



Micro-Ear



KTS4499801
3 1/4" (8cm)
Straight-Serrated
4.5 x 0.8mm



KTS4499802
3 1/4" (8cm)
Straight-Serrated
4.5 x 0.8mm

E



KTS0440310
3" (7.5cm)
With Stylus
1.0mm



KTS0440315
3" (7.5cm)
With Stylus
1.5mm



KTS0440320
3" (7.5cm)
With Stylus
2.0mm



Micro ENT/Bellucci



KTS4407001
C 128

Hartmann/French



Tuning Forks / Ear Syringes



Aluminum / Tuning Forks



KTS4406001
C 128



KTS4406002
C 256



KTS4406007
C 64 Low F



KTS4407100
C128-C1256-C2512
C31024-C42048

Hartmann/Set of 5

Reiner



KTS4409601
Tip Only



KTS4409602
Tip Only



KTS4409603
Shield

Reiner



KTS4409650
50CC Complete
With Tip & Shield

KTS4409710
100CC Complete
With Tip & Shield



Nasal Specula / Nasal Cutting Forceps




1/3

Killian




1/3

Modified Wien




1/4

KTS4611201
5 1/4" (13.25cm)
1" (2.5cm) Tip




KTS4611202
5 1/4" (13.25cm)
1 3/8" (3cm) Tip

KTS4612035
5" (13cm)
35mm




1/4



KTS4612050
5" (13cm)
50mm



KTS4612075
5" (13cm)
75mm




KTS4612090
5" (13cm)
90mm



1/3

Takahashi




1/4

KTS4631301
5" (13cm)
Straight
5 x 10mm




Nasal Cutting Forceps / Nasal Knives / Nasal Periosteal Elevator; Nasal Rasp



KTS4637415
6" (15cm)
8 x 8mm
blade size



Freer "D" Shaped



KTS4632319
7 1/2" (19cm)
Double Action
5 x 15mm tip



Jansen-Middleton



KTS4632519
7 1/2" (19cm)
Double Action
4 x 20mm tip



Jansen-Middleton



KTS4639004
8" (20cm)
Swivel
4mm




Ballenger




KTS4642516
6 1/4" (16cm)
4.5mm




Joseph



KTS4641518
7" (17.75cm)
Double Ended
Sharp/Blunt
4.5mm & 4.5mm
blade size



Freer




KTS4642322
8 1/2" (21cm)
Double Ended
7.5mm & 2.2mm




Cottle



Nasal Elevator / Nasal Rasp



Lewis



Maltz



E



KTS4652021
8 1/4" (21cm)
Double Ended



Fomon




KTS4652209
8 1/4" (21cm)
Double Ended
Tungsten Carbide



KTS4652207
8 1/4" (21cm)
Double Ended
Tungsten Carbide



KTS4652205
8 1/4" (21cm)
Double Ended
Tungsten Carbide



Nasal Forceps / ENT Morselizer / Nasal Fracture Elevator




Asch

KTS4657123
9" (23cm)



Rubin



KTS4625121
8" (20cm)



Boies

KTS4630418
7" (18 cm)



Cottle Bone




KTS4625800
Cottle Cartilage Crusher



F

- Cardiovascular
- Suction Tubes
- Intestinal / Rectal



F

Category F:

Cardio Vascular surgeries are generally considered high risk surgeries. We very responsibly manufacture rib spreaders, vascular clamps, shunt clamps, suction tubes, tissue seizing forceps and many other cardio vascular instruments.

Rib Spreader / Peripheral Vascular Clamps



Finnocchietto/Baby

KTS5621475
3" (7.5cm)
spread
18mm x 20mm x 75mm



Finochietto
KTS5621020
Finnocchietto
7 3/4" (20cm)
spread-Curved Arms
46mm x 55mm x 200mm
KTS5621025
10" (25cm)
spread-Curved Arms
70mm x 65mm x 250mm



Burford / Aluminum

KTS5624218
7" (18cm)
spread with 2 sets of Blades
65mm x 62mm
47mm x 62mm



KTS1342114
5 1/2" (14cm)
Angled, 32mm


Cooley / Pediatric



Vascular Clamps / Vena Cava Forceps / Shunt Clamps



1/2
Satinsky




1/4
KTS1325320
8" (20.5cm)
DeBakey
35 x 11mm

KTS1325324
9 3/4" (24cm)
DeBakey
45 x 15mm

KTS1325325
10" (25cm)
DeBakey
55 x 16mm


KTS1325326
10 1/2" (26.5cm)
DeBakey
65 x 16mm

KTS1325327
10 3/4" (27.5cm)
DeBakey
75 x 16mm



1/3
De Baky / Pediatric

1/4
KTS1311112
4 3/4" (12cm)
Angled, 43mm tip



1/2
Baby-Satinsky

1/4
KTS5632515
6" (15cm)
Cross serrated



1/2
Javid

1/4
KTS1334318
7" (18cm)
Small,
4mm Ring



Suction Tubes / Adapters



Yankauer

KTS0415027
11" (27cm)



1/4



Andrews-Pynchon

KTS0415222
9 1/4" (23.5cm)

1/4



Poole

KTS0413010
8 1/4" (21cm)
Straight, 30FR



1/3

Poole

KTS0414008
8 1/4" (21cm)
Curved, 23FR



Luer

KTS2238200
Luer lock adapter



1/3


Frazier




Fig.	Item	Description	W.L	Tip
A	KTS4029506	7" (18cm) Angled-30°	4.75"	6 French
B	KTS4029507	7" (18cm) Angled-30°	4.75"	7 French
C	KTS4029508	7" (18cm) Angled-30°	4.75"	8 French
D	KTS4029509	7" (18cm) Angled-30°	4.75"	9 French
E	KTS4029510	7" (18cm) Angled-30°	4.75"	10 French
F	KTS4029511	7" (18cm) Angled-30°	4.75"	11 French
G	KTS4029512	7" (18cm) Angled-30°	4.75"	12 French
H	KTS4029514	7" (18cm) Angled-30°	4.75"	14 French
I	KTS4029710	7" (18cm) Angled-75°	4.75"	10 French
J	KTS4029712	7" (18cm) Angled-75°	4.75"	12 French



Suction Tubes / Adapters




Tracheal Tubes, Speaking Tubes



Trachea Tube

outer tube window
inner tube window
KTS0088500



Fig.	Item	Ø	L
A	KTS0088501	7mm	60mm
A	KTS0088502	8mm	65mm
A	KTS0088503	9mm	70mm
A	KTS0088504	10mm	75mm
A	KTS0088505	11mm	80mm
A	KTS0088506	12mm	85mm
A	KTS0088507	13mm	90mm



Speaking Tube

outer tube closed
inner tube closed
KTS0088400



Fig.	Item	Ø	L	Radius
A	KTS0088401	8mm	95mm	50mm
A	KTS0088402	9mm	95mm	50mm
A	KTS0088403	10mm	95mm	50mm
A	KTS0088404	11mm	95mm	50mm
A	KTS0088405	12mm	95mm	50mm
A	KTS0088406	13mm	95mm	50mm



with sliding valve

F

Tissue Seizing Forceps



F



Rectal Speculas / Stomach & Intestinal Clamp Forceps



Babcock

KTS6416016
6 1/4" (16cm)
8mm

1/4




KTS6416018
7" (18cm)
10mm



KTS6416020
8" (20cm)
13mm



KTS6416024
9 1/2" (24cm)
14mm



Doyen
KTS6421123
8 3/4" (22cm)
Straight, grooved

KTS6421023
8 3/4" (22cm)
Curved, grooved

1/4

Rectal Speculas



Ives-Fansler
KTS6821022
3 1/4" x 7/8
80 x 22mm

1/4




Pratt
KTS6812100
8 1/2" (21.5cm)
85 x 25mm blade

1/3


Non-Traumatic Stomach & Intestinal Clamp Forceps



Kocher
22.0cm 8 3/4"
25.0cm 10"
28.0cm 11"



Brunner
KTS6426024
24.0cm 9 1/2"



Resano
KTS6429330
30.0cm 12"


F




Stomach & Intestinal Clamp Forceps



$\frac{1}{2}$
Kocher



$\frac{1}{4}$
KTS6422022
8 3/4" (22cm)
Straight



$\frac{1}{3}$
Fehland
KTS6424723
24.0cm 9 1/2"

F



Intestinal & Stomach Crushing Clamps, Anastomosis Clamps



Payr-Baby
15.0cm 6"



KTS6431015 without pin
KTS6431215 with pin




Payr



KTS6431021
21.0cm 8 1/4"
without pin
KTS6431221
21.0cm 8 1/4"
with pin



KTS6431029
29cm 11 1/2"
without pin
KTS6431229
29cm 11 1/2"
with pin



KTS6431035
35cm 13 3/4"
without pin
KTS6431235
35cm 13 3/4"
with pin



Rectal Specula



De Martel-Wolfson
KTS6428226
26.0cm 10 1/4"



Mathieu
KTS6813019
19.5cm 7 3/4"
90 x 15mm



Sims
KTS6811515
15.0cm 6"




Sims
KTS6811015
15.0cm 6"



G

- **Urology**
- **Gynecology**
- **Obstetric**

**Category G:**

This section covers catheters, dilators, bougies, Stone-Forceps Kidney, Vaginal speculum / retractors, and cervical forceps for gynecology and urology surgical needs.

Catheter / Dilators / Bougies / Stone Forceps-Kidney








Fig.	Item	Description	Dia.
A	KTS6620108	11" (28cm)	8 French
B	KTS6620110	11" (28cm)	10 French
C	KTS6620112	11" (28cm)	12 French
D	KTS6620114	11" (28cm)	14 French
E	KTS6620116	11" (28cm)	16 French
F	KTS6620118	11" (28cm)	18 French
G	KTS6620120	11" (28cm)	20 French
H	KTS6620122	11" (28cm)	22 French
I	KTS6620124	11" (28cm)	24 French
J	KTS6620126	11" (28cm)	26 French
K	KTS6620128	11" (28cm)	28 French
L	KTS6620130	11" (28cm)	30 French
M	KTS6620132	11" (28cm)	32 French
N	KTS6620134	11" (28cm)	34 French
O	KTS6620136	11" (28cm)	36 French




Vaginal Retractor / Vaginal Speculum



De Vilbiss

KTS9930141
4 x 1 1/4"
102mm x 32mm
Medium



Graves

KTS7011001
3 X 3/4" (75 x 20mm)
Small
KTS7011002
3 3/4" X 1 3/4" (95 x 35mm)
Medium
KTS7011003
4 1/2" X 1 3/8" (115 x 35mm)
Large



Pederson

KTS7012001
3 X 1/2" (75 x 13mm)
Small
KTS7012002
4" X 7/8" (100 x 22mm)
Medium
KTS7012003
4 3/4" X 1" (120 x 25mm)
Large



Auvard


KTS7030800
Long Blade
(102 x 44mm)
KTS7030900G
Fixed Weight
(86 x 32mm)




Tenaculum Forceps / Cervical Forceps



Kogan Gyne



Hulka Gyne



Schröder Gyne




Schröder



Pozzi Gyne



Uterine Dilators Sets



Pratt Gyne



1/3

Pratt Gyne

KTS7044708

Set of 8 Double ended
size 13/14 thru 41/43
4.5/5.0mm thru 13.6/14.3mm

KTS7043913	Fig. 13/15	4.5/5.0 mm
KTS7043917	Fig. 17/19	5.6/6.3 mm
KTS7043921	Fig. 21/23	7.0/7.6 mm
KTS7043925	Fig. 25/27	8.3/9.0 mm
KTS7044729	Fig. 29/31	9.6/10.3 mm
KTS7044733	Fig. 33/35	11.0/11.6 mm
KTS7044737	Fig. 37/39	12.3/13.0 mm
KTS7044741	Fig. 41/43	13.6/14.3 mm



1/3


Hegar

KTS7042608


Set of 8 dilator
3.0/4.0 mm - 17.0/18.0 mm

KTS7042613

Set of 13 dilator
1.0/2.0 mm - 25.0/26.0 mm



Hank



1/3

Hank Gyne

KTS7043608
Set of 6 dilators
4.5/5.0 mm - 9.5/10.0 mm

KTS7043010	Fig. 9/10	4.5/5.0 mm
KTS7043012	Fig. 11/12	5.5/6.0 mm
KTS7043014	Fig. 13/14	6.5/7.0 mm
KTS7043016	Fig. 15/16	7.5/8.0 mm
KTS7043018	Fig. 17/18	8.5/9.0 mm
KTS7043020	Fig. 19/20	9.5/10.0 mm




Uterine Dilators sets



KTS7042408
Set of 8 dilators
3.0/4.0 mm - 17.0/18.0 mm
KTS7042413
Set of 13 dilators
1.0/2.0 mm - 25.0/26.0 mm

Hegar Gyne

KTS7042001	1.0/2.0 mm
KTS7042003	3.0/4.0 mm
KTS7042005	5.0/6.0 mm
KTS7042007	7.0/8.0 mm
KTS7042009	9.0/10.0 mm
KTS7042011	11.0/12.0 mm
KTS7042013	13.0/14.0 mm
KTS7042015	15.0/16.0 mm
KTS7042017	17.0/18.0 mm
KTS7042019	19.0/20.0 mm
KTS7042021	21.0/22.0 mm
KTS7042023	23.0/24.0 mm
KTS7042025	25.0/26.0 mm



Hegar Gyne
KTS7042508
Set of 8 dilators
3.0/4.0 mm - 17.0/18.0 mm



KTS7041001	1.0 mm
KTS7041002	2.0 mm
KTS7041003	3.0 mm
KTS7041004	4.0 mm
KTS7041005	5.0 mm
KTS7041006	6.0 mm
KTS7041007	7.0 mm
KTS7041008	8.0 mm
KTS7041009	9.0 mm
KTS7041010	10. mm
KTS7041011	11.0 mm
KTS7041012	12.0 mm
KTS7041013	13.0 mm
KTS7041014	14.0 mm
KTS7041015	15.0 mm
KTS7041016	16.0 mm
KTS7041017	17.0 mm
KTS7041018	18.0 mm
KTS7041019	19.0 mm
KTS7041020	20.0 mm
KTS7041021	21.0 mm
KTS7041022	22.0 mm
KTS7041023	23.0 mm
KTS7041024	24.0 mm
KTS7041025	25.0 mm
KTS7041026	26.0 mm
KTS7041027	27.0 mm
KTS7041028	28.0 mm
KTS7041029	29.0 mm
KTS7041030	30.0 mm

KTS7040414
Set of 14 dilators
4.0 mm - 17.0 mm

KTS7040426
Set of 26 dilators
1.0 mm - 26.0 mm

KTS7041501	1.5 mm
KTS7041502	2.5 mm
KTS7041503	3.5 mm
KTS7041504	4.5 mm
KTS7041505	5.5 mm
KTS7041506	6.5 mm
KTS7041507	7.5 mm
KTS7041508	8.5 mm
KTS7041509	9.5 mm
KTS7041510	10.5 mm
KTS7041511	11.5 mm
KTS7041512	12.5 mm
KTS7041513	13.5 mm
KTS7041514	14.5 mm
KTS7041515	15.5 mm
KTS7041516	16.5 mm
KTS7041517	17.5 mm
KTS7041518	18.5 mm
KTS7041519	19.5 mm



Hegar Gyne

Hegar Gyne
KTS7040614
Set of 14 dilators
4.0 mm - 17.0 mm



Tenaculum / Curettes / Biopsy Forceps



KTS7058132
12 1/2" (32cm)
Graduated-Malleable

1/4
Sims



KTS7050425
10" (25cm)

1/3
Braun



1/3
Kevorkian



KTS7059412
12" (30cm)
1.2mm



KTS7059440
12" (30cm)
4.0 mm
Without Basket



KTS7059640
12" (30cm)
4.0mm
With Basket



1/3
Recamier /
Sims



KTS7071508
10 1/2" (26cm)
Malleable-sharp
8mm




KTS7071509
10 1/2" (26cm)
Malleable-sharp
9mm



KTS7071511
10 1/2" (26cm)
Malleable-sharp
11mm




Curettes / Biopsy Forceps / Cervical



Tischler-Kevorkian

KTS7062700
8 1/2" (22cm)
with latch
3 x 9.5mm




Tischler


KTS7063101
8 1/2" (22cm)



Biopsy Forceps / Cervical



KTS7063102
8" (20cm)
With Fixing Screw
2.3 x 4.2mm
KTS7063103
8" (20cm)
Without Fixing Screw
2.3 x 4.2mm



KTS7062101
8 1/2" (22cm)
3 x 7mm



KTS7061723
9" (23cm)
3 x 9.7mm
KTS7061704
10" (25cm)
3 x 9.7mm



Vesiculum Forceps / Delivery Forceps / Tenaculum Forceps



Jacobs
KTS7053221
8 1/4" (21cm)
Straight



Iowa
KTS7053926
10 1/4" (26cm)



Simpson
KTS7210230
12" (30cm)



Piper
KTS7215544
18" (44cm)



H

• *Titanium***Category H:**

In this category we have included the most valuable Titanium Surgical Instruments. You will find Hartmann, Bengolea, Halstead, Crile, Hesis, Doyen, Mayo-Hegar, Vannas, Jacobson, Potts, and other Titanium instruments in different sizes both straight & curved.





KTS7053222
4 3/4" (12cm)

ADSON



KTS7053223
4 3/4" (12cm)



Mc IND OE
KTS7053224
6" (15cm)



GILLIES
KTS7053225
6" (15cm)



CUSHING
KTS7053226
7 1/8" (18cm)



CUSHING
KTS7053227
8" (20cm)





BARRAYA
KTS7053228
8" (20cm)



GUILMET
KTS7053229
8" (20cm)








KTS7053242
8 1/2" (21.5cm)

O' SHAUGHNESSY



KTS7053243
8 1/4" (21cm)

HEISS



O' SHAUGHNESSY
KTS7053244/KTS7053253

KTS7053244
7 1/4" (18cm)

KTS7053245
7 1/4" (18cm)

KTS7053246
8" (20cm)

KTS7053247
8 3/4" (22cm)

KTS7053248
8" (20cm)

KTS7053249
8 3/4" (22cm)

KTS7053250
9 1/2" (24cm)






KTS7053251
9 1/2" (24cm)

KTS7053252
11" (28cm)





KTS7053253
11" (28cm)



H



Pinces a fil pour micro suture, a manches ronds
draadvoerders voor micro suture, met ronde grepen
micro suture tying forceps, with round handles
pinze a refe per micro suture con manici rotonde
pinzas a hilos para micro suture con mangos redondas



Micro ciseaux a manches plants
micro scharen met vlakke grepen
micro scissors with flat handles
micro forbici con manici plane
micro tijeras con mangos planos



YASARGIL
KTS7053287
7 1/4" (18.5cm)



KTS7053288



KTS7053289



KTS7053290



KTS7053291



KTS7053292



KTS7053293



POTTS
7 1/4" (18.5cm)




25° **KTS7053294**



45° **KTS7053295**



60° **KTS7053296**



90° **KTS7053297**



130° **KTS7053298**



KTS7053299



KTS7053300



KTS7053301




KTS7053302


4 3/4" (12cm)

4 3/4" (12cm)


Micro ciseaux a manches plans
micro scharen met vlakke grepen
micro scissors with flat handles
micro forbici con manici plane
micro tijeras con mangos planos




6" (15cm)




6" (15cm)




7 1/4" (18.5cm)



7 1/4" (18.5cm)




6 1/4" (16cm)



6 1/4" (16cm)



7 1/4" (18.5cm)




7 1/4" (18.5cm)






ecarteur avec quatre valves, pour
thorcoplastie en chirurgie neonatale
wondhaken met vier spatels, voor
thoracoplastie in neonatale chirurgie
retractor, with four blades, for
thoracoplasty in neonatal surgery
divaricatore, con quattro valve, per
thoracoplastia en neonatle chirurgia
separador, con cuatro valvas, para
thoracoplastia en neonatal cirugia



$\frac{3}{4}$


DENIS BROWNE
KTS7053319
6" (15cm)



KTS

$\frac{1}{2}$

DENIS BROWNE
KTS7053320
8" (20cm)



En chirurgie enfant
in kinder chirurgie
en bambino chirurgia
en nino cirugia

FINOCHIETTO
KTS7053321

$\frac{1}{2}$

ecarteur avec quatre valves, pour
thorcoplastie en chirurgie infantile
wondhaken met vier spatels, voor
thoracoplastie in kinderchirurgie
retractor, with four blades, for
thoracoplasty in infant surgery
divaricatore, con quattro valve, per
thoracoplastia en chirurgia pediatrica
separador, con cuatro valvas, para
thoracoplastia en cirugia infantil

Surgical Instruments Sets

QTY .	DESCRIPTION	QTY .	DESCRIPTION
MAJOR BASIC SET			
2 each	Mayo Dissecting Scissors, Straight. 6 3/4"	1 each	Dressing Forceps, Serrated 5 1/2"
1 each	Mayo Dissecting Scissors, Curved. 5 1/2"	1 each	Tissue Forceps 1 X 2 Teeth 4 1/2"
2 each	Mayo Dissecting Scissors, Curved. 6 3/4"	1 each	Tissue Forceps 1 X 2 Teeth 5 1/2"
12 each	Halsted Mosquito Forceps, Curved. 5"	2 each	Crile-wood T.c. Needle Holders 6"
6 each	Kelly Forceps, Straight. 5 1/2"	2 each	Mayo-hegar T.c.needle Holders 7"
18 each	Kelly Forceps, Curved. 5 1/2"	1 each	Metzenbaum Scissors, Curved 5 3/4"
18 each	Rochester-pen Forceps, Curved 6 1/4"	INTESTINAL INSTRUMENT SET	
6 each	Rochester-pen Forceps, Curved 7 1/4"	<i>(Used With Major Basic Set)</i>	
12 each	Rochester-ochsner Forceps, Straight 6 1/4"	1 each	Dressing Forceps 10"
2 each	Dressing Forceps, Serrated 5 1/2"	2 each	Mayo-hegar Neeldle Holder 7"
1 each	Dressing Forceps, Serrated 8"	1 each	Payr Pylorus Clamp. 5 3/4"
2 each	Tissue Forceps 1 X 2 Teeth 5 1/2"	6 each	Percy Intestinal Forceps. 6"
1 each	Tissue Forceps 1 X 2 Teeth 8"	6 each	Babcock Intestinal Forceps 7 3/4"
2 each	Crile-wood T.c. Needle Holders. 6"	3 each	Allis Tissu Forceps, 5 X 6 Teeth. 7 3/4"
2 each	Mayo-hegar T.c.needle Holders. 7"	2 each	Pean G.i. Forceps, Straight. 9"
2 each	Mayo-hegar T.c.needle Holders. 8"	2 each	Pean G.i. Forceps, Curved. 9"
1 each	Metzenbaum Scissors, Curved 5 3/4"	2 each	Allen Intestinal Clamps, Straight. 8"
1 each	Wire Scissors, Angular 4 3/4"	2 each	Doyen Forceps, Straight. 8 3/4"
1 each	Director And Tongue Tie, Probe End. 6"	2 each	Doyen Forceps.curved. 8 3/4"
1 each	Probe, W/ Eye 6"	1 each	Payr Pylorus Clamp. 8"
1 each	Ochsner Trocar, 20 Fr	1 each	Payr Pylorus Clamp. 11"
1 each	Poole Suction Tube, Angular, 23 Fr	1 each	Payr Pylorus Clamp. 13 3/4"
4 each	Poole Suction Tube, Straight, 30 Fr	1 each	Deaver Retractor. 2" X 12"
4 each	Babcock Forceps, 6 1/4"	2 each	Sponge Forceps, Straight, Serrated. 9 1/2"
12 each	Babcock Forceps, 7 3/4"	CAESAREAN SECTION SET	
4 each	Allis Forceps, 4 X 5 Teeth 6"	<i>(Used With Major Basic Set)</i>	
1 each	Allis Forceps, 6 X 7 Teeth 7 1/2"	1 each	Bandage Scissors. 7 1/4"
1 each	U.s. Army Retractor, Double Ended, Set Of 2.	8 each	Pennington Clamps, Straight. 5 3/4"
2 each	Parker Retractors, D.e, Set Of 2 . . . 5/8" X 5"	1 each	Delee O.b. Forceps, Short. 12"
6 each	Kelly Retractors, Medium	GALL BLADDER SET	
1 each	Richarson Retractors, Medium	<i>(Used With Major Basic Set)</i>	
1 each	Richarson-eastman Retractor, Small	2 each	Crile- Wood Needle Holder. 7"
1 each	Richarson-eastman Retractor, Large	1 each	Mayo Cystic Duct Scoop
1 each	Deaver Retractor 1" X 12"	1 each	Mayo Common Duct Scoop, Small
1 each	Deaver Retractor 2" X 12"	1 each	Mayo Common Duct Scoop, Medium
2 each	Deaver Retractor 3" X 12"	1 each	Mayo Common Duct Scoop, Large
2 each	Malleable Retractor 1 1/2" X 13"	4 each	Bakes Common Duct Dilators, 3mm-11 mm, Set Of 9
2 each	Volkman Retractor, Sharp Points, 4 Prog	4 each	Johns Hopkins Gall Bladder Forceps. 8"
1 each	Balfour Retractor	6 each	Lahey Gall Duct Forceps. 7 1/2"
12 each	Metzenbaum Scissors, Curved 7"	1 each	Lahey Gall Duct Forceps. 9"
8 each	Deaver Retractor. 1" X 13"	1 each	Grey Cystic Duct Forceps, Without Teeth, Set Of 2
1 each	Backhaus Towel Clamp 5 1/4"	1 each	Grey-cystic Duct Forceps, With Teeth, Set Of 2
1 each	Sponge Forceps, Straight, Serrated 9 1/2"	1 each	Blake Gall Stone Forceps. Curved. 8"
1 each	Deaver Retractor 1 1/2" X 12"	1 each	Desjardin Gall Stone Forceps. 9 1/2"
1 each	Deaver Retractor 1" X 9 12"	1 each	Ochsner Trocar. Size 16 (for 12 Fr. Catheter)
MAJOR BASIC SET		1 each	Ochsner Trocar.size 20 (for 16 Fr. Catheter)
1 each	Mayo Dissecting Scissors, Straight. 5"	1 each	Ochsner Trocar, Size 22 (for 18 Fr. Catheter)
1 each	Mayo Dissecting Scissors, Curved. 5 1/2"	3 each	Poole Suction Tube. Angular, 23 Fr.
12 each	Halsted Mosquito Forceps, Curved 5"		
6 each	Kelly Forceps, Straight 5 1/2"		
12 each	Kelly Forceps, Curved. 5 1/2"		
6 each	Rochester-ochsner Forceps, Straight 6 1/4"		
1 each	Dressing Forceps, Serrated 4 1/2"		



Surgical Instruments Sets

QTY .	DESCRIPTION	QTY .	DESCRIPTION
MAJOR BASIC SET			
8 each	Allis Forceps, 5*6 Teeth. 7 3/4"	2 each	Probes, W/eye, Silver 6"
1 each	Pean Forceps, Curved. 9"	12 each	Allis Forceps. 5x6 Teeth 7 3/4"
2 each	Nelson-metzenbaum Scissors, Curved 9"	1 each	Stille-liston Bone Forceps, Straight. 10 1/2"
4 each	Sponge Forceps. Straight. Serrated 9 1/2"	1 each	Stille-luer Bone Rongeur, Curved. 8 1/2"
1 each	Peanut Sponge Forceps, Straight, Serrated. 9 1/2"	1 each	Stille-luer Bone Rongeur, Angular. 9"
THYROIDECTOMY SET			
24 each	Providence Hospital Forceps, Straight. 5 1/2"	1 each	Gerzog Bone Mallet
24 each	Providence Hospital Forceps, Curved. 5 1/2"	1 each	Army Pattern Osteotome,6mm
6 each	Lahey Goiter Forceps. 6 1/4"	1 each	Army Pattern Osteotome,12mm
4 each	Mastin Muscle Clamps 6 3/4"	1 each	Army Pattern Osteotome,18mm
2 each	Lahey Retractor 7 3/4"	1 each	Langenbeck Elevator, Narrow 7 1/2"
1 each	Green Retractor 8 1/2"	1 each	Langenbeck Elevator, Broad 7 1/2"
VAGOTOMY SET			
2 each	Dressing Forceps, Serrated 12"	1 each	Sayre Elevator, Blunt, Double Ended
2 each	Masson Needle Holders. T.C Jaw 10 1/2"	1 each	Cushing Decompression Retractor. 8"
2 each	Lower Gall Duct Forceps. 7"	1 each	Cushing Nerve And Vein Retractor. 9"
2 each	Lahey Gall Duct Forceps. 7 1/2"	1 each	Weitlaner Retractor, Sharp 6 1/2"
6 each	Allis Forceps, 5*6 Teeth. 6"	24 each	Dandy Scalp Hemostats. 5 1/2"
1 each	Mayo Dissecting Scissors. Curved 9"	6 each	Adson Forceps, Straight. 7 1/4"
1 each	Nelson-metzenbaum Thoracic Scissors. Cvd. 9"	6 each	Adson Forceps, Curved 7 1/4"
1 each	Nelson-metzenbaum Thoracic Scissors. Str 9"	1 each	Adson Dura Hook, Sharp
1 each	Parker-mott Retractors. Double Ended, Set of 2	1 each	Frazier Brain Suction Tube, 6fr.
2 each	Deaver Retractor 2"*12"	4 each	Gigli Saw, Stille Pattern. 20"
2 each	Green Retractor 8 1/2"	1 each	Handles, Solid, For Gigli Saws, In Pairs
2 each	Sponge Forceps. Straight. Serrated. 9 1/2"	CARDIAC ARREST SET	
NEUROSURGICAL LAMINECTOMY INSTRUMENTS			
1 each	Cushing Brain Forceps, Serrated. 7"	1 each	Mayo Dissecting Scissors, Straight. 6 3/4"
1 each	Cushing Brain Forceps. W/mouse Teeth 7"	1 each	Mayo Dissecting Scissors, Curved. 6 3/4"
2 each	Adson Brain Forceps, W/mouse Teeth. 4 3/4"	6 each	Halstead Mosquito Forceps, Straight. 5"
2 each	Adson Brain Forceps, Serrated. 4 3/4"	6 each	Rochester-pean Forceps, Curved. 6 1/4"
1 each	Adson Forceps. Bayonet. Serrated. 7 1/2"	1 each	Dressing Forceps. 6"
1 each	Smithwick Clip Applying Forceps 10 1/2"	1 each	Tissue Forceps, 1x2 Teeth. 4 1/2"
1 each	Sweet Clip Applying Forceps. 9 1/4"	1 each	Mayo-heger Needle Holder, T.c. Jaw. 7"
1 each	Sweet Clip Applying Forceps 11"	2 each	Allis Intestinal Forceps, 4x5 Teeth. 6"
1 each	Strully Scissors, Curved. 8 3/4"	4 each	Backhaus Towel Clamp. 5 1/2"
2 each	Senn Retractors, Double Ended. Sharp. 6 1/4"	2 each	Sponge Forceps, Straight, Serrated. 9 1/2"
10 each	Backhaus Towel Clamp. 5 1/2"	THORACIC AND CARDIOVASCULAR INSTRUMENTS	
8 each	Sponge Forceps, Straight, Serrated. 9 1/2"	2 each	Yankauer Suction Tubes
1 each	Yabkauer Suction Tube	2 each	Mayo Dissecting Scissors, Straight. 6 3/4"
2 each	Operating Scissors, Straight, Sharp/blunt. 5 1/2"	2 each	Mayo Dissecting Scissors, Curved. 6 3/4"
2 each	Mayo Dissecting Scissors, Straight 5 1/2"	36 each	Halstead Mosquito Forceps, Straight. 5"
2 each	Mayo Dissecting Scissors, Straight. 6 3/4"	36 each	Halstead Mosquito Forceps, Curved.. . . . 5"
2 each	Mayo Dissecting Scissors, Curved 5 1/2"	36 each	Rochester-pean Forceps, Curved. 6 1/4"
2 each	Mayo Dissecting Scissors, Curved 6 3/4"	18 each	Rochester-ochsner Forceps, Straight. 6 1/4"
24 each	Crile Forceps, Curved 5 1/2"	12 each	Rochester-ochsner Forceps, Straight 8"
2 each	Crile-wood Needle Holders, T.c. Jaw 6"	36 each	Crile Forceps, Curved. 5 1/2"
1 each	Suture Wire Cutting Scissors, Angular	2 each	Dressing Forceps, Serrated. 5 1/2"
1 each	Grooved Director. 6"	2 each	Dressing Forceps, Serrated. 8"
		2 each	Dressing Forceps, Serrated. 10"
		2 each	Tissue Forceps 1x2 Teeth. 5 1/2"
		1 each	Russian Tissue Forceps. 6"
		2 each	Russian Tissue Forceps. 8"
		2 each	Crile-wood Needle Holder. 7"
		2 each	Mayo-hegar Needle Holder. 7"
		1 each	Hegar Needle Holder, T.c. Jaw. 8"



Surgical Instruments Sets

QTY .	DESCRIPTION	QTY .	DESCRIPTION
2 each	Masson Needle Holder, T.c. Jaw. 10 1/2"	3 each	Blades For Laforce Adenotome, Small
1 each	Suture Wire Cutting Scissors, Angular	1 each	Barnhill Adenoid Curette, Size 1
2 each	Grooved Director. 6"	1 each	Barnhill Adenoid Curette, Size 2
2 each	Probes W/eye, Silver. 6"	2 each	Laryngeal Mirrors, W/handle, Size 4
1 each	Ochsner Trocar, Size 20	1 each	Denhardt Mouth Gag
6 each	Babcock Forceps. 6 1/4"	1 each	Jennings Mouth Gag, Large
10 each	Babcock Forceps. 7 3/4"	1 each	Jennings Mouth Gag, Small
2 each	Babcock Forceps. 9 1/4"	1 each	Wieder Tongue Depressor, Small
12 each	Allis Forceps, 5x6 Teeth. 6"	3 each	Foerster Sponge Forceps, Serrated, Cvd. 7"
10 each	Allis Forceps, 5x6 Teeth. 10"	1 each	White Tonsil Forceps, Curved. 9"
6 each	Pean Forceps, Straight. 9"	2 each	Schmidt Tonsil Forceps, Half Curved Jaw. 7 1/2"
1 each	Stille-liston Bone Forceps, Curved. 10 1/2"	1 each	Boettcher Tonsil Scissors, Curved. 7"
1 each	Stille-horsley Bone Forceps, Bayonet. 10 1/2"	1 each	Metzenbaum Scissors, Curved. 7"
1 each	Stille-luer Bone Rongeur, Curved. 8 1/2"	1 each	Hurd Dissector And Retractor
1 each	Potts-smith Scissors, 45 Degree Angle. 7 1/2"	1 each	Abraham Double Edged Tonsil Kinfe 8"
1 each	Potts-smith Scissors, 60 Degree Angle. 7 1/4"	1 each	Fisher Tonsil Knife, Serrated Edge. 9"
2 each	Potts-smith Dressing Forceps, serrated. 9 1/2"	1 each	Beck- Schenck Tonsil Snare, Complete
1 each	Potts-smith Tissue Forceps, 1x2 Teeth. 7"	1 each	Tyding Tonsil Snare, Complete
1 each	Bethune Rib Shears. 14"	2 each	Backhaus Towel Clamp. 5 1/4"
1 each	Stille Rib Shears. 8 1/2"		
1 each	Lebsche Sternum Knife. 10"		
2 each	Mixer Clamps, Curved. 1 1/2"		

THORACIC AND CARDIOVASCULAR INSTRUMENTS

6 each	Duval Lung Grasping Forceps, 1"jaw. 8"
1 each	Matson Rib Stripper. 8 3/4"
1 each	Alexander Costal Periosteotome. 8 1/2"
1 each	Doyen Rib Raspatory And Elevator, Right. 6"
1 each	Doyen Rib Raspatory And Elevator, Left. 6"
1 each	Harrington Splanchnic Retractor, Large. 12 1/2"
2 each	Hibbs Retractor, Blade. 3" X 1"
1 each	Deaver Retractor. 1" X 12"
1 each	Deaver Retractor. 1" X 13"
2 each	Deaver Retractor. 1 1/2" X 12"
2 each	Ribbon Malleable Retractor. 1 1/2" X 13"
2 each	Volkman Rake Retractor, 4 Spong Sharp
18 each	Backhaus Towel Clamp. 5 1/2"
18 each	Sponge Forceps, Serrated, Straight. 9 1/2"
2 each	Deaver Retractor. 2" X 12"
2 each	Nelson-metzenbaum Scissors, Curved. 9"
2 each	Nelson-metzenbaum Scissors, Curved. 11"
2 each	Green Retractor. 8 1/2"
1 each	Parker Retractors, Set of 2. 5" X 3/4"
2 each	Langenbeck Retractor. 3/8" X 1 1/2"
2 each	Deaver Retractor, 1"blade. 9 1/2"
1 each	Adson Forceps, Mouse Teeth. 4 3/4"
1 each	Adson Forceps, Serrated. 4 3/4"
1 each	Satinsky Vena Cava Clamp. 10 1/4"
1 each	Potts-smith Scissors, 25 Degree Angle. 7 1/2"

TONSIL INSTRUMENT SET

1 each	Yankauer Suction Tube
1 each	Mayo-hegar Needle Holder, T.c. Jaw. 6"
1 each	Laforce Adenotome, Large
1 each	Laforce Adenotome, Medium
1 each	Laforce Adenotome, Small
3 each	Blades For Laforce Adenotome, Large
3 each	Blades For Laforce Adenotome, Medium

TRACHEOTOMY SET

1 each	Operating Scissors, Straight, S/b. 5 1/2"
1 each	Mayo Dissecting Scissors, Curved. 5 1/2"
6 each	Halstead Mosquito Forceps, Curved. 5"
2 each	Rochester-pean Forceps, Curved. 5 1/2"
10 each	Crile Forceps, Straight. 5 1/2"
2 each	Dressing Forceps, Serrated. 5"
2 each	Tissue Forceps, 2 X 3 Teeth. 5 1/2"
1 each	Crile-wood T.c. Jaw Needle Holder. 6"
2 each	Allis Forceps, 4 X 5 Teeth. 6"
1 each	Metzenbaum Scissors, Curved. 7"
1 each	Tracheal Retr. Sharp, 1 Prong
1 each	Laborde Dilator, Tri-valve
4 each	Backhaus Towel Clamp. 5 1/2"
4 each	Sponge Forceps, Straight, Serrated. 9 1/2"

OBSTETRICAL FORCEPS

2 each	Delee Forceps, Standard. 14"
1 each	Simpson Forceps, Short. 12"
1 each	Simpson Forceps, Long. 14"
2 each	Kieland Luikart Forceps. 15 3/4"
2 each	Simpson Luikart Forceps. 14"
1 each	Piper Forceps. 17 1/2"

OBSTETRICAL INSTRUMENTS

2 each	Operating Scissors, Straight, Blunt/blunt. 15 1/2"
6 each	Mayo Scissors, Straight. 6 3/4"
6 each	Mayo Scissors, Curved. 6 3/4"
12 each	Halstead Mosquito Forceps, Straight
1 each	Halstead Mosquito Forceps, Curved
12 each	Kelly Forceps, Straight
12 each	Kelly Forceps, Curved
12 each	Rochester-pean Forceps, Straight. 6 1/4"



Surgical Instruments Sets

QTY.	DESCRIPTION	QTY.	DESCRIPTION
12 each	Dressing Forceps. 6"	MISCELLANEOUS INSTRUMENTS	
12 each	Tissue Forceps, 2 X 3 Teeth. 6"	For Supplementing Sets or for Miscellaneous use Throughout the Hospital	
4 each	Mayo-hegar Needle Holders, T.c. Jaw. 7"	1 each	Finger Ring Cutter
4 each	Mayo-hegar Needle Holders, T.c. Jaw. 7"	1 each	Finger Nail Drill W/three Drills
2 each	Probe W/eye, Silver. 6"	6 each	Yankauer Suction Tubes
12 each	Allis Tissue Forceps, 5 X 6 Teeth. 6"	MISCELLANEOUS INSTRUMENTS	
24 each	Backhaus Towel Clamp. 5 1/4"	6 each	Sterilizer Utility Forceps, Straight 11"
24 each	Sponge Forceps, Straight, Serrated. 9 1/2"	12 each	Lister Bandage Scissors 5 1/2"
2 each	Mayo Scissors, Curve. 9"	6 each	Lister Bandage Scissors. 7 1/4"
2 each	Jackson Vaginal Retractors, Medium	6 each	Operating Scissors, Str., Sharp/sharp. 5 1/2"
2 each	Jackson Vaginal Retractors, Small	12 each	Operating Scissors, Str., Sharp/blunt. 5 1/2"
1 each	Auvard Weighted Speculum, Medium	6 each	Operating Scissors, Str., Blunt/blunt. 5 1/2"
2 each	Gelpi Retractor. 7"	6 each	Operating Scissors, Str., Sharp/blunt. 6 1/2"
2 each	Graves Vaginal Speculum, Small	6 each	Mayo Dissecting Scissors, Straight. 5 1/2"
OBSTETRICAL INSTRUMENTS		12 each	Mayo Dissecting Scissors, Straight. 6 3/4"
3 each	Graves Vaginal Speculum, Medium	12 each	Mayo Dissecting Scissors, Curved. 5 1/2"
2 each	Graves Vaginal Speculum, Large	12 each	Mayo Dissecting Scissors, Curved. 6 3/4"
1 each	Pederson Vaginal Speculum, Narrow	24 each	Halstead Mosquito Forceps, Straight. 5"
4 each	Bozeman Uterine Dressing Forceps, Str. 10"	48 each	Halstead Mosquito Forceps, Curved. 5"
2 each	Schroeder Uterine Tenaculum	24 each	Rochester-pean Forceps, Curved. 7 1/4"
2 each	Iowa Membrane Puncturing Forceps	24 each	Crile Forceps, Straight. 5 1/2"
1 each	Collyer Pelvimeter	12 each	Crile Forceps, Curved. 5 1/2"
DILATON AND CURETTAGE SET		1 each	Dressing Forceps, Serrated. 5 1/2"
1 each	Operating Scissors, Str., Sharp/blunt. 5 1/2"	12 each	Virtus Splinter Forceps. 6"
6 each	Backhaus Towel Clamp. 5 1/2"	12 each	Baumgartner Needle Holders, T.c. Jaw. 5 1/2"
4 each	Sponge Forceps., Serrated. 9 1/2"	6 each	Mayo-hegar Needle Holder, T.c. Jaw. 6"
1 each	Jackson Retractor, Large Blade. 4" X 1 1/2"	6 each	Mayo-hegar Needle Holder, T.c. Jaw. 7"
1 each	Jackson Retractor, Medium Blade. 3 1/2" X 1 1/2"	12 each	Mayo-hegar Needle Holder, T.c. Jaw. 8"
1 each	Jackson Retractor Small, Blade. 3" X 1 1/2"	12 each	Metzenbaum Scissors, Curved. 5 3/4"
1 each	Sims Retractor, Double Ended, No. 4	2 each	Adson Forceps, Serrated. 4 3/4"
1 each	Sims Retractor, Double Ended, No. 6	4 each	Malleable Retractor. 1" X 13"
1 each	Auvard Weighted Speculum, Medium	2 each	Malleable Retractor. 1 3/4" X 13"
1 each	Graves Speculum, Medium	2 each	Malleable Retractor. 2" X 13"
1 each	Goodell Uterine Dilator, Large. 13"	2 each	Volkman Retractor, Blunt, 2 Prong
1 each	Goodell Uterine Dilator, Small. 11"	2 each	Volkman Retractor, Blunt, 4 Prong
1 each	Hegar Dilators, Set Of 8	2 each	Volkman Retractor, Blunt, 6 Prong
1 each	Sims Sound, Plain	6 each	Strabismus Scissors, Straight. 4 1/4"
1 each	Simpson Sound, Graduated	6 each	Strabismus Scissors, Curved. 4 1/4"
1 each	Sims Curette, Size 1, Sharp	6 each	Iris Scissors, Curved. 4 1/2"
1 each	Sims Curette, Size 2	6 each	Iris Scissors, Straight. 4 1/2"
1 each	Sims Curette, Size 3	1 each	Coakley Curette, Set Of 6
1 each	Sims Curette, Size 4	1 set	Van Buren Sounds, Set Of 12, sizes 8-30 Fr.
1 each	Sims Curette, Size 5	1 each	Jones Towel Clamp, Cross Action. 3"
1 each	Sims Curette, Size 6	12 each	Sponge Forceps, Straight, Serrated. 9 1/2"
1 each	Thomas Uterine Curette, Size 1, Blunt	12 each	Sponge Forceps, Straight, Serrated. 7"
1 each	Thomas Uterine Curette, Size 2	6 each	Graves Speculum, Small
1 each	Thomas Uterine Curette, Size 3	1 each	Graves Speculum, Medium
1 each	Thomas Uterine Curette, Size 4	6 each	Graves Speculum, Large
1 each	Thomas Uterine Curette, Size 5	6 each	Suture Wire Cutting Scissors, Angular
1 each	Thomas Uterine Curette, Size 6	2 each	Ochsner Gall Bladder Trocar, 14 Fr.
1 each	Thoms-gaylor Biopsy Forceps, Improved Jaw. 8 1/2"	2 each	Ochsner Gall Bladder Trocar, 20 Fr.
1 each	Bozeman Uterine Dressing Forceps, Cvd. 10"	2 each	Ochsner Gall Bladder Trocar, 24 Fr.
2 each	Schroeder Tenaculum. 9"	24 each	Backhaus Towel Clamp. 3 1/2"
		36 each	Backhaus Towel Clamp. 5 1/2"
		12 each	Adson Forceps, Mouse Teeth. 4 3/4"
		2 each	Pederson Speculum, Small. 3 1/2" X 7 1/2"
		6 each	Bozeman Uterine Dressing Forceps, Cvd. 10"



Surgical Instruments Sets

QTY.	DESCRIPTION	QTY.	DESCRIPTION
PELVIC SET (Used With Major Basic Set)			
1 each	Mayo Dissecting Scissors, Straight. 6 3/4"	1 each	Liston Bone Cutting Forceps, Straight. 7 3/4"
1 each	Mayo Dissecting Scissors, Curved. 6 3/4"	1 each	Stille-liston Bone Forceps, Straight. 10 1/2"
12 each	Kelly Forceps, Curved. 5 1/2"	1 each	Stille-liston Bone Forceps, Curved. 10 1/2"
10 each	Ochsner Forceps, Curved. 8"	1 each	Stille-luer Bone Rongeur, Straight. 8 1/2"
6 each	Ochsner Forceps, Straight. 6 1/4"	1 each	Stille-luer Bone Rongeur, Angular. 9"
6 each	Ochsner Forceps, Straight. 8"	2 each	Gigli Bone Saws, Stille Pattern. 12"
1 each	Dressing Forceps. 6"	2 each	Gigli Bone Saws, Stille Pattern. 20"
1 each	Tissue Forceps, 1 X 2 Teeth. 6"	2 each	Gigli Bone Saws, Stille Pattern. 30"
2 each	Masson Needle Holder, T.c. Jaw. 10 1/2"	1 pair	Handles For Gigli Saw, Solid
6 each	Pean Forceps, Curved. 9"	1 each	Sayre Elevator, Sharp Double Ended 6 1/2"
6 each	Pean Forceps, Curved. 10"	1 each	Murphy Bone Skid, double Ended. 12"
4 each	Backhaus Towel Clamp. 5 1/2"	1 each	Bone File. 9 1/2"
2 each	Sponge Forceps, Straight, Serrated. 9 1/2"	2 each	Hibbs Retractor, Blade 3" X 1" Length. 9 1/4"
6 each	Heaney-ballantine Hysterectomy Forceps, Cvd. . . 8"	2 each	Green Retractor. 8 1/2"
2 each	Heaney Needle Holders, T.c. Jaws. 8 1/4"	1 each	Mayo-collins Retractor, Set of 2
1 each	Schroeder Tenaculum. 9"	1 each	Parker Retractor, Set Of 2
1 each	Schroeder Vulsellum Forceps, Straight. 9"	2 each	Langenbeck Retractor, Medium
PLASTIC SET			
1 each	Mayo Dissecting Scissors, Straight. 5 1/2"	2 each	Volkman Retractor, Blunt, 4 Prong. 8 1/2"
1 each	Mayo Dissecting Scissors, Curved. 5 1/2"	1 each	Lucae Bayonet Forceps. 5 1/2"
6 each	Halstead Mosquito Forceps, Straight. 5"	1 each	Metzenbaum Scissors, Curved. 7"
12 each	Halstead Mosquito Forceps, Curved. 5"	10 each	Backhaus Tower Clamp. 5 1/2"
4 each	Rochester-ochsner Forceps, Straight. 5 1/2"	8 each	Sponge Forceps, Str., Serrated. 9 1/2"
2 each	Crile-wood Needle Holder. 6"	1 each	Volkman Bone Curettes, Oval, # 000. 6 3/4"
2 each	Baumgartner Needle Holders. 5"	1 each	Volkman Bone Curettes, Oval, # 00. 6 3/4"
4 each	Allis Forceps, 4 X 5 Teeth. 6"	1 each	Volkman Bone Curettes, Oval, # 0. 6 3/4"
4 each	Senn Retractor, Sharp	1 each	Volkman Bone Curettes, Oval, # 1. 6 3/4"
1 each	Stevens Teotomoy Scissors, Cvd., Blunt	1 each	Volkman Bone Curettes, Oval, # 2. 6 3/4"
1 each	Stevens Teotomoy Scissors, Str., Blunt	1 each	Volkman Bone Curettes, Oval, # 3. 6 3/4"
1 each	Vienna Nasal Speculum, Adult	1 each	Volkman Bone Curettes, Oval, # 4. 6 3/4"
6 each	Backhaus Towel Clamp. 3 1/2"	1 each	Volkman Bone Curettes, Oval, # 5. 6 3/4"
BASIC ORTHOPEDIC SET			
2 each	Mayo Dissecting Scissors, Straight. 5 1/2"	1 each	Volkman Bone Curettes, Oval, # 6. 6 3/4"
1 each	Mayo Dissecting Scissors, Straight. 6 3/4"	POST MORTEM INSTRUMENTS	
2 each	Mayo Dissecting Scissors, Curved. 5 1/2"	1 each	Sterilizer Utility Forceps, Straight. 11"
1 each	Mayo Dissecting Scissors, Curved. 6 3/4"	1 each	Operating Scissors, Straight, S/s. 5"
12 each	Rochester-pean Forceps, Curved. 6 1/4"	1 each	Mayo Dissecting Scissors, Straight. 6 3/4"
10 each	Rochester-ochsner Forceps, Straight. 6 1/4"	1 each	Mayo Dissecting Scissors, Curved. 6 3/4"
18 each	Crile Forceps, Curved. 5 1/2"	1 each	Dressing Forceps, Serrated. 6"
2 each	Dressing Forceps, Serrated. 5 1/2"	1 each	Tissue Forceps, 1 X 2 Teeth. 6"
2 each	Tissue Forceps, 1 X 2 Teeth. 5 1/2"	1 each	Tissue Forceps, 1 X 2 Teeth. 10"
1 each	Crile-wood Needle Holder, T.c. Jaw. 6"	1 each	Probe, W/eye. 10"
2 each	Mayo-heagar Needle Holder, T.c. Jaw. 7"	1 each	Langeneack Metacarpal Saw. 4 1/2"
1 each	Suture Wire Cutting Scissors, Angular. 4 3/4"	1 each	Gluck Rib Shears. 8 3/4"
1 each	Grooved Director. 6"	1 each	Mayo Dissecting Scissors, Straight. 9"
2 each	Probes, W/eye, Silver. 6"	1 each	Cartilage Knife, Blade. 6 1/2"
10 each	Allis Forceps, 5 X 6 Teeth. 6"	1 each	Virchow Knife, Blade. 3 1/2"
1 each	Lane Screw Driver	1 each	Suture Wire Scissors, Angular. 4 3/4"
1 each	Van Buren Sequester Forceps, Bayonet . . . 8 1/2"	1 each	Director And Tongue Tie. 6"
1 each	Lane Bone Holding Forceps, W/ratchet. 13"	8 each	Allis Forceps, 4 X 5 Teeth. 6"
1 each	Kern Bone Holding Forceps, W/ratchet. 8 1/2"	2 each	Green Retractor. 8 1/2"
1 each	Lambotte Bone Forceps. 8 3/4"	1 each	Parker Retractor, D.e., Set Of 2. 7 1/8" X 7"
1 each	Lowman Clamp, Large	2 each	Richardson Retractor, Small
1 each	Lowman Clamp, Large	2 each	Volkman Retractors, Sharp, 4 Prong
		2 each	Senn D.e., Retractor. 6 1/4"
		1 each	Alm Retractor. 2 3/4"
		1 each	Metzenbaum Scissors, Curved. 5 3/4"
		12 each	Backhaus Towel Clamp. 5 1/4"
		8 each	Sponge Forceps, Straight Serrated 9 1/2"



Surgical Instruments Sets

QTY .	DESCRIPTION	QTY .	DESCRIPTION	
NEPHRECTOMY SET (Used With Major Basic Set)		NORMAL DELIVERY SET		
8 each	Rochester-pean Forceps, Curved. 8"	1 each	Large Bowel	
8 each	Rochester-ochsner Forceps, Straight. 8"	2 each	Sim,s Spacullm (small, Large)	
2 each	Dressing Forceps, Serrated. 10"	1 each	Episotomy Scissor	
2 each	Masson Needle Holder, T.c. Jaw. 10 1/2"	1 each	Non Toothed Forceps	
2 each	Lower Gall Duct Forceps. 7"	2 each	Artery Forceps Straight (small, Large)	
1 set	Grey Cystic Duct Forceps, W/o Teeth, Set Of 2 .	2 each	Artery Forceps Curved	
1 each	Blake Gall Stone Forceps, Curved. 8 1/4"	3 each	Spong Holding Forceps	
1 each	Desjardin Gall Stone Forceps. 9 1/4"	1 each	Needle Holder Mediumsize	
1 each	Randall Kidney Stone Forceps, No.1 .	1 each	Kidney Tray Big Size	
1 each	Randall Kidney Stone Forceps, No.2 .	2 each	Kidney Tray Small Size	
1 each	Randall Kidney Stone Forceps, No.3 .			
1 each	Randall Kidney Stone Forceps, No.4 .			
2 each	Herrick Kidney Clamp. 9"			
1 each	Ochsner Trocar, Size 20 Fr. .			
2 each	Deaver Retractor. 2" X 12"			
1 each	Metzenbaum Scissors, Curved. 7"			
1 each	Bozeman Dressing Forceps, Curved. 10"			
VENOUS SECTION SET		D&C SET		
1 each	Operating Scissors, Str., Sharp/blunt. 5 1/2"	1 each	Instrument Tray	
2 each	Halstead Mosquito Forceps. Straight. 5"	2 each	Bowel Small and Medium	
2 each	Halstead Mosquito Forceps. Curved. 5"	2 each	Sim,s Spaculum	
1 each	Dressing Forceps, Serrated. 5 1/2"	3 each	Spong Holder	
1 each	Crile-wood Needle Holder, T.c. Jaw. 6"	12 each	Dilator Set In All Size	
1 each	Adson Tissue Forceps, W/mouse Teeth .	1 each	Matalic Catheter	
4 each	Backhaus Towel Clamp. 5 1/4"	1 each	B.p Handle No 3	
		1 each	Tissue Scissor	
		3 each	Suction Curette	
		3 each	Curett With Handle Small, Medium, Large	
		2 each	Toothed And Non Toothed Forceps	
		2 each	Complete Set Of Curette	
		1 each	Wall, Spaculum Forceps	
RECTAL INSTRUMENT SET		MICRO SURGERY SET (MAJOR SET)		
1 each	Operating Scissors, Str., Sharp/blunt. 5 1/2"	1 each	Micro Tying Forceps	4 3/8" 0.3mm tip
1 each	Mayo Dissecting Scissors, Straight. 6 3/4"	1 each	Barraquer Needle Holder	5" curved, with lock
1 each	Mayo Dissecting Scissors, Curved. 6 3/4"	1 each	Micro Scissors	5" straight 1x2 teeth
8 each	Rochester-pean Forceps, Curved. 6 1/2"	1 each	Graefe Iris Forceps	4" straight 1x2 teeth
8 each	Rochester-ochsner Forceps, Curved. 6 1/2"	1 each	Micro Clamps Applicator	5 1/2"
6 each	Pennington Clamp, Straight. 5 3/4"	1 each	Approximation Clamp	without suture holding fra
8 each	Crile Forceps, Straight. 5 1/2"	1 each	Micro Scissors (reill)	5 1/2" curved, serrated
2 each	Dressing Forceps, Serrated. 6"	1 each	Platform Tying Forceps	5 1/2" 0.7mm tip
2 each	Tissue Forceps, 1 X 2 Teeth. 6"	1 each	Hepp-scheidel Dissecting Scissors	5 1/2" curved
2 each	Mayo-hegar Needle Holder, T.c. Jaw. 7"	1 each	Clamp Box	for micro clamps
1 each	Suture Wire Scissors, Angular. 4 1/4"	1 each	Sterilization And Instrument Case	for eight piece micro set
1 each	Grooved Director. 6"	1 each	Instruments Rack	for eight piece micro set
2 each	Probe, W/eye, Silver. 6"			
4 each	Allis Forceps, 5 X 6 Teeth. 7 3/4"			
1 each	Buie Pile Clamp. 8 1/4"			
6 each	Pratt "t" Shaped Forceps. 5 1/2"			
1 each	Smith Anal Retractor .			
1 each	Sim Speculum, Fenestrated Blad.e .			
1 each	Hirschman Anascope, Small .			
1 each	Hirschman Anascope, Medium .			
1 each	Hirschman Anascope, Large .			
6 each	Backhaus Towel Clamp. 5 1/2"			
6 each	Foerster Sponge Forceps, Straight, Serr.. . . . 9 1/2"			
		MICRO SURGERY SET (MINI SET)		
		1 each	Platform Tying Forceps	5 1/2" 0.7mm tip
		1 each	Barraquer Needle Holder	5" curved with lock
		1 each	Micro Scissors	5" curved, long blade
		1 each	Sterilization and Instrument Case	for three piece micro set



Item Name	Page #
A	
Adson / Micro-Adson	18
Adson Titanium	113
Adson Forceps	12,18,69
Adson-Brown	18
Allis / Adair	96
Allis	96
Alm	50
Anaesthophon	59
Anderson-Adson	51
Andrews-Pynchon Suction Tube	93
Asch	89
Atrauma Debakey	06
Auvard	104
B	
Babcock	96,97
Baby Crile-Wood	39
Baby Senn-Miller / Double Ended	45
Baby-Adson Forceps	12
Baby-Mixer Forceps	13
Baby-Satinsky	92
Backhaus Towel Clamps	14
Balfour	53,54
Balfour/Baby	52
Balfour/US Type	52
Ball & Socket	15
Ballenger	87
Bane	69
Barraquer	43,50
Barraquer/with Catch	43
Barraya Titanium	113
Baumgartner	39
Beckmann	51
Beckmann/Adson	52
Beckmann/Eaton	51
Bengolea Titanium	114
Bennett	66
Bergmann	33,66
Beyer	70
Bishop-Harman	22
Blount Knee Retractor	73
Böhler	72
Boies	89
Bone Reduction Forceps	63
Bone Hook	73
Bonney Dissecting	19
Bowman-Lacrima / Double Ended-Malleable	57
Boy-Allis	96
Bozemann	23
Braun	108
Bridge Forceps	13
Brown	19,39
Brunner	98
Bruns	68
Buck Hammar	59
Burford / Aluminum	91
Byrd	94
C	
Caliper	74
Caplan	31
Carmalt Splinter	21
Carroll	41,73,74
Castroviejo / 1:2	21
Castroviejo	43,75
Castroviejo-Epker	75
Catheter	103
Chandler	65
Christophe Forceps	08

Item Name	Page #
Coller Forceps	09
Collin Forceps	34,36,54,66
Colver	82
Cone	52
Converse	42
Cooley / Pediatric	91
Cooley	94
Corwin Round	34
Cottle	87,66
Cottle Bone	89
Cottle Cartilage Crusher	89
Coupland	94
Crafoord Titanium	116
Crile Forceps	09
Crile Titanium	114
Crile-Rankin Forceps	09
Crile-Wood	40
Cushing Titanium	113
Cushing	19,45,78
Cuticle Nipper	57
D	
Dressing tip standard & USA Model	17
Davis-Boyle	80
De Bakey / Lightweight	21
De Bakey / Pediatric	92
De Bakey Titanium	116
De Martel-Wolfson	101
De Vilbiss	104
Dean-Shellcross	14
Deaver	49
DeBakey	21
Debakey-Vascular/Extra Light	41
Dejerine	59
Denis Browne Titanium	120
Derf	42
Diagnostic Sets	60,61
Diamant Pin Cutter	63
Dieffenbach	35,36
Diethrich/Micro Coronary	06
Doyen Titanium	117
Doyen	49,56,97
Duplex	59
Duplex-Baby	59
E	
Ear syringes / Reiner	85
Echlin / Stille-Luer	71
ENT-Dean	30
Esmarch	33
Eufate Pasque	41
Eves	82
Ewald / Hudson	19
Eye Scissors	29
F	
Fehland	99
Fergusson	25
Ferris-Kerrison	77
Ferris-Smith	19
Finger Rake	45
Finnocchietto / Baby	91
Finnocchietto	91
Finocchietto Titanium	120
Flat Nosed Parallel Pliers	65
Foerster-Ballenger	23
Fomon	88
Foregger	62
Fox	58
Frazier	93
Frazier/Dura Hook	44

Item Name	Page #
Freer "D" Shaped	87
Freer	44,87
Friedman & Mini-Friedman	69
G	
Gelpi	52
Gemini Forceps	13
Gerald	20
Gerster-Lowmann	66
Gerzog	66
Gillies Titanium	113
Gillies	20
Gillies/Skin Hook	44
Gil-Vernet	45
Goelet	50
Gorney Septal Scissors	31
Graduated Metal Ruler	74
Graduated Tape Measurer	74
Graves	104
Green	47
Guedel-Negus / Miller	61
Guilmet Titanium	113
Guthrie/Skin Hooks O'Conner	44
H	
Hajek	67
Halsey	39
Halstead Titanium	114
Halsted-Mosquito	07
Handle T-form	56
Hank Gyne	106
Hank	106
Harrington/Sweetheart	48
Hartmann / Mosquito	07
Hartmann	82,83,85
Hartmann/French	84
Harvey	30
Heaney / Ballantine Hysterectomy Forceps	14
Heaney / Rezak	14
Heaney	14
Heaney-Gyne	39
Hegar Gyne	107
Hegar	106
Hegemann-Diethrich	28
HEISS Titanium	115
Heiss	50
Hennig	32
Heywood-Smith	23
Hibbs-Spratt / Angled	68
Hibbs-Spratt / Straight	68
Hohmann	47
Hulka Gyne	105
Hunter	20
I	
Inge	77
Intestinal	20
Intracardiac Titanium	116
Iowa	111
Iris Scissors	29
Israel	46
Ives-Fansler	97
J	
Jackson	81
Jacobs	111
Jacobson Forceps	10,43



Item Name	Page #
Jameson.....	74
Jansen.....	50
Jansen-Gruenwald.....	21
Jansen-Middleton.....	87
Javid.....	92
Jeweler's Forceps.....	22
John-Hopkins.....	06
Jones / Perforating.....	14
Jorgenson.....	25
Joseph.....	87
Joseph/Skin Hooks.....	44
K	
Kantrowitz Forceps.....	13
Kelly Forceps.....	08
Kelly-Rankin Forceps.....	08
Kevorkian.....	110,108
Key.....	66
Killian.....	86
Kleinert-Kutz.....	70
Kleinert-Kutz/Skin.....	44
Knapp.....	45
Knight.....	28
Knives.....	35,36,37
Knife Handle / Long Angled #3.....	34
Knife Handle #3.....	34
Knife Handle/Graduated #3.....	34
Knife Handle/Long #3.....	34
Knife Handle #4.....	34
Knife Handle #7.....	34
Kocher / Delicate Forceps.....	11
Kocher.....	46,98,99
Koerte.....	46
Kogan Gyne.....	105
Knife Blades.....	34
L	
La Borde.....	81
Lahey Forceps.....	13,47,81
Lambotte.....	65
Larry-Rectal / Grooved Director.....	57
Laryngeal Mirror with Handle.....	83
Laryngoscope.....	61-62
Lead Hand.....	73
Leksell-Stille.....	71
Lempert.....	69
Leriche Titanium.....	114
Lewin / Bolt / Pin Cutter.....	64
Lewin.....	63
Lewis.....	88
Lightweight-Adson.....	18
Lincoln.....	28
Lister Bandage.....	33
Liston.....	71,72
Liston-Key (Horsley).....	72
Littauer / Liston.....	71
Littler.....	29
Lockhart-Mummery.....	56
Locking Pliers.....	65
Lorna.....	15
Love-Gruenwald.....	78
Lucae / ENT.....	18
Luer.....	69,93
M	
Magill.....	16

Item Name	Page #
Malleable Ribbone Retrater.....	54
Maltz.....	88
Markham/Meyerding.....	52
Mathieu Rectal Specula.....	101
Mayfield.....	71
Mayo / Develed Blades.....	26
Mayo Dissecting Scissors.....	25,26
Mayo-hegar Titanium.....	116
Mayo-Hegar.....	40
Mayo-Noble / Beveled Blades.....	25
Mayo-Russian.....	20
Mayo-Stille / Rounded Blades.....	26
Mc Indoe Titanium.....	113
McIntosh.....	62
Mcivor.....	80
Meeker Forceps.....	12
Metzenbaum / Baby / Delicate.....	26,27
Metzenbaum.....	27
Metzenbaum-Jamison.....	27
Meyerding.....	49
Micro ENT/Bellucci.....	84
Micro Scissor.....	31
Micro Scissors Titanium.....	118,119
Micro Suture Tying Forceps Titanium.....	117
Micro-Ear Scissors.....	83
Micro-Ear.....	84
Micro-Mosquito / Petit Point.....	07
Microvascular Titanium.....	116
Miller.....	23,69
Milligan.....	50
Mini Hohmann.....	46
Mini Lambotte.....	65
Mini-Gemini Forceps.....	13
Miskimon.....	50
Mixer Forceps.....	13
Model USA / Dressing.....	17,63
Modified Bonn.....	22
Modified Wien.....	86
Modified-USA.....	28
Molt.....	80
Möltgen.....	74
N	
Nail Nipper.....	57,58
Nabatoff.....	56
Nail Clipper.....	57
Narrow Tip Dressing Forceps.....	17
Nasal Specula.....	86
Needle Nose Pliers.....	64
Needle Nosed.....	73
Nelaton Prob.....	56
Nelson.....	20
New Orleans.....	41
News.....	81
Non-Perforating.....	15
O	
O'Conner.....	44
O' SHAUGHNESSY Titanium.....	115
O'Sullivan/O'Connor.....	53
Oehler.....	19
Olsen-Hegar/With Scissor.....	41
Ombredanne.....	66
Operating Fine Scissors.....	28
Ophthalmic Iris.....	29
Ortho / With Plastic Caps.....	67
Ortho / Mallets.....	63,67
P	
Penfield.....	77
Parker-Langenbeck/Army-Navy.....	49

Item Name	Page #
Partsch.....	67
Payr.....	100
Payr-Baby.....	100
Pean - Delicate Forceps.....	09
Pederson Speculam.....	104
Peers-Bertman.....	15
Pennington / Siezing.....	06
Piper.....	111
Plaster-Bruns.....	63
Poole Suction Tube.....	93
Potts Titanium.....	118
Potts-Smith.....	27
Pozzi Gyne.....	105
Pratt Gyne.....	106
Pratt.....	14,97
Presbyterian.....	16
Probe With Eye / Double Ended.....	57
Probe with Eye.....	57
Providence / Hospital Forceps.....	10
Putti.....	69
R	
Ragnell-Davis / Double Ended.....	45
Ralk.....	20
Randall.....	103
Recamier / Sims Curette.....	108
Rectal Specula.....	97,101
Reiner / Ear Syringe.....	85
Resano.....	98
Reverdin Needle.....	43
Reynolds Tenotomy.....	27
Ribban Retrator Malleable.....	54
Richardson/Kelly.....	48
Richardson-Eastman.....	47
Ring Cutter.....	58
Rochester / Carmalt Forceps.....	12
Rochester / Ochsner / Kocher Forceps.....	11
Rochester / Pean Forceps.....	10
Roeder.....	15
Rollet.....	45
Rubin.....	89
Ruskin.....	70
Ruskin-Liston.....	72
Ruskin-Mini.....	70
Russian-Modell Grasping.....	20
Ryder French-Eye/Standard/Micro.....	42
Ryder Titanium.....	116
S	
Saw Blade for Ring Cutter.....	58
Sarot Forceps.....	12,42
Satinsky.....	92
Self Retaining Retractors.....	51,52,53
Schlesinger.....	79
Schmidt / Sawtell Forceps.....	12
Schröder Gyne.....	105
Schröder.....	105
Semb.....	70
Senn-Miller / Double Ended.....	45
Senn-Miller / Baby.....	45
Shield for Ear Syringe.....	85
Sims Curettes / recamier.....	108
Sims Uterine.....	108
Simpson.....	111
Sims.....	101,108
Singley-Tuttle.....	20
Slip Joint.....	72
Speaking Tube.....	95
Spencer / Delicate.....	30
Spencer.....	30
Spencer-Wells Forceps.....	11



Item Name	Page #
Splinter	21
Spreader	91
Spratt / Hollow Handle	67
Spurling	79
S-Shaped Retractors.....	49
Standard / Operating Scissors	25
Standard / Thumb	17
Stille	32,70
Stille-Listo.....	72
Stille-Luer	71
Stille-Ruskin	71
Strabismus	29
St. Marks	47
T	
Takahashi	86
Taylor Hammer	59
Tendon	73
Tenotomy.....	30
Tenotomy-Stevens	30
Thomas-Allis	96
Tischler	109,110
Tischler-Kevorkian.....	109
Tischler-Morgan.....	110
Tohoku.....	15
Townley	74
Traction Cable	56
Trachea Tube	95
Trousseau	81
Tubing Clamp	16
Tuning Forks Aluminum.....	85
Tyding.....	82
U	
Universal	33,69
Uterine dilators	106,107
Utility Shears	31
V	
Van Buren	103
Vanderbilt / Vessel Clamp	08
Vaginal Speculams.....	104
Vernier	74
Verres.....	58
Volkman	46,67,68
W	
Waller Cutter	58
Walton-Liston	72
Wangensteen	41
Wartenberg.....	59
Webster	39
Weider	80
Weitlaner	51
Westcott	33
White	81
Williger	68
Wire Cutter	63,64
Wire Twisting.....	64
Wires/Loops	82
Wire-Universal Scissors	33
Wolff	32
Y	
Yankauer	93
Yasargil Titanium / Micro Sissors.....	118
Z	
Zaufal-Jansen	70

Item Name	Page #
-----------	--------

Item Name	Page #
-----------	--------



Part No.	Page #	Part No.	Page #	Part No.	Page #	Part No.	Page #
KTS0023401	74	KTS0422002	62	KTS0800411	28	KTS0834111	29
KTS0023403	74	KTS0422003	62	KTS0800511	28	KTS0834211	29
KTS0023405	74	KTS0422004	62	KTS0801914	25	KTS0834312	29
KTS0023407	74	KTS0422300	62	KTS0801918	25	KTS0834411	29
KTS0023409	74	KTS0422500	61	KTS0802718	27	KTS0834511	29
KTS0023411	74	KTS0422501	61	KTS0802723	27	KTS0836111	30
KTS0023413	75	KTS0422502	61	KTS0803012	28	KTS0836210	30
KTS0023415	75	KTS0422503	61	KTS0803025	28	KTS0836211	30
KTS0023417	75	KTS0422504	61	KTS0803045	28	KTS0836310	30
KTS0023419	75	KTS0422600	61	KTS0803060	28	KTS0836311	30
KTS0023421	75	KTS0440310	84	KTS0803090	28	KTS0836314	30
KTS0023423	75	KTS0440315	84	KTS0806109	30	KTS0837111	29
KTS0025710	60	KTS0440320	84	KTS0809711	29	KTS0837511	29
KTS0025712	60	KTS0610116	34	KTS0809714	29	KTS0837602	29
KTS0025714	60	KTS0610300	34	KTS0810214	25	KTS0837701	29
KTS0025716	60	KTS0610301	34	KTS0810215	25	KTS0837702	29
KTS0025718	60	KTS0610302	34	KTS0810216	25	KTS0838410	29
KTS0025720	61	KTS0610303	34	KTS0811014	25	KTS0838510	29
KTS0025722	61	KTS0610400	34	KTS0811414	25	KTS0840918	28
KTS0088400	95	KTS0610700	34	KTS0812014	25	KTS0842717	30
KTS0088401	95	KTS0610901	34	KTS0816014	26	KTS0848025	27
KTS0088402	95	KTS0639801	35	KTS0816015	26	KTS0848040	27
KTS0088403	95	KTS0639802	35	KTS0816017	26	KTS0848060	27
KTS0088404	95	KTS0639803	35	KTS0816019	26	KTS0855010	28
KTS0088405	95	KTS0639804	35	KTS0816023	26	KTS0859411	29
KTS0088406	95	KTS0639805	35	KTS0816028	26	KTS0865022	25
KTS0088500	95	KTS0639806	35	KTS0816114	26	KTS0870722	31
KTS0088501	95	KTS0639807	35	KTS0816117	26	KTS0870622	31
KTS0088502	95	KTS0639808	35	KTS0816119	26	KTS0880111	30
KTS0088503	95	KTS0639809	35	KTS0816128	26	KTS0880209	30
KTS0088504	95	KTS0639815	35	KTS0816414	26	KTS0883112	33
KTS0088505	95	KTS0639816	35	KTS0816417	26	KTS0883512	33
KTS0088506	95	KTS0640001	36	KTS0816423	26	KTS0885421	31
KTS0088507	95	KTS0640002	36	KTS0816514	26	KTS0886212	30
KTS0228817	59	KTS0640003	36	KTS0816517	26	KTS0890109	33
KTS0228820	59	KTS0640004	36	KTS0816523	26	KTS0890111	33
KTS0230420	59	KTS0640005	36	KTS0817017	26	KTS0890114	33
KTS0230818	59	KTS0640006	36	KTS0817117	26	KTS0890115	33
KTS0231528	59	KTS0640007	36	KTS0817415	26	KTS0890118	33
KTS0231529	59	KTS0640008	36	KTS0817417	26	KTS0890120	33
KTS0231530	59	KTS0640009	36	KTS0817517	26	KTS0890121	33
KTS0234015	74	KTS0640010	36	KTS0818017	25	KTS0890518	33
KTS0235015	74	KTS0640401	37	KTS0818117	25	KTS0891123	33
KTS0235200	59	KTS0640403	37	KTS0828011	26	KTS0891320	33
KTS0413010	93	KTS0640405	37	KTS0828015	27	KTS0891322	33
KTS0414008	93	KTS0640407	37	KTS0828018	27	KTS0893723	63
KTS0415027	93	KTS0640409	37	KTS0828111	26	KTS0895014	33
KTS0415222	93	KTS0640410	37	KTS0828115	27	KTS0895015	33
KTS0416500	94	KTS0640412	37	KTS0828118	27	KTS0895016	33
KTS0416715	94	KTS0640414	37	KTS0828120	27	KTS0895018	33
KTS0416725	94	KTS0640416	37	KTS0828123	27	KTS0895245	31
KTS0416735	94	KTS0640418	37	KTS0828315	27	KTS0910211	33
KTS0417035	94	KTS0640601	36	KTS0828318	27	KTS0910311	33
KTS0417133	94	KTS0640602	36	KTS0828328	27	KTS0910411	33
KTS0417530	94	KTS0640603	36	KTS0828415	27	KTS0920418	31
KTS0417731	94	KTS0640604	36	KTS0828515	27	KTS1000112	18
KTS0419017	16	KTS0640605	36	KTS0828518	27	KTS1000212	18
KTS0419020	16	KTS0640606	36	KTS0828523	27	KTS1000215	18
KTS0419025	16	KTS0642028	34	KTS0828530	27	KTS1010213	17
KTS0420300	62	KTS0716014	25	KTS0828615	27	KTS1010214	17
KTS0420400	62	KTS0716017	25	KTS0828715	27	KTS1010216	17
KTS0421000	62	KTS0716117	25	KTS0828718	27	KTS1010225	17
KTS0421001	62	KTS0728011	26	KTS0828723	27	KTS1010414	17
KTS0421002	62	KTS0728115	27	KTS0828914	27	KTS1010415	17
KTS0421003	62	KTS0728118	27	KTS0829311	27	KTS1010418	17
KTS0421004	62	KTS0728123	27	KTS0832115	27	KTS1010425	17
KTS0422001	62	KTS0734110	29	KTS0832118	27	KTS1011014	17
		KTS0734111	29	KTS0834009	29	KTS1011016	17
		KTS0736311	30	KTS0834010	29	KTS1011020	17
		KTS0800117	28	KTS0834011	29	KTS1012014	17
		KTS0800217	28	KTS0834109	29	KTS1012016	17
		KTS0800315	27	KTS0834110	29	KTS1012025	17



Part No.	Page #	Part No.	Page #	Part No.	Page #	Part No.	Page #
KTS1012214	17	KTS1130816	21	KTS1232026	11	KTS1617025	23
KTS1012215	17	KTS1130820	21	KTS1232030	11	KTS1617118	23
KTS1012225	17	KTS1130824	21	KTS1232116	11	KTS1617120	23
KTS1013014	17	KTS1130830	21	KTS1232118	11	KTS1617125	23
KTS1013016	17	KTS1131020	21	KTS1232120	11	KTS1617225	23
KTS1013020	17	KTS1131024	21	KTS1232213	11	KTS1617325	23
KTS1014212	18	KTS1131030	21	KTS1232214	11	KTS1623126	23
KTS1014412	18	KTS1132218	19	KTS1232216	11	KTS1800416	52
KTS1015212	19	KTS1133220	21	KTS1232314	11	KTS1800513	52
KTS1015412	19	KTS1207001	06	KTS1233115	11	KTS1800516	52
KTS1016812	18	KTS1207101	06	KTS1234116	12	KTS1800520	52
KTS1016912	18	KTS1207108	06	KTS1234120	12	KTS1800631	51
KTS1017012	18	KTS1207112	06	KTS1235815	06	KTS1800731	52
KTS1017312	18	KTS1207208	06	KTS1235820	06	KTS1801031	52
KTS1017412	18	KTS1207301	06	KTS1236015	14	KTS1801045	52
KTS1017612	18	KTS1209307	06	KTS1241118	12	KTS1813121	49
KTS1018212	18	KTS1211000	06	KTS1241514	12	KTS1813519	50
KTS1018612	18	KTS1218113	08	KTS1241518	12	KTS1814000	47
KTS1018812	18	KTS1218115	08	KTS1243113	13	KTS1814001	47
KTS1018915	19	KTS1221010	07	KTS1243118	13	KTS1814002	47
KTS1019912	18	KTS1221109	07	KTS1243120	13	KTS1815001	48
KTS1022017	19	KTS1221110	07	KTS1243318	13	KTS1815002	48
KTS1022502	19	KTS1222012	07	KTS1243520	13	KTS1815003	48
KTS1023017	19	KTS1222014	07	KTS1243524	13	KTS1815004	48
KTS1024618	19	KTS1222112	07	KTS1243528	13	KTS1815231	48
KTS1024818	19	KTS1222114	07	KTS1244721	08	KTS1815251	48
KTS1025615	20	KTS1222210	07	KTS1245124	12	KTS1815626	48
KTS1025620	20	KTS1222212	07	KTS1245918	12	KTS1815851	48
KTS1025625	20	KTS1222310	07	KTS1246128	13	KTS1815876	48
KTS1026114	19	KTS1222312	07	KTS1247116	13	KTS1816017	45
KTS1026116	19	KTS1223012	07	KTS1247119	13	KTS1816117	45
KTS1026118	19	KTS1223014	07	KTS1247123	13	KTS1816215	45
KTS1026120	19	KTS1223112	07	KTS1247514	13	KTS1816315	45
KTS1026125	19	KTS1223114	07	KTS1247922	13	KTS1816414	45
KTS1026130	19	KTS1224014	08	KTS1249119	13	KTS1816812	54
KTS1027215	20	KTS1224114	08	KTS1259120	14	KTS1816817	54
KTS1030218	20	KTS1224216	08	KTS1259420	14	KTS1816820	54
KTS1031018	20	KTS1224316	08	KTS1259520	14	KTS1816825	54
KTS1032023	20	KTS1225014	09	KTS1259618	14	KTS1817025	54
KTS1032215	20	KTS1225114	09	KTS1259921	14	KTS1817032	54
KTS1045118	19	KTS1226014	09	KTS1260021	14	KTS1817040	54
KTS1045218	19	KTS1226114	09	KTS1270416	16	KTS1817050	54
KTS1045615	20	KTS1226216	09	KTS1270418	16	KTS1817612	49
KTS1045623	20	KTS1226316	09	KTS1271018	16	KTS1818325	49
KTS1060301	23	KTS1227014	10	KTS1271020	16	KTS1818750	49
KTS1060302	23	KTS1227114	10	KTS1271218	16	KTS1819140	48
KTS1080209	21	KTS1228014	09	KTS1299215	73	KTS1819162	48
KTS1080211	21	KTS1228016	09	KTS1311112	92	KTS1819164	48
KTS1081210	20	KTS1228114	09	KTS1325320	92	KTS1821201	44
KTS1081310	20	KTS1228116	09	KTS1325324	92	KTS1821202	44
KTS1081610	21	KTS1229518	10	KTS1325325	92	KTS1821301	44
KTS1081710	21	KTS1231016	10	KTS1325326	92	KTS1821302	44
KTS1082815	20	KTS1231018	10	KTS1325327	92	KTS1821303	44
KTS1082915	20	KTS1231020	10	KTS1334318	92	KTS1821413	44
KTS1083010	22	KTS1231030	10	KTS1342114	91	KTS1821601	44
KTS1083012	22	KTS1231116	10	KTS1392310	15	KTS1821602	44
KTS1083015	22	KTS1231118	10	KTS1392313	15	KTS1821605	44
KTS1083040	22	KTS1231120	10	KTS1410309	14	KTS1821607	44
KTS1083050	22	KTS1231122	10	KTS1411109	14	KTS1821610	44
KTS1083060	22	KTS1231124	10	KTS1411111	14	KTS1822001	45
KTS1083070	22	KTS1231126	10	KTS1411113	14	KTS1822816	44
KTS1083080	22	KTS1231130	10	KTS1411115	14	KTS1823421	45
KTS1083090	22	KTS1231214	09	KTS1411713	15	KTS1826513	45
KTS1084015	21	KTS1231216	09	KTS1412110	15	KTS1826515	45
KTS1130416	21	KTS1231314	09	KTS1412111	15	KTS1828001	46
KTS1130420	21	KTS1231316	09	KTS1412113	15	KTS1828002	46
KTS1130424	21	KTS1232014	11	KTS1412313	15	KTS1828003	46
KTS1130430	21	KTS1232016	11	KTS1412511	15		
KTS1130613	21	KTS1232018	11	KTS1412513	15		
KTS1130616	21	KTS1232020	11	KTS1416514	15		
KTS1130620	21	KTS1232022	11	KTS1617018	23		
KTS1130624	21	KTS1232024	11	KTS1617020	23		



Part No.	Page #	Part No.	Page #	Part No.	Page #	Part No.	Page #
KTS1828004	46	KTS2010616	57	KTS2434418	39	KTS3241113	63
KTS1828006	46	KTS2011114	56	KTS2441414	41	KTS3244120	66
KTS1828101	46	KTS2011716	56	KTS2441417	41	KTS3245001	73
KTS1828102	46	KTS2011815	57	KTS2441419	41	KTS3246996	67
KTS1828103	46	KTS2012101	56	KTS2442014	41	KTS3246997	67
KTS1828104	46	KTS2225010	58	KTS2442017	41	KTS3246998	67
KTS1828106	46	KTS2225012	58	KTS2442019	41	KTS3247000	67
KTS1830419	47	KTS2225015	58	KTS2466114	43	KTS3247001	67
KTS1833304	46	KTS2238200	93	KTS2466314	43	KTS3247002	67
KTS1833306	46	KTS2367911	73	KTS2467414	43	KTS3247003	67
KTS1834822	47	KTS2410013	39	KTS2467614	43	KTS3247004	67
KTS1836440	46	KTS2410413	39	KTS2467618	43	KTS3247005	67
KTS1836608	46	KTS2411412	42	KTS2467714	43	KTS3247006	67
KTS1837825	49	KTS2411612	42	KTS2469421	43	KTS3247197	68
KTS1852075	49	KTS2413013	39	KTS2469521	43	KTS3247198	68
KTS1852301G	47	KTS2413413	39	KTS2480818	64	KTS3247200	68
KTS1852302G	47	KTS2414013	39	KTS2482214	43	KTS3247201	68
KTS1852303G	47	KTS2414413	39	KTS2482414	43	KTS3247202	68
KTS1862007	50	KTS2415014	39	KTS2482514	43	KTS3247203	68
KTS1862107	50	KTS2415412	39	KTS2620119	43	KTS3247204	68
KTS1862410	50	KTS2415414	39	KTS2620123	43	KTS3247205	68
KTS1862910	50	KTS2416015	40	KTS2620219	43	KTS3247206	68
KTS1863410	50	KTS2416018	40	KTS2620223	43	KTS3247495	67
KTS1863510	50	KTS2416115	39	KTS2620319	43	KTS3247496	67
KTS1864826	51	KTS2416215	39	KTS2620323	43	KTS3247497	67
KTS1864926	51	KTS2416415	40	KTS2620419	43	KTS3247498	67
KTS1865010	51	KTS2416418	40	KTS2620423	43	KTS3247500	67
KTS1865110	51	KTS2416420	40	KTS3010320	32	KTS3247501	67
KTS1865413	51	KTS2416423	40	KTS3010323	32	KTS3247502	67
KTS1865513	51	KTS2416430	40	KTS3010326	32	KTS3247503	67
KTS1867013	51	KTS2416615	40	KTS3010337	32	KTS3247504	67
KTS1867016	51	KTS2418014	40	KTS3035028	32	KTS3247505	67
KTS1867020	51	KTS2418016	40	KTS3036623	63	KTS3249213	68
KTS1867025	51	KTS2418018	40	KTS3039518	32	KTS3249216	68
KTS1867113	51	KTS2418020	40	KTS3039724	32	KTS3249220	68
KTS1867116	51	KTS2418024	40	KTS3200218	64	KTS3249600	68
KTS1867120	51	KTS2418026	40	KTS3210112	69	KTS3261200	68
KTS1867125	51	KTS2418030	40	KTS3210114	69	KTS3261201	68
KTS1867820	51	KTS2418414	40	KTS3213515	69	KTS3261202	68
KTS1867920	51	KTS2418416	40	KTS3214220	69	KTS3261203	68
KTS1869031	51	KTS2418418	40	KTS3215116	69	KTS3261204	68
KTS1869132	51	KTS2418420	40	KTS3215418	69	KTS3261205	68
KTS1872802	50	KTS2418424	40	KTS3215518	69	KTS3261206	68
KTS1875400	53	KTS2418426	40	KTS3215613	69	KTS3261297	68
KTS1875401	53	KTS2418430	40	KTS3216717	71	KTS3261298	68
KTS1875402	53	KTS2419615	41	KTS3217115	70	KTS3261300	68
KTS1875403	53	KTS2419818	41	KTS3217523	70	KTS3261301	68
KTS1876418	53	KTS2421418	41	KTS3220118	70	KTS3261302	68
KTS1876425	54	KTS2421420	41	KTS3220518	70	KTS3261303	68
KTS1876545	53	KTS2421423	41	KTS3221115	72	KTS3261304	68
KTS1876570	54	KTS2421426	41	KTS3221315	72	KTS3261305	68
KTS1876690	52	KTS2421431	41	KTS3221715	70	KTS3261306	68
KTS1877022	50	KTS2422213	42	KTS3221718	70	KTS3261397	68
KTS1880209	52	KTS2422215	42	KTS3223022	71	KTS3261398	68
KTS1880211	52	KTS2422218	42	KTS3223122	71	KTS3265303	65
KTS1880214	52	KTS2422413	42	KTS3223523	71	KTS3265304	65
KTS1880218	52	KTS2422415	42	KTS3223723	71	KTS3265305	65
KTS1881210	54	KTS2422418	42	KTS3224523	71	KTS3265306	65
KTS1881260	54	KTS2422420	42	KTS3224723	71	KTS3265308	65
KTS1881280	54	KTS2422423	42	KTS3225902	71	KTS3265310	65
KTS1882018	52	KTS2422426	42	KTS3226122	70	KTS3265312	65
KTS1882118	52	KTS2422614	41	KTS3227414	71	KTS3265315	65
KTS1891800	52	KTS2424121	39	KTS3227419	71	KTS3265319	65
KTS2010211	57	KTS2424521	39	KTS3228415	71	KTS3265403	65
KTS2010213	57	KTS2424526	39	KTS3228618	72	KTS3265404	65
KTS2010214	57	KTS2426027	41	KTS3229027	72	KTS3265406	65
KTS2010216	57	KTS2426427	41	KTS3230015	72	KTS3265408	65
		KTS2429818	43	KTS3230215	72	KTS3265410	65
		KTS2429821	43	KTS3230327	72	KTS3265412	65
		KTS2432418	42	KTS3230523	72	KTS3265415	65
		KTS2432426	42	KTS3240217	63	KTS3265706	65
		KTS2434413	39	KTS3240913	63	KTS3265712	65



Part No.	Page #	Part No.	Page #	Part No.	Page #	Part No.	Page #
KTS3265719	65	KTS4029712	93	KTS4611201	86	KTS2015002	56
KTS3265725	65	KTS4040517	77	KTS4611202	86	KTS2015003	56
KTS3265732	65	KTS4042001	78	KTS4612035	86	KTS2015006	56
KTS3265738	65	KTS4042201	78	KTS4612050	86	KTS2015009	56
KTS3265806	65	KTS4042202	78	KTS4612075	86	KTS2015012	56
KTS3265812	65	KTS4042203	78	KTS4612090	86	KTS2015015	56
KTS3265819	65	KTS4043201	78	KTS4615416	21	KTS2015080	56
KTS3265825	65	KTS4043202	78	KTS4615418	21	KTS2015090	56
KTS3265832	65	KTS4043301	78	KTS4615720	21	KTS6410214	96
KTS3265838	65	KTS4043401	78	KTS4620915	44	KTS6411015	96
KTS3268625	66	KTS4044101	79	KTS4625121	89	KTS6412015	96
KTS3269045	67	KTS4044201	79	KTS4625800	89	KTS6412019	96
KTS3269145	67	KTS4044301	79	KTS4630418	89	KTS6412024	96
KTS3269090	67	KTS4045602	79	KTS4631301	86	KTS6412416	96
KTS3269222	66	KTS4045603	79	KTS4632319	87	KTS6412420	96
KTS3274101	66	KTS4045604	79	KTS4632519	87	KTS6412425	96
KTS3274224	66	KTS4051201	77	KTS4637415	87	KTS6412430	96
KTS3274624	66	KTS4051202	77	KTS4639004	87	KTS6413015	96
KTS3274822	67	KTS4051203	77	KTS4641518	87	KTS6414815	96
KTS3275903	66	KTS4051204	77	KTS4642322	87	KTS6415020	96
KTS3275906	66	KTS4051205	77	KTS4642516	87	KTS6416016	97
KTS3275909	66	KTS4051801	77	KTS4650618	88	KTS6416018	97
KTS3275912	66	KTS4051802	77	KTS4651616	88	KTS6416020	97
KTS3275915	66	KTS4051803	77	KTS4652021	88	KTS6416024	97
KTS3275919	66	KTS4051804	77	KTS4652205	88	KTS6416816	96
KTS3275925	66	KTS4051805	77	KTS4652207	88	KTS6416820	96
KTS3276118	65	KTS4210103	50	KTS4652209	88	KTS6421023	97
KTS3276120	65	KTS4237613	45	KTS4656018	44	KTS6421123	97
KTS3276123	65	KTS4237813	45	KTS4656218	44	KTS6422022	99
KTS3276125	65	KTS4261599	57	KTS4657123	89	KTS6422622	98
KTS3277171	66	KTS4261600	57	KTS4811501	80	KTS6422625	98
KTS3277214	67	KTS4261612	57	KTS4811502	80	KTS6422628	98
KTS3277215	67	KTS4261634	57	KTS4826514	80	KTS6422722	98
KTS3277620	67	KTS4261656	57	KTS4831000	80	KTS6422725	98
KTS3280110	47	KTS4261678	57	KTS4831010	80	KTS6422728	98
KTS3280118	47	KTS4294005	22	KTS4832000	80	KTS6424723	99
KTS3280724	47	KTS4294008	22	KTS4832010	80	KTS6426024	98
KTS3281106	46	KTS4294205	22	KTS5018319	12	KTS6428226	101
KTS3281108	46	KTS4294208	22	KTS5020219	82	KTS6429330	98
KTS3282717	73	KTS4294407	22	KTS5020319	82	KTS6431015	100
KTS3283030	73	KTS4294607	22	KTS5020518	81	KTS6431021	100
KTS3297327	69	KTS4315001	21	KTS5025228	82	KTS6431029	100
KTS3297530	69	KTS4315003	21	KTS5025903	82	KTS6431035	100
KTS3304056	64	KTS4315005	21	KTS5025955	82	KTS6431215	100
KTS3306019	72	KTS4401000	83	KTS5026000	82	KTS6431221	100
KTS3317300	58	KTS4401001	83	KTS5227208	83	KTS6431229	100
KTS3317301	58	KTS4401030	83	KTS5227210	83	KTS6431235	100
KTS3317500	58	KTS4401035	83	KTS5227212	83	KTS6620100	103
KTS3317501	58	KTS4401040	83	KTS5227214	83	KTS6620108	103
KTS3352219	64	KTS4401050	83	KTS5227216	83	KTS6620110	103
KTS3353217	65	KTS4401065	83	KTS5227218	83	KTS6620112	103
KTS3353318	65	KTS4401075	83	KTS5227220	83	KTS6620114	103
KTS3354016	63	KTS4406001	85	KTS5227222	83	KTS6620116	103
KTS3354215	63	KTS4406002	85	KTS5227224	83	KTS6620118	103
KTS3354517	64	KTS4406007	85	KTS5227226	83	KTS6620120	103
KTS3354522	64	KTS4407001	84	KTS5411111	81	KTS6620122	103
KTS3356614	64	KTS4407100	85	KTS5411114	81	KTS6620124	103
KTS3357218	64	KTS4409601	85	KTS5411314	81	KTS6620126	103
KTS3390214	73	KTS4409602	85	KTS5411615	81	KTS6620128	103
KTS4016101	77	KTS4409603	85	KTS5413115	81	KTS6620130	103
KTS4016102	77	KTS4409650	85	KTS5452015	81	KTS6620132	103
KTS4016103	77	KTS4409710	85	KTS5452020	81	KTS6620134	103
KTS4016104	77	KTS4420614	18	KTS5621020	91	KTS6620136	103
KTS4029506	93	KTS4420616	18	KTS5621025	91	KTS6621314	103
KTS4029507	93	KTS4426108	82	KTS5621475	91	KTS6622300	58
KTS4029508	93	KTS4426114	82	KTS5624218	91	KTS6811015	101
KTS4029509	93	KTS4431040	83	KTS5632515	92	KTS6811515	101
KTS4029510	93	KTS4499401	84	KTS6223101G	103		
KTS4029511	93	KTS4499601	84	KTS6223102G	103		
KTS4029512	93	KTS4499602	84	KTS6223103G	103		
KTS4029514	93	KTS4499801	84	KTS6223104G	103		
KTS4029710	93	KTS4499802	84	KTS2015000	56		



Part No.	Page #	Part No.	Page #	Part No.	Page #	Part No.	Page #
KTS6812100	97	KTS7042005	107	KTS7053255	116	KTS7061704	110
KTS6813019	101	KTS7042007	107	KTS7053256	116	KTS7061723	110
KTS6821022	97	KTS7042009	107	KTS7053257	116	KTS7062101	110
KTS6840021	23	KTS7042011	107	KTS7053258	116	KTS7062700	109
KTS7011001	104	KTS7042013	107	KTS7053259	116	KTS7063101	109
KTS7011002	104	KTS7042015	107	KTS7053260	116	KTS7063102	110
KTS7011003	104	KTS7042017	107	KTS7053261	116	KTS7063103	110
KTS7012001	104	KTS7042019	107	KTS7053262	116	KTS7071508	108
KTS7012002	104	KTS7042021	107	KTS7053263	116	KTS7071509	108
KTS7012003	104	KTS7042023	107	KTS7053264	116	KTS7071511	108
KTS7030800	104	KTS7042025	107	KTS7053265	116	KTS7210230	111
KTS7030900G	104	KTS7042408	107	KTS7053266	116	KTS7215544	111
KTS7035603	105	KTS7042413	107	KTS7053267	116	KTS7412804	58
KTS7040414	107	KTS7042508	107	KTS7053268	116	KTS7412806	58
KTS7040426	107	KTS7042608	106	KTS7053269	116	KTS7433114	57
KTS7040614	107	KTS7042613	106	KTS7053270	116	KTS7433614	57
KTS7041001	107	KTS7043010	106	KTS7053271	116	KTS7434113	58
KTS7041002	107	KTS7043012	106	KTS7053272	116	KTS7436013	57
KTS7041003	107	KTS7043014	106	KTS7053273	116	KTS7436314	57
KTS7041004	107	KTS7043016	106	KTS7053274	116	KTS7436511	57
KTS7041005	107	KTS7043018	106	KTS7053275	116	KTS7436815	58
KTS7041006	107	KTS7043020	106	KTS7053276	116	KTS9900004	80
KTS7041007	107	KTS7043608	106	KTS7053277	117	KTS9900006	69
KTS7041008	107	KTS7043913	106	KTS7053278	117	KTS9911050	69
KTS7041009	107	KTS7043917	106	KTS7053279	117	KTS9930141	104
KTS7041010	107	KTS7043921	106	KTS7053280	117		
KTS7041011	107	KTS7043925	106	KTS7053281	117		
KTS7041012	107	KTS7044708	106	KTS7053282	117		
KTS7041013	107	KTS7044729	106	KTS7053283	117		
KTS7041014	107	KTS7044733	106	KTS7053284	117		
KTS7041015	107	KTS7044737	106	KTS7053285	117		
KTS7041016	107	KTS7044741	106	KTS7053286	117		
KTS7041017	107	KTS7050225	105	KTS7053287	118		
KTS7041018	107	KTS7050425	108	KTS7053288	118		
KTS7041019	107	KTS7050927	105	KTS7053289	118		
KTS7041020	107	KTS7051424	105	KTS7053290	118		
KTS7041021	107	KTS7051524	105	KTS7053291	118		
KTS7041022	107	KTS7051825	105	KTS7053292	118		
KTS7041023	107	KTS7053221	111	KTS7053293	118		
KTS7041024	107	KTS7053222	113	KTS7053294	118		
KTS7041025	107	KTS7053223	113	KTS7053295	118		
KTS7041026	107	KTS7053224	113	KTS7053296	118		
KTS7041027	107	KTS7053225	113	KTS7053297	118		
KTS7041028	107	KTS7053226	113	KTS7053298	118		
KTS7041029	107	KTS7053227	113	KTS7053299	118		
KTS7041030	107	KTS7053228	113	KTS7053300	118		
KTS7041501	107	KTS7053229	113	KTS7053301	118		
KTS7041502	107	KTS7053230	114	KTS7053302	118		
KTS7041503	107	KTS7053231	114	KTS7053303	119		
KTS7041504	107	KTS7053232	114	KTS7053304	119		
KTS7041505	107	KTS7053233	114	KTS7053305	119		
KTS7041506	107	KTS7053234	114	KTS7053306	119		
KTS7041507	107	KTS7053235	114	KTS7053307	119		
KTS7041508	107	KTS7053236	114	KTS7053308	119		
KTS7041509	107	KTS7053237	114	KTS7053309	119		
KTS7041510	107	KTS7053238	114	KTS7053310	119		
KTS7041511	107	KTS7053239	114	KTS7053311	119		
KTS7041512	107	KTS7053240	114	KTS7053312	119		
KTS7041513	107	KTS7053241	114	KTS7053313	119		
KTS7041514	107	KTS7053242	115	KTS7053314	119		
KTS7041515	107	KTS7053243	115	KTS7053315	119		
KTS7041516	107	KTS7053244	115	KTS7053316	119		
KTS7041517	107	KTS7053245	115	KTS7053317	119		
KTS7041518	107	KTS7053246	115	KTS7053318	119		
KTS7041519	107	KTS7053247	115	KTS7053319	120		
KTS7042001	107	KTS7053248	115	KTS7053320	120		
KTS7042003	107	KTS7053249	115	KTS7053321	120		
		KTS7053250	115	KTS7053926	111		
		KTS7053251	115	KTS7058132	108		
		KTS7053252	115	KTS7059412	108		
		KTS7053253	115	KTS7059440	108		
		KTS7053254	116	KTS7059640	108		



Aesculap	K.T Surgico	Aesculap	K.T Surgico	Aesculap	K.T Surgico	Aesculap	K.T Surgico	Aesculap	K.T Surgico
AA050R	KTS0610301	BC277R	KTS0828723	BD512R	KTS1014412	BH205R	KTS1241118	BM124R	KTS2441414
AA804R	KTS0234015	BC289R	KTS0828530	BD520R	KTS1018615	BH226R	KTS1245124	BM201R	KTS2413013
AA821C	KTS0235015	BC314R	KTS0810214	BD537R	KTS1013014	BH231R	KTS1246128	BM202R	KTS2410013
AC023C	KTS0228817	BC315R	KTS0810215	BD539R	KTS1013016	BH332R	KTS1218113	BM205R	KTS2414013
AC023C	KTS0228820	BC316R	KTS0810216	BD557R	KTS1012014	BH335R	KTS1233115	BM215R	KTS2415014
AN379R	KTS0419017	BC324R	KTS0811014	BD557R	KTS1012214	BH335R	KTS1218115	BM218R	KTS2416115
AN380R	KTS0419020	BC341R	KTS0800411	BD559R	KTS1012215	BH422R	KTS1231213	BM220R	KTS2416015
AN381R	KTS0419025	BC344R	KTS0812014	BD559R	KTS1012016	BH442R	KTS1231016	BM220R	KTS2416018
AN404R	KTS0421000	BC544R	KTS0816014	BD561R	KTS1013020	BH443R	KTS1231316	BM235R	KTS2418014
AN405R	KTS0424119	BC544W	KTS0716014	BD562R	KTS1012225	BH443R	KTS1231116	BM237R	KTS2418024
AN406R	KTS0422501	BC545R	KTS0816015	BD562R	KTS1012025	BH444R	KTS1231018	BM245R	KTS2418016
AN407R	KTS0422502	BC547R	KTS0816017	BD610R	KTS1024618	BH445R	KTS1231118	BM247R	KTS2418018
AN410R	KTS0423901	BC547W	KTS0716017	BD611R	KTS1024818	BH446R	KTS1231020	BM249R	KTS2418020
AN447	KTS0425001	BC554R	KTS0816114	BD660R	KTS1027215	BH447R	KTS1231120	BM259R	KTS2424121
AN447R	KTS0423902	BC557R	KTS0816117	BD662R	KTS1031018	BH449R	KTS1231122	BM272R	KTS2426027
AN460R	KTS0424128	BC557W	KTS0716117	BD693R	KTS1045118	BH451R	KTS1231124	BM290R	KTS2442014
AN490R	KTS0421100	BC558R	KTS0816423	BD693R	KTS1132218	BH473R	KTS1231126	BM291R	KTS2442017
AN491R	KTS0421001	BC559R	KTS0816523	BD700R	KTS1018812	BH474R	KTS1231030	BM562R	KTS2466114
AN492R	KTS0421002	BC560R	KTS0816023	BD700R	KTS1019912	BH475R	KTS1231130	BM563R	KTS2482214
AN493R	KTS0421003	BC562R	KTS0816019	BD701R	KTS1018915	BH618R	KTS1232014	BM563R	KTS2482414
AN494R	KTS0421004	BC563R	KTS0816119	BD733R	KTS1032215	BH630R	KTS1232216	BM564R	KTS2482514
AN560R	KTS0423801	BC566R	KTS0818017	BD755R	KTS1025615	BH642R	KTS1232016	BM606R	KTS2422415
AN590R	KTS0420401	BC567R	KTS0816128	BD755R	KTS1045615	BH643R	KTS1232116	BMO19R	KTS2416420
AN591R	KTS0421101	BC568R	KTS0816028	BD760R	KTS1025620	BH644R	KTS1232018	BMO34R	KTS2421420
AN592R	KTS0421102	BC569R	KTS0818117	BD760R	KTS1045623	BH645R	KTS1232118	BMO44R	KTS2422614
AN593R	KTS0421103	BC584R	KTS0817017	BD801R	KTS1032023	BH646R	KTS1232020	BN034R	KTS2010214
AN594R	KTS0421104	BC587R	KTS0817117	BD878R	KTS4420614	BH647R	KTS1232120	BN036R	KTS2010216
AN710R	KTS6840021	BC601R	KTS0828015	BD880R	KTS4420616	BH648R	KTS1232022	BN136R	KTS2010616
AN909R	KTS1271218	BC602R	KTS0828018	BD883R	KTS4615720	BH650R	KTS1232024	BT004R	KTS1816414
BB070R	KTS0610116	BC605R	KTS0828115	BD952R	KTS1080209	BH652R	KTS1232026	BT006R	KTS1816215
BB073R	KTS0610300	BC605R	KTS0828914	BD980R	KTS1081210	BH652R	KTS1232030	BT007R	KTS1816315
BB075R	KTS0610302	BC605R	KTS0828111	BD980R	KTS1081610	BH801R	KTS1234116	BT008R	KTS1816017
BB077R	KTS0610700	BC605W	KTS0728115	BD980R	KTS1080211	BH805R	KTS1234120	BT041R	KTS1813121
BB079R	KTS0610303	BC606R	KTS0828118	BD997R	KTS1082815	BH902R	KTS1235815	BT0700R	KTS1816117
BB084R	KTS0610400	BC606W	KTS0728118	BF112R	KTS1617018	BH959R	KTS5018319	BT111R	KTS1822001
BC023R	KTS0828315	BC607R	KTS0828120	BF113R	KTS1617118	BJ012R	KTS1247116	BT182R	KTS1823421
BC050R	KTS0838410	BC615R	KTS0828123	BF120R	KTS1617225	BJ013R	KTS1241514	BT202R	KTS1826515
BC056R	KTS0834009	BC615W	KTS0728123	BF121R	KTS1617325	BJ014R	KTS1241518	BT240R	KTS3283030
BC057R	KTS0834109	BC618R	KTS0828328	BF122R	KTS1617025	BJ057R	KTS1243520	BT241R	KTS1828101
BC060R	KTS0834010	BC631R	KTS0828318	BF122R	KTS1617020	BJ058R	KTS1243524	BT242R	KTS1828002
BC061R	KTS0834110	BC647R	KTS0848025	BF123R	KTS1617125	BJ059R	KTS1243528	BT243R	KTS1828003
BC061R	KTS0838510	BC648R	KTS0848040	BF123R	KTS1617120	BJ061R	KTS1247922	BT244R	KTS1828004
BC061W	KTS0734110	BC652R	KTS0848060	BF310R	KTS5452015	BJ100R	KTS1243113	BT246R	KTS1828006
BC064R	KTS0834011	BC654R	KTS0803012	BF407R	KTS1410309	BJ102R	KTS1243318	BT252R	KTS1828102
BC065R	KTS0834111	BC662R	KTS0803025	BF426R	KTS1414016	BJ524R	KTS1260021	BT253R	KTS1828103
BC065W	KTS0734111	BC663R	KTS0803045	BF431R	KTS1411109	BJ532R	KTS1242321	BT254R	KTS1828104
BC106R	KTS0859411	BC664R	KTS0803060	BF432R	KTS1411111	BM001R	KTS2467414	BT256R	KTS1828106
BC106R	KTS0834211	BC801R	KTS0880209	BF433R	KTS1411113	BM002R	KTS2467614	BT264R	KTS1833304
BC107R	KTS0834312	BC802R	KTS0880111	BF433R	KTS1411115	BM004R	KTS2467618	BT266R	KTS1833306
BC112R	KTS0837602	BC815R	KTS0886212	BF448R	KTS1411713	BM008R	KTS2411412	BT279R	KTS1836608
BC113R	KTS0837702	BC848R	KTS0890115	BF458R	KTS1412511	BM010R	KTS2413413	BT284R	KTS1836440
BC131R	KTS0800511	BC860R	KTS0890111	BF459R	KTS1412513	BM011R	KTS2411612	BT408R	KTS1834822
BC163R	KTS0809711	BC861R	KTS0890114	BF460R	KTS1416514	BM012R	KTS2414413	BT420R	KTS3282717
BC163R	KTS0809714	BC862R	KTS0890118	BF461R	KTS1412110	BM013R	KTS2416215	BT470R	KTS1814000
BC165R	KTS0837111	BC863R	KTS0890120	BF461R	KTS1412111	BM014R	KTS2410413	BT471R	KTS1814001
BC171R	KTS0836111	BC881R	KTS0895018	BF463R	KTS1412113	BM015R	KTS2415412	BT472R	KTS1814002
BC177R	KTS0800315	BC881R	KTS0895016	BF464R	KTS1412313	BM016R	KTS2416415	BT475R	KTS1815231
BC180R	KTS0836314	BC885R	KTS0895014	BF508R	KTS1325325	BM017R	KTS2416418	BT475R	KTS1815001
BC180R	KTS0832115	BC885R	KTS0895015	BH104R	KTS1222210	BM020R	KTS2432418	BT476R	KTS1815251
BC181R	KTS0832118	BD027R	KTS1011014	BH104R	KTS1221010	BM021R	KTS2416430	BT476R	KTS1815002
BC194R	KTS0836210	BD029R	KTS1011016	BH105R	KTS1222310	BM023R	KTS2432426	BT477R	KTS1815003
BC194R	KTS0836211	BD045R	KTS1010213	BH105R	KTS1221109	BM029R	KTS2416615	BT478R	KTS1815004
BC195R	KTS0836310	BD045R	KTS1010414	BH105R	KTS1221110	BM035R	KTS2421418	BT481R	KTS1815851
BC195R	KTS0836311	BD047R	KTS1010214	BH108R	KTS1222212	BM036R	KTS2421423	BT482R	KTS1815626
BC195W	KTS0736311	BD047R	KTS1010415	BH109R	KTS1222312	BM037R	KTS2421426	BT483R	KTS1815876
BC224R	KTS0811414	BD049R	KTS1010216	BH110R	KTS1222012	BM038R	KTS2421431	BT486R	KTS1852075
BC240R	KTS0817415	BD050R	KTS1010418	BH110R	KTS1222014	BM042R	KTS2426427	BT551R	KTS1819140
BC242R	KTS0817417	BD051R	KTS1011020	BH111R	KTS1222112	BM044R	KTS2422413	BT552R	KTS1819162
BC243R	KTS0817517	BD052R	KTS1010225	BH111R	KTS1222114	BM045R	KTS2422418	BT553R	KTS1819164
BC250R	KTS0816414	BD052R	KTS1010425	BH120R	KTS1223012	BM046R	KTS2422420	BT602R	KTS1818325
BC251R	KTS0816417	BD151R	KTS1000112	BH120R	KTS1223014	BM047R	KTS2422426	BT605R	KTS1818750
BC251R	KTS0716517	BD151R	KTS1017412	BH121R	KTS1223112	BM054R	KTS2422213	BT667R	KTS1837825
BC252R	KTS0816514	BD160R	KTS1022502	BH121R	KTS1223114	BM055R	KTS2422215	BT751R	KTS1816812
BC253R	KTS0816517	BD220R	KTS1016812	BH130R	KTS1227014	BM056R	KTS2422218	BT752R	KTS1816817
BC257R	KTS0829311	BD220R	KTS1017612	BH131R	KTS1227114	BM065R	KTS2418414	BT755R	KTS1816825
BC257R	KTS0837511	BD222R	KTS1017012	BH134R	KTS1224014	BM065R	KTS2418416	BT756R	KTS1817025
BC258R	KTS0828615	BD223R	KTS1017615	BH135R	KTS1224114	BM066R	KTS2418418	BT761R	KTS1817032
BC259R	KTS0828715	BD225R	KTS1015212	BH144R	KTS1225014	BM066R	KTS2419418	BT762R	KTS1817040
BC260R	KTS0828415	BD238R	KTS1022017	BH145R	KTS1225114	BM067R	KTS2418420	BT763R	KTS1817050
BC261R	KTS0828515	BD285R	KTS1030218	BH154R	KTS1226014	BM067R	KTS2418424	BV010R	KTS1862007
BC262R	KTS0800217	BD505R	KTS4294607	BH155R	KTS1226114	BM069R	KTS2418426		
BC263R	KTS0800117	BD506R	KTS1016912	BH164R	KTS1224216	BM069R	KTS2418430		
BC263R	KTS0828518	BD507R	KTS1017312	BH165R	KTS1224316	BM076R	KTS2422423		
BC267R	KTS0828523	BD510R	KTS1018612	BH166R	KTS1226216	BM081R	KTS2424521		
BC271R	KTS0828718	BD511R	KTS1018012	BH167R	KTS1226316	BM084R	KTS2424526		



Aesculap	K.T Surgico	Aesculap	K.T Surgico	Aesculap	K.T Surgico	Aesculap	K.T Surgico	Aesculap	K.T Surgico
BV010R	KTS1862107	EF809R	KTS6620126	FK180R	KTS3277171	LX164R	KTS3304056	MD594R	KTS4652207
BV015R	KTS1863410	EF810R	KTS6620128	FK185R	KTS3280110	LX172R	KTS3352219	MD595R	KTS4652205
BV016R	KTS1863510	EF811R	KTS6620130	FK310R	KTS4642516	LX174R	KTS3353318	MD596R	KTS4652209
BV030R	KTS1862410	EF812R	KTS6620132	FK410R	KTS3275906	LX189R	KTS3353217	MD597R	KTS4625800
BV032R	KTS1862910	EF813R	KTS6620134	FK411R	KTS3275909	LX550R	KTS0891320	MD598R	KTS4657123
BV070R	KTS1865010	EF814R	KTS6620136	FK412R	KTS3275912	LX552R	KTS0891322	MD613R	KTS4630418
BV071R	KTS1867013	EL192R	KTS7012002	FK413R	KTS3275919	LX553R	KTS0891123	MD615R	KTS1821301
BV072R	KTS1867016	EL193R	KTS7012003	FK414R	KTS3275925	LX558R	KTS0893723	MD617R	KTS1821302
BV073R	KTS1865110	EL194R	KTS7011001	FK505R	KTS3297327	LX605R	KTS3036623	MD652R	KTS1821303
BV074R	KTS1867113	EL195R	KTS7011002	FK505R	KTS3297530	MB213R	KTS2469521	MD653R	KTS4656018
BV075R	KTS1867116	EL196R	KTS7011003	FK630R	KTS3246996	MB223R	KTS0890121	OA334R	KTS4656218
BV085R	KTS1800416	EM244R	KTS7042408	FK631R	KTS3247497	MB236R	KTS4615416	OA335R	KTS1822816
BV085R	KTS1800513	EO012R	KTS7058132	FK631R	KTS3246997	MB238R	KTS1334318	OA337R	KTS1821601
BV085R	KTS1800516	EO108R	KTS7050425	FK632R	KTS3247498	MB275R	KTS3276118	OA351R	KTS4620915
BV085R	KTS1800520	EO110R	KTS7050225	FK632R	KTS3246998	MB276R	KTS3276120	OB510	KTS1821602
BV101R	KTS1865413	EO125R	KTS7051825	FK633R	KTS3247500	MB277R	KTS3276123	OB511	KTS1821605
BV104R	KTS1865513	EO150R	KTS7051424	FK633R	KTS3247000	MB278R	KTS3276125	OB512	KTS1821607
BV200R	KTS1867020	EO152R	KTS7051524	FK634R	KTS3247501	MB291R	KTS1015412	OB513	KTS4826514
BV201R	KTS1867120	EO202R	KTS7053221	FK634R	KTS3247001	MB576R	KTS3306019	OB514	KTS5227208
BV206R	KTS1867025	ER068R	KTS7062700	FK635R	KTS3247502	MB603R	KTS2010213	OB515	KTS5227210
BV207R	KTS1867125	ER068R	KTS7062101	FK635R	KTS3247002	MB641R	KTS3265706	OC051R	KTS5227212
BV208R	KTS1864826	ER080R	KTS7061704	FK636R	KTS3247503	MB642R	KTS3265712	OC093R	KTS5227214
BV209R	KTS1864926	ER080R	KTS7061723	FK636R	KTS3247003	MB643R	KTS3265719	OC094R	KTS5227216
BV250R	KTS1867820	ET107R	KTS7053926	FK637R	KTS3247504	MB644R	KTS3265725	OC094R	KTS5227218
BV250R	KTS1867920	ET152R	KTS0855010	FK637R	KTS3247004	MB645R	KTS3265732	OF028	KTS5227220
BV271R	KTS1801031	ET220R	KTS7210230	FK638R	KTS3247505	MB646R	KTS3265738	OF120C	KTS5227222
BV271R	KTS1801045	ET237R	KTS7215544	FK638R	KTS3247005	MB651R	KTS3265806	OF411R	KTS5227224
BV272R	KTS1800731	FA005R	KTS5411114	FK639R	KTS3247006	MB652R	KTS3265812	OG331R	KTS5227226
BV272R	KTS1869031	FA007R	KTS5411314	FK814R	KTS3249600	MB653R	KTS3265819	OG332R	KTS4832000
BV273R	KTS1800631	FB400R	KTS1130816	FK845R	KTS3249213	MB654R	KTS3265825	OG350R	KTS4832010
BV276R	KTS1869132	FB402R	KTS1130820	FK846R	KTS3249216	MB655R	KTS3265832	OK060R	KTS4831000
BV297R	KTS4040517	FB404R	KTS1130824	FK848R	KTS3249220	MB656R	KTS3265838	OK061R	KTS4831010
BV601R	KTS1876690	FB405R	KTS1130830	FL036R	KTS3277214	MB736R	KTS5411615	OK091R	KTS4811502
BV605R	KTS1876418	FB406R	KTS1131020	FL037	KTS3277201	MB768R	KTS1826513	OK092R	KTS5020319
BV606R	KTS1876545	FB407R	KTS1131024	FL038R	KTS3268625	MB785R	KTS1830419	OK093R	KTS5020518
BV610R	KTS1876425	FB408R	KTS1131030	FL038R	KTS3274624	MD076R	KTS1817612	OK094R	KTS0415027
BV611R	KTS1876570	FB411R	KTS1130616	FL044R	KTS3274822	MD093R	KTS4029506	OK350R	KTS0842717
BV632R	KTS1891800	FB411R	KTS1130613	FL066R	KTS3274224	MD094R	KTS4029507	OK366R	KTS5025903
BV996R	KTS1880214	FB412R	KTS1130620	FL066R	KTS3269045	MD095R	KTS4029508	OK369R	KTS5025955
BV997R	KTS1880218	FB414R	KTS1130416	FO115R	KTS3240913	MD096R	KTS4029509	OK400R	KTS5025228
DF114R	KTS3277620	FB415R	KTS1130420	FO123R	KTS3241113	MD097R	KTS4029510	OK405R	KTS5026000
DO208R	KTS0834411	FB415R	KTS1133220	FO191R	KTS3244120	MD098R	KTS4029511	OK525R	KTS1875400
DO209R	KTS0834511	FB416R	KTS1130424	FO407R	KTS3210112	MD099R	KTS4029512	OL120R	KTS1025625
DO771R	KTS9900006	FB417R	KTS1130430	FO407R	KTS3210114	MD127R	KTS4029710	OL152R	KTS7043012
DP031R	KTS3356614	FB418R	KTS1130624	FO410R	KTS3215613	MD129R	KTS4029712	OL165R	KTS7043014
DP551R	KTS0883112	FB421R	KTS1209307	FO411R	KTS3213515	MD452R	KTS1271018	OL354R	KTS7043016
DP552R	KTS0883512	FB504R	KTS1325327	FO450R	KTS3215418	MD453R	KTS1271020	OL406R	KTS7043018
EA010R	KTS6414815	FB506R	KTS1325320	FO451R	KTS3215518	MD455R	KTS3240217	OL407R	KTS7043020
EA012R	KTS6410214	FB507R	KTS1325324	FO457R	KTS9911050	MD456R	KTS2419818	OL409R	KTS0415222
EA015R	KTS6411015	FB540R	KTS1310212	FO500R	KTS3221715	MD456R	KTS2480818	OL427R	KTS3214220
EA016R	KTS6412015	FB543R	KTS1311112	FO500R	KTS3221315	MD463R	KTS1023017	OL429R	KTS3215116
EA017R	KTS6412019	FB600R	KTS5632515	FO511R	KTS3216717	MD464R	KTS1228014	OL429R	KTS1270416
EA018R	KTS6412024	FB704R	KTS1342114	FO514R	KTS3220118	MD465R	KTS1228114	OL460R	KTS1270418
EA020R	KTS6415020	FB803R	KTS5621020	FO515R	KTS3220518	MD466R	KTS1228016	OL470R	KTS1821201
EA030R	KTS6416016	FB804R	KTS5621025	FO518R	KTS3221718	MD467R	KTS1228116	OL477R	KTS1821202
EA031R	KTS6416018	FB830R	KTS5621475	FO522R	KTS3225902	MD470R	KTS6416024	OL601R	KTS4237613
EA032R	KTS6416020	FC411R	KTS1084015	FO527R	KTS3224723	MD472R	KTS7059440	OL602R	KTS1821413
EA037R	KTS1235820	FD046R	KTS2429818	FO530R	KTS3217523	MD473R	KTS3261397	OL603R	KTS4261599
EA050R	KTS6416816	FD233R	KTS2466314	FO532R	KTS3224523	MD474R	KTS3261398	OL611R	KTS4261600
EA052R	KTS6416820	FF291R	KTS4016101	FO534R	KTS3223723	MD475R	KTS3261300	OL612R	KTS4261612
EA100R	KTS6412430	FF292R	KTS4016102	FO550R	KTS3223022	MD476R	KTS3261301	OL619R	KTS4261634
EA127R	KTS6421123	FF293R	KTS4016103	FO551R	KTS3223122	MD477R	KTS3261302	OL620R	KTS4261656
EA150R	KTS6422022	FF294R	KTS4016104	FO555R	KTS3226122	MD478R	KTS3261303	OL621R	KTS4261678
EA202R	KTS6422622	FF677R	KTS4051202	FO620R	KTS3227414	MD479R	KTS3261304	OL622R	KTS4294005
EA203R	KTS6422625	FF678R	KTS4051203	FO622R	KTS3228415	MD480R	KTS3261305	OL623R	KTS4315001
EA204R	KTS6422628	FF679R	KTS4051204	FO624R	KTS3227419	MD481R	KTS3261306	OL624R	KTS4315005
EA831R	KTS6812100	FF680R	KTS4051205	FO640R	KTS3228618	MD482R	KTS3261297	OM038R	KTS4315003
EA880R	KTS6821022	FF715R	KTS4051201	FO645R	KTS3230523	MD483R	KTS3261298	OM079R	KTS4407001
ED0082R	KTS1623126	FF737R	KTS4051802	FO646R	KTS3229027	MD484R	KTS3261200	OM080R	KTS4407100
EF016R	KTS1259618	FF739R	KTS4051803	FO649R	KTS3230327	MD485R	KTS3261201	OM081R	KTS4426108
EF314C	KTS6621314	FF745R	KTS4051804	FO900R	KTS3317500	MD486R	KTS3261202	OM082R	KTS4499801
EF800R	KTS6620108	FF751R	KTS4051805	FO901R	KTS3317501	MD487R	KTS3261203	OM083R	KTS4499802
EF801R	KTS6620110	FF799R	KTS4042001	FO920	KTS3245001	MD488R	KTS3261204	OM084R	KTS4431040
EF802R	KTS6620112	FF803R	KTS4042201	FO950R	KTS1299215	MD489R	KTS3261205	OM085R	KTS4611201
EF803R	KTS6620114	FF805R	KTS4043201	GF860C	KTS0414008	MD490R	KTS3261206	OM086R	KTS4611202
EF804R	KTS6620116	FF806R	KTS4043202	GF862C	KTS0413010	MD491R	KTS0865022	OM087R	KTS4612035
EF805R	KTS6620118	FF807R	KTS4044201	GF935R	KTS0440310	MD497R	KTS1813519	OM088R	KTS4612050
EF806R	KTS6620120	FH116R	KTS0642028	GF936R	KTS0440315	MD499R	KTS2416423	OM100R	KTS4612075
EF807R	KTS6620122	FK117R	KTS3281106	GF936R	KTS0440320	MD568R	KTS6413015	OM104R	KTS4612090
EF808R	KTS6620124	FK118R	KTS3281108	HF213R	KTS7433114	MD570R	KTS1245918	OM110R	KTS0840918
		FK146R	KTS3276409	HF234R	KTS7433614	MD579R	KTS7043010	OM116R	KTS0806109
		FK147R	KTS3276413	HF243R	KTS7434113	MD580R	KTS4294205	OM237R	KTS0837701
		FK148R	KTS3276419	LX150R	KTS3354215	MD581R	KTS1821610	OM611R	KTS4632519
		FK171R	KTS3280118	LX157R	KTS3354522	MD582R	KTS1236015	OM612R	KTS4632319
		FK176R	KTS3280724	LX158R	KTS3357218	MD591R	KTS1259420	OM671R	KTS4631301
				LX159R	KTS3354517	MD592R	KTS1259520	OM700R	KTS4637415
				LX163R	KTS3354016	MD593R	KTS7043506	OM750R	KTS4642322
								OM751R	KTS4641518
								OM761R	KTS4639004
								OM781R	KTS4652021
								US100R	KTS4650618
								XG006R	KTS4651616



Codman	K.T Surgico	Codman	K.T Surgico	Codman	K.T Surgico	Codman	K.T Surgico	Codman	K.T Surgico
36-1002	KTS2429818	38-6531	KTS1821201	50-4124	KTS1815231	54-6630	KTS0910211	50-4303	KTS1818325
36-1001	KTS2429821	38-6530	KTS1821202	50-4131	KTS1815251	54-6513	KTS0910411	50-4307	KTS1818750
64-6502	KTS2466114	38-1040	KTS1821413	50-4115	KTS1815851	36-6000	KTS1000112	50-8100	KTS1819140
18-1839	KTS2466314	38-6642	KTS1821601	50-4117	KTS1815876	30-4133	KTS1010213	50-8102	KTS1819162
36-1008	KTS2469421	38-6645	KTS1821602	50-2608	KTS1816414	30-4134	KTS1010214	50-1020	KTS1823421
36-1001	KTS2469521	38-6646	KTS1821605	50-4040	KTS1830419	30-4135	KTS1010216	50-4220CO	KTS1828002
18-1828	KTS2482214	38-6647	KTS1821607	50-4246	KTS1836440	30-4138	KTS1010225	50-4221CO	KTS1828003
18-1832	KTS2482514	38-6648	KTS1821610	50-1066	KTS1852075	30-4152	KTS1012014	50-4222CO	KTS1828004
54-8017	KTS0803025	38-1038	KTS1822001	23-1043	KTS3261200	30-4153	KTS1012016	50-4223CO	KTS1828006
54-8018	KTS0803045	50-4265	KTS1862007	23-1044	KTS3261201	30-4156	KTS1012025	50-4225CO	KTS1828102
54-8019	KTS0803060	50-4265	KTS1862107	23-1045	KTS3261202	30-1185	KTS1017012	504226CO	KTS1828103
54-8021	KTS0803090	50-1226	KTS1862410	23-1046	KTS3261203	30-4214	KTS1017612	50-4227CO	KTS1828104
36-5050	KTS0816414	50-1227	KTS1862910	23-1047	KTS3261204	30-1186	KTS1018012	50-4228CO	KTS1828106
36-5051	KTS0816417	50-1186	KTS1867820	23-1048	KTS3261205	30-1187	KTS1018212	50-4222	KTS1833304
36-5052	KTS0816423	50-1165	KTS1869031	23-1049	KTS3261206	30-4218	KTS1018612	50-4223	KTS1833306
36-5060	KTS0816514	50-1163	KTS1869132	23-1041	KTS3261297	30-1189	KTS1018812	50-4297	KTS1834822
36-5061	KTS0816517	36-1010	KTS2411412	23-1042	KTS3261298	30-1151	KTS1022017	50-5646	KTS1865010
36-5062	KTS0816523	64-6558	KTS2413013	23-1053	KTS3261300	36-6005	KTS1022502	50-5647	KTS1865110
34-4092	KTS0817017	36-1016	KTS2413413	23-1054	KTS3261301	30-1142	KTS1023017	50-1194	KTS1867013
34-4097	KTS0817117	36-1015	KTS2414413	23-1055	KTS3261302	30-4152	KTS1026114	50-1196	KTS1867016
36-5050	KTS0817415	64-2509	KTS2416115	23-1056	KTS3261303	30-4153	KTS1026116	50-1198	KTS1867020
36-5051	KTS0817417	36-2000	KTS2416415	23-1057	KTS3261304	30-4155	KTS1026120	50-1200	KTS1867025
36-5061	KTS0817517	36-2001	KTS2416418	23-1058	KTS3261305	30-4156	KTS1026125	50-1195	KTS1867113
36-5021	KTS0828615	36-2002	KTS2416420	23-1059	KTS3261306	30-4157	KTS1026130	50-1197	KTS1867116
54-6520	KTS0836210	36-2003	KTS2416423	23-1051	KTS3261397	30-1214	KTS1030218	50-1199	KTS1867120
54-6522	KTS0836211	36-3010	KTS2421418	23-1052	KTS3261398	30-4255	KTS1081610	50-1201	KTS1867125
54-6521	KTS0836310	36-3011	KTS2421423	40-6771	KTS4639004	32-4000	KTS1222012	50-4400	KTS1875400
54-6523	KTS0836311	36-3000	KTS2422413	65-6620	KTS4641518	32-4001	KTS1222112	50-4500	KTS1876418
54-6517	KTS0837111	36-3001	KTS2422415	52-1005	KTS4642516	30-4105	KTS1223012	50-4516	KTS1876545
54-8001	KTS0848025	36-3002	KTS2422418	88-4000	KTS0228820	30-4106	KTS1223112	50-4522	KTS1876570
54-8002	KTS0848040	36-3003	KTS2422420	70-4001	KTS0413010	32-4020	KTS1224014	50-5651	KTS1876690
54-8003	KTS0848060	36-3004	KTS2422423	70-4000	KTS0414008	32-4021	KTS1224114	50-1363	KTS1877022
54-1075	KTS0859411	36-3005	KTS2422426	70-6560	KTS0415222	32-4050	KTS1224216	50-1212	KTS1880214
30-4153	KTS1012215	36-3020	KTS2432418	30-4340	KTS0419025	32-4051	KTS1224316	50-1210	KTS1880218
30-4156	KTS1012225	36-3021	KTS2432426	37-4107	KTS0440315	32-4030	KTS1225014	50-1205	KTS1882018
30-4152	KTS1013014	36-4000	KTS2441414	54-1075	KTS0800411	32-4031	KTS1225114	50-1206	KTS1882118
30-4153	KTS1013016	36-1006	KTS2467414	54-1076	KTS0800511	32-4040	KTS1226216	50-4509	KTS1891800
30-4215	KTS1015212	36-1000	KTS2467614	54-4085	KTS0801918	32-4041	KTS1226316	49-4030	KTS2010211
30-4216	KTS1015412	36-1002	KTS2467618	54-5573	KTS0802718	32-4070	KTS1231016	49-4031	KTS2010213
30-4220	KTS1025615	181832ATC	KTS2467714	36-5034	KTS0802723	32-4080	KTS1231018	49-4032	KTS2010214
30-4221	KTS1025620	53-1133	KTS3217523	54-4028	KTS0810214	32-4090	KTS1231020	49-4033	KTS2010216
30-4222	KTS1025625	53-1100	KTS3220118	54-4034	KTS0810216	32-4071	KTS1231116	37-4029	KTS2410413
30-4517	KTS1032023	53-1095	KTS3220518	54-4029	KTS0811014	32-4081	KTS1231118	64-2511	KTS2415014
30-1129	KTS1045118	30-1310	KTS3227414	36-5040	KTS0811414	32-4091	KTS1231120	36-2010	KTS2415412
30-4255	KTS1081210	30-1290	KTS3228618	54-4030	KTS0812014	30-4565	KTS1231124	36-2005	KTS2416430
30-4258	KTS1082815	79-6675	KTS4611201	54-4098	KTS0816014	30-4573	KTS1231126	64-2512	KTS2418016
37-1000	KTS1130416	79-6676	KTS4611202	54-4090	KTS0816017	30-8093	KTS1231130	64-2513	KTS2418018
37-1001	KTS1130420	79-7011	KTS4612050	54-5535	KTS0816023	32-4110	KTS1232016	64-2503	KTS2418020
37-1002	KTS1130424	79-7013	KTS4612075	54-4099	KTS0816114	32-4120	KTS1232018	36-2015	KTS2418414
37-1003	KTS1130430	32-7000	KTS6411015	54-4095	KTS0816117	32-4130	KTS1232020	36-2016	KTS2418416
30-8131	KTS1130616	32-7001	KTS6412015	54-4075	KTS0818017	30-4066	KTS1232024	36-2017	KTS2418418
30-8132	KTS1130620	32-7010	KTS6412019	54-4076	KTS0818117	30-8090	KTS1232026	36-2018	KTS2418420
30-8133	KTS1130624	30-4542	KTS6413015	54-1020	KTS0828015	32-4111	KTS1232116	36-2019	KTS2418426
37-1010	KTS1130816	32-7020	KTS6416016	54-1026	KTS0828018	32-4121	KTS1232118	36-2020	KTS2418430
37-1011	KTS1130820	30-7025	KTS6416020	54-1023	KTS0828111	32-4131	KTS1232120	36-3023	KTS2422213
37-1012	KTS1130824	30-7030	KTS6416024	54-1021	KTS0828115	30-4095	KTS1236015	36-3014	KTS2422215
37-1013	KTS1130830	30-5605	KTS7053221	54-1027	KTS0828118	30-1011	KTS1241118	36-3012	KTS2422218
37-1001	KTS1131020	11-5530	KTS0610300	54-1034	KTS0828123	37-1521	KTS1242321	64-2595	KTS2424121
37-1002	KTS1131024	11-5531	KTS0610301	54-5573	KTS0828318	19-5579	KTS1243113	36-4010	KTS2424521
37-1003	KTS1131030	11-5535	KTS0610302	36-5010	KTS0828415	19-5580	KTS1243318	36-4020	KTS2426427
37-1044	KTS1209307	11-5531	KTS0610400	36-5015	KTS0828515	30-4080	KTS1247116	56-2510	KTS3200218
30-4471	KTS1222210	11-5534	KTS0610700	36-5016	KTS0828518	30-4485	KTS1247614	53-1040	KTS3214220
30-4473	KTS1222212	37-1070	KTS1310212	36-5017	KTS0828523	19-8010	KTS1270416	53-1108	KTS3215418
30-4472	KTS1222310	37-1072	KTS1311112	36-5031	KTS0828715	19-8010	KTS1270418	53-1311	KTS3221718
30-4474	KTS1222312	37-1202	KTS1325320	36-5033	KTS0828718	30-2540	KTS1299215	53-1130	KTS3223022
32-4011	KTS1227114	37-1201	KTS1325324	36-5035	KTS0828723	33-5016	KTS1412110	53-1131	KTS3223122
30-4035	KTS1228116	37-1200	KTS1325325	54-6500	KTS0834010	33-5015	KTS1412113	53-1125	KTS3224523
30-4462	KTS1231316	37-1193	KTS1325326	34-6500	KTS0834011	19-4146	KTS1416514	53-1127	KTS3224723
32-4141	KTS1234116	37-1194	KTS1325327	54-6503	KTS0834110	33-6030	KTS1617020	53-1045	KTS3225902
32-4151	KTS1234120	37-4019	KTS1334318	34-6503	KTS0834111	33-6031	KTS1617120	30-1295	KTS3229027
19-4000	KTS1235815	36-5050	KTS0716014	36-5000	KTS0834411	30-5561	KTS1623126	30-1298	KTS3230327
19-4001	KTS1235820	36-5051	KTS0716017	36-5001	KTS0834511	50-1167	KTS1800416	23-1041	KTS3246997
25-836	KTS1243520	36-5061	KTS0716117	36-5006	KTS0837511	50-1166	KTS1800516	23-1042	KTS3246998
19-8145	KTS1243524	36-5061	KTS0716517	54-6524	KTS0837602	50-4071	KTS1813121	23-1043	KTS3247000
25-838	KTS1243528	36-5010	KTS0728115	54-6527	KTS0837701	50-4120	KTS1815001	23-1044	KTS3247001
30-8059	KTS1245124	36-5011	KTS0728118	54-6525	KTS0837702	50-4121	KTS1815002	23-1045	KTS3247002
30-4700	KTS1246128	36-5017	KTS0728123	54-6666	KTS0840918	50-4122	KTS1815003	23-1046	KTS3247003
30-5800	KTS1260021	36-5001	KTS0734111	54-5545	KTS0855010	50-4111	KTS1815626	23-1047	KTS3247004
19-1016	KTS1410309	30-4200	KTS1018915	54-6624	KTS0880209	50-4260	KTS1816017	23-1048	KTS3247005
33-5000	KTS1411109	30-1132	KTS1024618	54-4110	KTS0883112	50-4261	KTS1816117	23-1049	KTS3247006
33-5005	KTS1411113	30-1132	KTS1024818	54-4012	KTS0890115	50-5637	KTS1816812		
19-4012	KTS1411713	30-1210	KTS1031018	54-4013	KTS0890118	50-5638	KTS1816817		
30-985	KTS1414016	30-8112	KTS1045623	54-4192	KTS0890121	50-4210	KTS1817025		
50-1160	KTS1801031	50-4095	KTS1813519	54-4192	KTS0891320	50-4212	KTS1817040		
50-1162	KTS1801045	50-4135	KTS1814000	34-4013	KTS0895018	50-4214	KTS1817050		



Codman	K.T Surgico	Codman	K.T Surgico
88-1006	KTS3269045	25-4000	KTS5411114
65-2510	KTS3275906	50-5625	KTS5621475
65-2511	KTS3275909	30-4554	KTS6415020
65-2512	KTS3275912	53-1290	KTS4042001
65-2513	KTS3275919	53-1230	KTS4042201
65-2514	KTS3275925	53-1232	KTS4042202
37-4102	KTS3353217	53-1233	KTS4042203
37-4101	KTS3353318	53-1235	KTS4043201
65-1015	KTS4016101	53-1236	KTS4043202
65-1016	KTS4016102	53-1210	KTS4044101
65-1017	KTS4016103	53-1240	KTS4044201
65-1018	KTS4016104	53-1190	KTS4045602
61-2516	KTS4040517	53-1191	KTS4045603
50-6536	KTS4237613	53-1192	KTS4045604
49-6500	KTS4261599	37-4030	KTS0804119
49-6501	KTS4261600		
49-6502	KTS4261612		
49-6503	KTS4261634		
49-6504	KTS4261656		
49-6505	KTS4261678		
30-6557	KTS4294005		
30-6555	KTS4294205		
30-6898	KTS4420614		
30-6951	KTS4426108		
30-6995	KTS4615416		
30-7150	KTS4632319		
38-2520	KTS4656018		
30-4290	KTS5018319		
57-6546	KTS5026000		
25-4002	KTS5411314		
30-4088	KTS5452015		
50-8023	KTS5621020		
50-8024	KTS5621025		
30-4306	KTS6421123		
30-4300	KTS6422628		
79-5550	KTS6812100		
79-5555	KTS6821022		
79-5514	KTS7011002		
79-5515	KTS7011003		
79-5521	KTS7012002		
79-5522	KTS7012003		
25-5536	KTS7042408		
25-5540	KTS7043010		
25-5541	KTS7043012		
25-5542	KTS7043014		
25-5543	KTS7043016		
25-5544	KTS7043018		
25-5545	KTS7043020		
25-5549	KTS7043506		
25-5550	KTS7044708		
30-5580	KTS7050225		
30-5595	KTS7051424		
23-5574	KTS7059440		
23-5535	KTS7071508		
23-5536	KTS7071509		
23-5537	KTS7071511		
30-5655	KTS7210230		
30-5680	KTS7215544		
36-5080	KTS0883512		
34-4013	KTS0890518		
19-8090	KTS1207108		
19-8091	KTS1207112		
19-8092	KTS1207208		
30-8018	KTS1244721		
30-5801	KTS1259420		
19-5505	KTS1259520		
19-4133	KTS1271218		
30-2540	KTS2367911		
56-2507	KTS3354016		
37-4100	KTS3390214		
53-1500	KTS4051202		
53-1402	KTS4051203		
53-1405	KTS4051204		
53-1408	KTS4051205		
53-1524	KTS4051802		
53-1414	KTS4051803		
53-1417	KTS4051804		
53-1420	KTS4051805		
24-1010	KTS4406001		
30-7315	KTS5020518		



Jarit	K.T Surgico	Jarit	K.T Surgico	Jarit	K.T Surgico	Jarit	K.T Surgico	Jarit	K.T Surgico
099-262	KTS0828118	106-131	KTS1224114	121-165	KTS2422426	140-329	KTS1243113	205-161	KTS1867120
100-110	KTS0812014	106-140	KTS1226216	121-170	KTS2422213	150-160	KTS5452015	205-165	KTS1867025
100-130	KTS0811014	106-141	KTS1226316	121-172	KTS2422215	170-100	KTS0414008	205-166	KTS1867125
100-131	KTS0810214	106-150	KTS1224216	121-174	KTS2422218	170-102	KTS0413010	205-210	KTS1869031
100-138	KTS0810216	106-151	KTS1224316	121-185	KTS2421418	170-115	KTS0415222	205-226	KTS1801045
100-200	KTS0816014	106-170	KTS1231016	121-186	KTS2421420	170-120	KTS0415027	205-255	KTS1882118
100-202	KTS0816114	106-171	KTS1231116	121-187	KTS2421423	175-120	KTS0419025	205-293	KTS1876570
100-210	KTS0816017	106-180	KTS1231018	121-188	KTS2421426	175-121	KTS0419017	205-298	KTS1876545
100-212	KTS0816117	106-181	KTS1231118	121-189	KTS2421431	175-121	KTS0419020	205-310	KTS1876690
100-220	KTS0817017	106-190	KTS1231020	121-200	KTS2432418	180-100	KTS0228817	205-315	KTS1876418
100-222	KTS0817117	106-191	KTS1231120	121-205	KTS2432426	180-105	KTS0228820	205-316	KTS1891800
100-224	KTS0818017	106-196	KTS1231122	121-209	KTS2424526	180-200	KTS2010211	205-325	KTS1875400
100-225	KTS0818117	106-196	KTS1231124	121-210	KTS2424521	180-201	KTS2010213	220-100	KTS0891320
100-233	KTS0801918	106-200	KTS1232016	121-229	KTS2467618	180-202	KTS2010214	220-105	KTS0893723
100-250	KTS0828015	106-202	KTS1232116	121-230	KTS2469421	180-203	KTS2010216	220-180	KTS3036623
100-251	KTS0828115	106-210	KTS1232018	121-231	KTS2469521	180-225	KTS2010616	225-137	KTS3241113
100-260	KTS0828018	106-212	KTS1232118	121-260	KTS2480818	185-135	KTS3317300	225-139	KTS3240913
100-261	KTS0828318	106-220	KTS1232020	122-500	KTS1000212	185-135	KTS3317500	225-175	KTS3240217
100-272	KTS0828315	106-222	KTS1232120	122-502	KTS1000215	185-136	KTS3317301	225-210	KTS3244120
100-295	KTS0832115	106-230	KTS1226014	122-520	KTS1132218	185-181	KTS7412804	225-402	KTS3277171
100-297	KTS0832118	106-232	KTS1226114	122-530	KTS1022502	185-183	KTS7412806	225-405	KTS3282717
100-300	KTS0800315	110-102/		122-560	KTS1026116	185-192	KTS7434113	225-440	KTS3281106
100-512	KTS0890111	111-102/		122-600	KTS1026120	185-195	KTS7436013	225-441	KTS3281108
100-522	KTS0890114	109-102	KTS1411113	125-170	KTS0883112	185-200	KTS7433614	225-452	KTS3280118
100-526	KTS0890118	110-105	KTS1411713	125-180	KTS0880209	185-205	KTS7433114	225-453	KTS3280724
100-528	KTS0890120	110-110	KTS1410309	129-260	KTS1018812	185-230	KTS7436511	225-520	KTS3245001
100-539	KTS0895014	110-117	KTS1412511	130-105	KTS1010213	189-110	KTS1816215	230-098	KTS3210112
100-540	KTS0895018	110-118	KTS1412513	130-110	KTS1010214	189-112	KTS1816315	230-099	KTS3210114
101-130	KTS0811414	110-119	KTS1412110	130-115	KTS1010216	190-110	KTS1816017	230-124	KTS9911050
101-215	KTS0816414	110-120	KTS1412113	130-116	KTS1010415	190-112	KTS1816117	230-131	KTS3215116
101-217	KTS0816514	110-121	KTS1416514	130-125	KTS1010225	190-118	KTS4237613	230-135	KTS3215418
101-218	KTS0817415	110-150	KTS1270416	130-160	KTS1012014	190-150	KTS1822001	230-136	KTS3215518
101-220	KTS0816417	110-157	KTS1271018	130-161	KTS1012214	200-105	KTS1813121	230-150	KTS3214220
101-220	KTS0817417	110-159	KTS1271020	130-165	KTS1012016	200-130	KTS1813519	230-196	KTS3221315
101-222	KTS0816517	110-161	KTS1271218	130-166	KTS1012215	200-132	KTS1834822	230-205	KTS3220118
101-222	KTS0817517	110-165	KTS0610301	130-175	KTS1012025	200-134	KTS1830419	230-207	KTS3220518
101-225	KTS0816423	110-167	KTS0610400	130-230	KTS1017012	200-135	KTS1823421	230-211	KTS3221718
101-226	KTS0816523	110-169	KTS0610700	130-231	KTS1017612	200-150	KTS1814000	230-230	KTS3223022
101-250	KTS0828615	110-170	KTS0610302	130-232	KTS1014212	200-151	KTS1814001	230-235	KTS3223122
101-251	KTS0828715	110-171	KTS0610303	130-234	KTS1018012	200-152	KTS1814002	230-240	KTS3217523
101-258	KTS0828518	110-175	KTS0610300	130-235	KTS1018612	200-175	KTS1819140	230-245	KTS3225902
101-261	KTS0828718	110-182	KTS0610116	130-236	KTS1000112	200-176	KTS1819162	230-266	KTS3224523
101-267	KTS0828523	111-100	KTS1411109	130-240	KTS1018212	200-178	KTS1819164	230-267	KTS3224723
101-272	KTS0828530	116-102	KTS1617018	130-245	KTS1014412	200-191	KTS1852075	230-305	KTS3228618
101-276	KTS0828723	116-105	KTS1617225	130-250	KTS1015212	200-202	KTS1818325	230-315	KTS3228415
101-300	KTS0834411	116-107	KTS1617025	130-255	KTS1015412	200-206	KTS1818750	230-320	KTS3229027
101-301	KTS0834511	116-112	KTS1617118	130-260	KTS1019912	200-217	KTS1816825	230-325	KTS3230327
101-311	KTS0837511	116-115	KTS1617325	130-265	KTS1018915	200-218	KTS1816817	230-330	KTS3227414
101-350	KTS0883512	116-117	KTS1617125	130-286	KTS1027215	200-219	KTS1816820	230-332	KTS3227419
102-100	KTS0716014	120-095	KTS2434413	130-289	KTS1024618	200-221	KTS1817025	240-099	KTS3246996
102-110	KTS0716017	120-112	KTS2415014	130-290	KTS1024818	200-222	KTS1817032	240-100	KTS3246997
102-112	KTS0716117	120-118	KTS2410013	130-293	KTS1045118	200-223	KTS1817040	240-101	KTS3246998
102-120	KTS0728115	120-120	KTS2413013	130-294	KTS1045218	200-225	KTS1817050	240-102	KTS3247000
102-122	KTS0728118	120-130	KTS2416115	130-390	KTS1025615	200-270	KTS1828002	240-103	KTS3247001
102-125	KTS0728123	120-130	KTS2416015	130-392	KTS1025620	200-271	KTS1828003	240-104	KTS3247002
102-131	KTS0734111	120-131	KTS2416018	130-395	KTS1025625	200-272	KTS1828004	240-105	KTS3247003
103-250	KTS0828415	120-135	KTS2418016	130-397	KTS1045623	200-273	KTS1828006	240-106	KTS3247004
105-090	KTS1222210	120-140	KTS2418018	130-425	KTS1022017	200-280	KTS1828102	240-107	KTS3247005
105-091	KTS1222310	120-145	KTS2418020	130-430	KTS1023017	200-281	KTS1828103	240-108	KTS3247006
105-095	KTS1222212	120-147	KTS2418030	130-500	KTS1080209	200-282	KTS1828104	240-120	KTS3261297
105-096	KTS1222312	120-150	KTS2442014	130-520	KTS1081210	200-283	KTS1828106	240-121	KTS3261298
105-110	KTS1223012	120-155	KTS2442017	130-520	KTS1081610	200-315	KTS1836440	240-122	KTS3261200
105-112	KTS1223112	120-210	KTS2424121	130-525	KTS1081310	201-156	KTS1815231	240-123	KTS3261201
105-114	KTS1221010	121-100	KTS2411412	130-525	KTS1081710	201-160	KTS1815251	240-124	KTS3261202
105-115	KTS1221110	121-105	KTS2415412	135-095	KTS6410214	201-165	KTS1815851	240-125	KTS3261203
105-142	KTS1228014	121-118	KTS2410413	135-100/		201-167	KTS1815876	240-126	KTS3261204
105-143	KTS1228016	121-120	KTS2413413	136-100	KTS6411015	202-155	KTS1815001	240-127	KTS3261205
105-144	KTS1228114	121-121	KTS2414413	135-105/		202-156	KTS1815002	240-128	KTS3261206
105-145	KTS1228116	121-125	KTS2419818	136-105	KTS6412015	202-157	KTS1815003	240-140	KTS3261397
105-160	KTS1227014	121-129	KTS2416430	135-107	KTS6412416	202-160	KTS1815004	240-141	KTS3261398
105-161	KTS1227114	121-130/		135-110/		202-166	KTS1815626	240-142	KTS3261300
105-200	KTS1232216	123-130	KTS2416415	136-110	KTS6412019	205-100	KTS1862107	240-143	KTS3261301
105-223	KTS1232024	121-133	KTS2416423	135-125	KTS6415020	205-110	KTS1863410	240-144	KTS3261302
105-225	KTS1232026	121-135	KTS2418416	135-127	KTS6414815	205-111	KTS1863510	240-145	KTS3261303
105-242	KTS1234116	121-140	KTS2418418	135-130	KTS6413015	205-120	KTS1880218	240-146	KTS3261304
105-245	KTS1234120	121-145	KTS2418420	135-181	KTS6416016	205-123	KTS1880209	240-147	KTS3261305
105-255	KTS1236015	121-146	KTS2418424	135-185	KTS6416018	205-124	KTS1880211	240-148	KTS3261306
105-320	KTS1235815	121-147	KTS2418430	135-190	KTS6416020	205-125	KTS1880214	240-197	KTS3283030
106-100/		121-150	KTS2441414	135-195	KTS6416024	205-140	KTS1865010	240-200	KTS3249213
105-100	KTS1222012	121-155	KTS2441417	135-250	KTS1032023	205-141	KTS1865110	240-203	KTS3249216
106-102/		121-160	KTS2422413	135-261	KTS1231126	205-150	KTS1867013		
105-102	KTS1222112	121-161	KTS2422415	135-355	KTS6421123	205-151	KTS1867113		
106-120	KTS1225014	121-162	KTS2422418	140-224	KTS1247116	205-155	KTS1867016		
106-121	KTS1225114	121-163	KTS2422420	140-225	KTS1247614	205-156	KTS1867116		
106-130	KTS1224014	121-164	KTS2422423	140-231	KTS1245918	205-160	KTS1867020		



Jarit	K.T Surgico	Jarit	K.T Surgico	Jarit	K.T Surgico	Jarit	K.T Surgico
240-205	KTS3249220	285-455	KTS4029511	410-150	KTS0836210	505-160	KTS5411114
250-250	KTS3265706	285-456	KTS4029512	410-151	KTS0836310	505-200	KTS5411111
250-251	KTS3265712	285-457	KTS4029514	410-152	KTS0806109	505-205	KTS5411314
250-252	KTS3265719	285-480	KTS0440310	410-153	KTS0910211	505-210	KTS5411615
250-253	KTS3265725	285-482	KTS0440315	410-154	KTS0234015	505-225	KTS5413115
250-254	KTS3265732	285-484	KTS0440320	410-166	KTS4261599	360-140	KTS6812100
250-255	KTS3265738	290-100	KTS1862410	410-170	KTS4261600	360-141	KTS6821022
250-260	KTS3265806	290-101	KTS1862910	410-220	KTS4261612	360-143	KTS2011815
250-261	KTS3265812	290-122	KTS1867820	410-335	KTS4261634	360-150	KTS6620108
250-262	KTS3265819	290-125	KTS1877022	415-145	KTS4261656	360-260	KTS6620110
250-263	KTS3265825	290-152	KTS4040517	420-135	KTS4261678	360-280	KTS6620112
250-264	KTS3265832	300-105	KTS5621020	420-145	KTS4420614	360-281	KTS6620114
250-265	KTS3265838	300-115	KTS5621025	420-150	KTS4615416	360-282	KTS6620116
250-398	KTS3277214	300-200	KTS5621475	420-170	KTS4426108	360-283	KTS6620118
250-405	KTS3269045	305-101	KTS0828123	420-171	KTS4426114	360-284	KTS6620120
250-406	KTS3269090	305-120	KTS0865022	420-172	KTS4406001	360-285	KTS6620122
250-410	KTS3277201	305-140	KTS0816023	440-304	KTS4406002	380-100	KTS6620124
250-420	KTS3275903	305-145	KTS0816028	450-106	KTS4409650	380-105	KTS6620126
250-421	KTS3275906	305-146	KTS0816128	450-110	KTS4409710	380-150	KTS6620128
250-422	KTS3275909	305-200	KTS1235820	450-140	KTS4409601	380-155	KTS6620130
250-423	KTS3275912	305-230	KTS6412024	450-141	KTS4499801	380-165	KTS6620132
250-424	KTS3275915	305-295	KTS1243520	450-150	KTS4499802	380-166	KTS6620134
250-425	KTS3275919	305-296	KTS1243524	450-151	KTS4431040	380-360	KTS6620136
250-426	KTS3275925	305-297	KTS1243528	450-200	KTS4611201	380-362	KTS6621314
250-480	KTS3276118	305-310	KTS1246128	450-201	KTS4611202	380-366	KTS7011002
250-481	KTS3276120	305-355	KTS1244721	450-210	KTS4612035	385-100	KTS7011003
250-482	KTS3276123	305-365	KTS1245124	450-310	KTS4612050	385-104	KTS7012002
250-483	KTS3276125	305-370	KTS1247922	450-450	KTS4612075	385-200	KTS7012003
260-200	KTS3297530	310-204	KTS5632515	450-501	KTS4612090	400-100	KTS7043506
260-201	KTS3297327	310-310	KTS0848025	450-510	KTS4615418	400-101	KTS7043010
260-205	KTS9900006	310-315	KTS0848040	450-524	KTS4615720	400-130	KTS7043012
270-185	KTS2367911	310-320	KTS0848060	450-529	KTS0840918	400-130	KTS7043014
270-195	KTS1299215	310-350	KTS0803025	460-100	KTS4637415	400-132	KTS7043016
275-155	KTS3354016	310-351	KTS0803045	460-101	KTS4639004	400-134	KTS7043018
275-510	KTS3354517	310-352	KTS0803060	460-102	KTS4642322	400-150	KTS7043020
275-515	KTS3354522	310-353	KTS0803090	460-103	KTS4630418	400-152	KTS7044708
275-525	KTS3357218	310-354	KTS0803012	460-104	KTS4642516	400-180	KTS7058132
275-550	KTS3304056	310-369	KTS0800217	460-105	KTS4632519	400-255	KTS7071508
275-555	KTS3306019	310-370	KTS0800117	460-106	KTS4632319	400-261	KTS7071509
275-560	KTS3353217	320-101	KTS1131020	460-107	KTS4657123	400-277	KTS7071511
275-565	KTS3356614	320-102	KTS1131024	460-400	KTS4625121	400-279	KTS7059440
275-565	KTS3390214	320-103	KTS1131030	460-401	KTS4625800	400-281	KTS7059640
275-589	KTS3353318	320-109	KTS1130830	460-405	KTS1821301	405-150	KTS7061704
277-141	KTS2429818	320-110	KTS1130816	460-420	KTS1821302	405-153	KTS7063101
280-041	KTS4051201	320-111	KTS1130820	460-430	KTS1821303	405-180	KTS7063103
280-060	KTS4051202	320-112	KTS1130824	475-100	KTS1821601	405-200	KTS1623126
280-061	KTS4051203	320-113	KTS1130616	475-130	KTS1821602	405-225	KTS7050225
280-062	KTS4051204	320-114	KTS1130620	475-337	KTS1821605	410-135	KTS7050425
280-063	KTS4051205	320-115	KTS1130624	475-550	KTS1821607	410-136	KTS1414016
280-386	KTS4044101	320-152	KTS1209307	475-551	KTS1821610	410-137	KTS7050927
280-400	KTS4042001	320-170	KTS1207108	475-552	KTS4656218	505-240	KTS7051424
280-405	KTS4042201	320-171	KTS1207208	475-553	KTS4656018	505-241	KTS7051524
280-406	KTS4042202	320-172	KTS1207112	475-554	KTS1816414	505-250	KTS7053221
280-407	KTS4042203	320-180	KTS1310212	475-555	KTS3274624	505-315	KTS1260021
280-415	KTS4043201	320-185	KTS1311112	475-556	KTS0809711	505-319	KTS1259420
280-416	KTS4043202	320-390	KTS1325325	475-557	KTS4650618	505-320	KTS1259520
280-425	KTS4044201	320-391	KTS1325324	475-558	KTS4651616	510-120	KTS7053926
280-430	KTS4045602	320-392	KTS1325320	475-559	KTS4652021	510-155	KTS0855010
280-431	KTS4045603	320-412	KTS1342114	475-560	KTS4652209	510-220	KTS7210230
280-432	KTS4045604	320-550	KTS1334318	475-561	KTS4652207	510-300	KTS7215544
281-041	KTS4051801	350-094	KTS4210103	475-562	KTS4652205	600-455	KTS1817612
281-060	KTS4051802	350-360	KTS1822816	475-563	KTS4631301		
281-061	KTS4051803	350-370	KTS1821201	475-564	KTS4826514		
281-062	KTS4051804	350-371	KTS1821202	475-635	KTS4831010		
281-063	KTS4051805	355-120	KTS4294005	500-105	KTS4832000		
281-415	KTS4043401	355-123	KTS4294008	500-110	KTS4832010		
285-106	KTS1241118	355-143	KTS4294607	500-120	KTS4811501		
285-120	KTS1030218	355-150	KTS4315005	500-125	KTS4811502		
285-125	KTS1031018	355-404	KTS2482514	500-230	KTS5020219		
285-200	KTS1821413	355-405	KTS2482214	500-232	KTS5020319		
285-365	KTS4016101	355-412	KTS2467414	500-233	KTS5020518		
285-366	KTS4016102	355-413	KTS2467614	500-234	KTS5018319		
285-367	KTS4016103	355-415	KTS2466314	500-235	KTS0842717		
285-368	KTS4016104	355-416	KTS2466114	500-236	KTS5025228		
285-380	KTS4641518	360-100	KTS0834011	500-237	KTS5026000		
285-450	KTS4029506	360-101	KTS0834111	500-240	KTS5025903		
285-451	KTS4029507	360-104	KTS0834009	500-252	KTS5025955		
285-452	KTS4029508	360-105	KTS0834109	500-302	KTS5227212		
285-453	KTS4029509	360-110	KTS0834010	500-303	KTS5227214		
285-454	KTS4029510	360-111	KTS0834110	500-305	KTS5227216		
		360-117	KTS0838410	500-335	KTS5227218		
		360-118	KTS0838510	500-336	KTS5227220		
		360-131	KTS0837111	505-140	KTS5227222		
		360-138	KTS0836211	505-145	KTS5227224		
		360-139	KTS0836311	505-146	KTS5227226		



V Mueller	K.T Surgico	V Mueller	K.T Surgico	V Mueller	K.T Surgico	V Mueller	K.T Surgico	V Mueller	K.T Surgico
AU11355-02	KTS3210114	CH5902	KTS1130820	M01601	KTS0828518	OP5668	KTS0910211	OS6085	KTS3036623
AU11370	KTS3215418	CH5902-1	KTS1130420	M01601-001	KTS0828718	OP5680	KTS0836210	OS910-032	KTS3244120
AU11385	KTS3215518	CH5904	KTS1130824	M01650	KTS0842717	OP5690	KTS0836310	OS911-003	KTS3241113
AU12220	KTS4431040	CH5904-1	KTS1130424	M01700	KTS5018319	OP5692	KTS0806109	OS911-008	KTS3240913
AU14122	KTS4499802	CH5906	KTS1130830	M01760	KTS1617018	OP5693	KTS0836311	OS928	KTS3282717
AU15030	KTS0440315	CH6076	KTS1207112	M01780	KTS1617118	OP5694	KTS0836211	OS930	KTS1816414
AU15032	KTS0440310	CH6085	KTS1209307	M01970	KTS5025228	OP5710	KTS0837111	OS993-001	KTS3261297
AU15034	KTS0440320	CH6095	KTS1310212	M01995	KTS5026000	OP5740	KTS0837511	OS993-002	KTS3261298
AU5000	KTS4406007	CH6097	KTS1311112	NL1003	KTS1823421	OP5780	KTS0880209	OS993-003	KTS3261200
AU5001	KTS4406001	CH6950-003	KTS1342114	NL1090	KTS4016101	OP7020-012	KTS4261612	OS993-004	KTS3261201
AU5002	KTS4406002	CH7296	KTS1325326	NL1091	KTS4016102	OP7020-034	KTS4261634	OS993-005	KTS3261202
AU6120	KTS4420614	CH7298	KTS1325327	NL1092	KTS4016103	OP7020-056	KTS4261656	OS993-006	KTS3261203
AU6410	KTS1862910	CH7300	KTS1325320	NL1093	KTS4016104	OP7020-078	KTS4261678	OS993-007	KTS3261204
AU6600	KTS3247497	CH7302	KTS1325324	NL1400	KTS1014412	OP7020-300	KTS4261600	OS993-008	KTS3261205
AU6601	KTS3247498	CH7304	KTS1325325	NL1401	KTS1017612	OP7020-700	KTS4261599	OS993-009	KTS3261206
AU6602	KTS3247500	CH7950	KTS1334318	NL1410	KTS1017012	OP7380	KTS2467714	RH100	KTS4612090
AU6603	KTS3247501	CH8591-001	KTS2469421	NL1411	KTS1018612	OP7399	KTS2467614	RH100-1	KTS4612075
AU6604	KTS3247502	CH8592	KTS2429818	NL1415	KTS1000112	OP7412-002	KTS2482214	RH1010	KTS4642516
AU6605	KTS3247503	DE10718	KTS0228817	NL1440	KTS1031018	OP7413	KTS2482514	RH102	KTS4612050
AU6606	KTS3247504	DE10720	KTS0228820	NL1450	KTS1030218	OP7418-001	KTS2466314	RH1045	KTS4630418
AU6607	KTS3247505	GL1000	KTS7051424	NL1460	KTS1022502	OP7418-002	KTS2466114	RH11	KTS4611202
AU6720	KTS9911050	GL1020	KTS7051524	NL1470	KTS1023017	OP7511	KTS2414413	RH1104	KTS1821301
AU6780	KTS3221718	GL11	KTS7011002	NL1480	KTS1022017	OP7512	KTS2411412	RH1105	KTS1821302
AU6784-002	KTS3221715	GL12	KTS7011003	NL1900	KTS4029507	OP7513	KTS2413413	RH1106	KTS1821303
AU6791	KTS3220118	GL1320	KTS7058132	NL1901	KTS4029508	OP7550	KTS1221010	RH1115	KTS1821601
BE390	KTS5411615	GL1451	KTS7044708	NL1902	KTS4029510	OP7551	KTS1221110	RH1120	KTS1821602
BE395	KTS5413115	GL1460	KTS7043506	NL1903	KTS4029512	OS1060-006	KTS3265706	RH1125	KTS1821605
BE460	KTS5411114	GL1461-001	KTS7043010	NL2440	KTS1821413	OS1060-013	KTS3265712	RH1130	KTS1821607
BE465	KTS5411314	GL1461-002	KTS7043012	NL3530	KTS1241118	OS1060-019	KTS3265719	RH1135	KTS1821610
BE470	KTS0415222	GL1461-003	KTS7043014	NL4271-072	KTS4051202	OS1060-025	KTS3265725	RH12	KTS4611201
CH1000	KTS1819140	GL1461-004	KTS7043016	NL4P71-073	KTS4051203	OS1060-032	KTS3265732	RH1540	KTS3274624
CH1005	KTS1819162	GL1461-005	KTS7043018	NL4271-074	KTS4051204	OS1060-038	KTS3265738	RH1545	KTS4625800
CH1010	KTS1819164	GL1461-006	KTS7043020	NL4271-075	KTS4051205	OS1061-006	KTS3265806	RH1634	KTS0809711
CH1520	KTS1045623	GL1600-002	KTS7071508	NL4271-082	KTS4051802	OS1061-013	KTS3265812	RH1640	KTS0828015
CH1560	KTS6412024	GL1600-003	KTS7071509	NL4271-083	KTS4051803	OS1061-019	KTS3265819	RH1650	KTS0828115
CH1570	KTS6416024	GL1600-004	KTS7071511	NL4271-084	KTS4051804	OS1061-025	KTS3265825	RH1650-S	KTS0728115
CH1590	KTS1235820	GL1775	KTS7059440	NL4271-085	KTS4051805	OS1061-032	KTS3265832	RH1651	KTS0828515
CH1653	KTS0828415	GL1925	KTS7061704	NL500	KTS3214220	OS1061-038	KTS3265838	RH1651-001	KTS0828715
CH1682	KTS1245124	GL1942	KTS7063101	NL5131	KTS1852075	OS1079-001	KTS3275903	RH1653-001	KTS0828615
CH1687	KTS1244721	GL1943	KTS7063103	NL5382	KTS1882118	OS1079-002	KTS3275906	RH1770	KTS4657123
CH1721	KTS1243520	GL31	KTS7012002	NL5383	KTS1882018	OS1079-003	KTS3275909	RH1790	KTS4625121
CH1722	KTS1243524	GL3101	KTS1260021	NL600	KTS3217523	OS1079-004	KTS3275912	RH2010-003	KTS4652021
CH1723	KTS1243528	GL3120	KTS1259520	NL6090	KTS4042001	OS1079-005	KTS3275915	RH2050	KTS4650618
CH1724-001	KTS1246128	GL3121	KTS1259420	NL610	KTS3223022	OS1079-006	KTS3275919	RH2081	KTS4651616
CH1725	KTS1247922	GL32	KTS7012003	NL6146	KTS4042201	OS1079-007	KTS3275925	RH2410	KTS4632519
CH2006	KTS0828523	GL3210	KTS1231122	NL6147	KTS4042202	OS1110	KTS3225902	RH2430	KTS4632319
CH2006-001	KTS0828723	GL3211	KTS1231126	NL6148	KTS4042203	OS1116	KTS3210112	RH2560	KTS2410013
CH2032	KTS0828123	GL3550	KTS0855010	NL6150	KTS4043201	OS1316	KTS4040517	RH2561	KTS2410413
CH2032-S	KTS0728123	GL500	KTS1880218	NL6151	KTS4043202	OS1512	KTS3268625	RH2636	KTS4656018
CH2070	KTS0865022	GL501	KTS1880214	NL6160	KTS4044101	OS1517	KTS3274101	RH2637	KTS4656218
CH2080	KTS0816028	GL5260	KTS7053926	NL6165	KTS4044201	OS3030	KTS3304056	RH3915	KTS4631301
CH2085	KTS0816023	GL5301-4	KTS7210230	NL630	KTS3224523	OS3040	KTS3357218	RH4700-002	KTS5227210
CH2090	KTS0816128	GL5500-4	KTS7215544	NL631	KTS3224723	OS3046	KTS3354215	RH4700-003	KTS5227212
CH2400	KTS2416018	GL600	KTS1623126	NL6760-001	KTS3261397	OS3052-001	KTS3354517	RH4700-004	KTS5227214
CH2411	KTS2432426	GL650	KTS1617025	NL6760-002	KTS3261398	OS3052-002	KTS3354522	RH4700-005	KTS5227216
CH2416	KTS2432418	GL651	KTS1617225	NL6760-003	KTS3261300	OS3062	KTS3354016	RH4700-006	KTS5227218
CH2419	KTS2416615	GL660	KTS1617125	NL6760-004	KTS3261301	OS3065	KTS3200218	RH4700-007	KTS5227222
CH2425	KTS2416423	GL661	KTS1617325	NL6760-005	KTS3261302	053070	KTS3306019	RH4700-008	KTS5227224
CH2428	KTS2416430	GL850	KTS7050425	NL6760-006	KTS3261303	053071	KTS3390214	RH4700-009	KTS5227226
CH2433	KTS2426427	GL91	KTS9930141	NL6760-007	KTS3261304	OS3077	KTS3353217	RH600	KTS0840918
CH2435	KTS2418430	GL950	KTS7053221	NL6760-008	KTS3261305	OS3090-002	KTS3353318	RH750	KTS4641518
CH2452	KTS2467414	GU4010-008	KTS6620108	NL6760-009	KTS3261306	OS4080-003	KTS3283030	RH760	KTS4639004
CH2456	KTS2467618	GU4010-010	KTS6620110	NL6960	KTS3240217	OS4130-002	KTS3277171	RH980	KTS4642322
CH2481	KTS2421418	GU4010-012	KTS6620112	NL960	KTS1877022	OS4150-001	KTS3276118	SU10500	KTS1247116
CH2482	KTS2421423	GU4010-014	KTS6620114	OP2321	KTS1822816	OS4150-002	KTS3276120	SU10530	KTS1243113
CH2483	KTS2421426	GU4010-016	KTS6620116	OP26-001	KTS4210103	OS4150-003	KTS3276123	SU10560	KTS1245918
CH2484	KTS2421431	GU4010-018	KTS6620118	OP2720	KTS0234015	OS4150-004	KTS3276125	SU10810450	KTS2010211
CH2506	KTS2422213	GU4010-020	KTS6620120	OP3302	KTS4315005	OS4160-002	KTS3281106	SU10810500	KTS2010213
CH2507	KTS2422215	GU4010-022	KTS6620122	OP3303	KTS4315003	OS4160-003	KTS3281108	SU10810550	KTS2010214
CH2508	KTS2422218	GU4010-024	KTS6620124	OP3410	KTS4294208	OS4160-004	KTS3280118	SU10810600	KTS2010216
CH470	KTS3228618	GU4010-026	KTS6620126	OP3412	KTS4294205	OS4160-011	KTS3280724	SU110	KTS6812100
CH5660	KTS0848060	GU4010-028	KTS6620128	OP3413	KTS4294008	OS4450	KTS3297327	SU13000	KTS0413010
CH5662	KTS0848040	GU4010-030	KTS6620130	OP3505	KTS4294407	OS4455	KTS3297530	SU13002	KTS0414008
CH5665	KTS0848025	GU4010-032	KTS6620132	OP3515	KTS4294607	OS4464-002	KTS9900006	SU1403-001	KTS0610301
CH5675	KTS0832115	GU4010-034	KTS6620134	OP5000	KTS0834009	OS4610	KTS3228415	SU1403-002	KTS0610302
CH5676	KTS0832118	GU4010-036	KTS6620136	OP5005	KTS0834109	OS4640	KTS3227414	SU1403-003	KTS0610303
CH5720	KTS0803025	GU4230	KTS6621314	OP5525	KTS0834511	OS4700	KTS3229027	SU1404-001	KTS0610400
CH5725	KTS0803045	GU8772	KTS1826513	OP5526	KTS0834411	OS4720	KTS3230327	SU1407	KTS0610700
CH5730	KTS0803060	GU8773	KTS1826515	OP5530	KTS0834010	OS4730	KTS3230215	SU15085	KTS7436511
CH5735	KTS0803090	M0110	KTS4826514	OP5530-S	KTS0734111	OS5003	KTS3245001	SU15100	KTS7436013
CH5893	KTS1130616	M0140	KTS4832000	OP5531	KTS0837602	OS5510	KTS1299215		
CH5895	KTS1130620	M0141	KTS4832010	OP5540	KTS0834110	OS5521	KTS2367911		
CH5897	KTS1130624	M01410	KTS5020518	OP5541	KTS0837702	OS5531	KTS1865110		
CH5900	KTS1130816	M01420	KTS5020219	OP5542	KTS0837701	OS5532	KTS1865010		
CH5900-1	KTS1130416	M01600-S	KTS0728118	OP5666	KTS0910411	OS6075	KTS0893723		



V Mueller	K.T Surgico	V Mueller	K.T Surgico	V Mueller	K.T Surgico
SU15120	KTS7433614	SU2542	KTS1045118	SU3554-01	KTS1828106
SU15130	KTS7434113	SU2601	KTS1080209	SU3560-01	KTS1828002
SU15330004	KTS7412804	SU2610	KTS1081610	SU3561-01	KTS1828003
SU15330006	KTS7412806	SU2611	KTS1081710	SU3562-01	KTS1828004
SU16005	KTS2416415	SU2630	KTS1082815	SU3564-01	KTS1828006
SU16020	KTS2418016	SU2698	KTS1222212	SU3580	KTS1836440
SU16025	KTS2418018	SU2699	KTS1222312	SU3660	KTS1813121
SU16030	KTS2418020	SU2700	KTS1222012	SU3675	KTS1813519
SU16060	KTS2418416	SU2702	KTS1222112	SU3760	KTS1834822
SU16061	KTS2418418	SU2705	KTS1223012	SU3770	KTS1830419
SU16062	KTS2418420	SU2707	KTS1223112	SU3785	KTS1816017
SU16063	KTS2418424	SU2710	KTS1227014	SU3786	KTS1816117
SU16080	KTS2415014	SU2712	KTS1227114	SU3960	KTS5452015
SU16081	KTS2415412	SU2720	KTS1224014	SU4054	KTS6411015
SU16091	KTS2419818	SU2722	KTS1224114	SU4055	KTS6412015
SU16120	KTS2424121	SU2723	KTS1228014	SU4056	KTS6412019
SU16122	KTS2424521	SU2724	KTS1228114	SU4060	KTS6414815
SU16124	KTS2424526	SU2725	KTS1228016	SU4070	KTS6415020
SU16145	KTS2422415	SU2726	KTS1228116	SU5000	KTS6416016
SU16146	KTS2422418	SU2730	KTS1225014	SU5001	KTS6416020
SU16147	KTS2422420	SU2732	KTS1226216	SU5002	KTS6416018
SU16148	KTS2422426	SU2735	KTS1225114	SU5075	KTS1032023
SU16150	KTS2426027	SU2750	KTS1231016	SU6082	KTS6421123
SU16252	KTS2441414	SU2752	KTS1231018	SU740	KTS2011815
SU16252-001	KTS2441417	SU2754	KTS1231020	SU960	KTS6413015
SU16253	KTS2441419	SU2758	KTS1231030	SU970	KTS1236015
SU1702	KTS0811014	SU2760	KTS1231116	SU980	KTS1235815
SU1712	KTS0812014	SU2762	KTS1231118		
SU1722	KTS0810214	SU2764	KTS1231120		
SU1723	KTS0810215	SU2768	KTS1231130		
SU1755	KTS0811414	SU2775	KTS1234116		
SU1791	KTS0859411	SU2777	KTS1234120		
SU1800	KTS0816014	SU2800	KTS1232016		
SU1800-S	KTS0716014	SU2802	KTS1232018		
SU1801	KTS0816017	SU2804	KTS1232020		
SU1801-S	KTS0716017	SU2806	KTS1232026		
SU1803	KTS0816414	SU2810	KTS1232116		
SU1803-002	KTS0817415	SU2812	KTS1232118		
SU1804	KTS0816417	SU2814	KTS1232120		
SU1804-002	KTS0817417	SU2854	KTS1271018		
SU1810	KTS0816114	SU2855	KTS1271218		
SU1811	KTS0816117	SU2856	KTS1271020		
SU1811-S	KTS0716117	SU2900	KTS1411109		
SU1813	KTS0816514	SU2905	KTS1411113		
SU1814	KTS0816517	SU2910	KTS1411713		
SU1814-002	KTS0817517	SU2920	KTS1410309		
SU1821	KTS0817017	SU2935	KTS1412511		
SU1826	KTS0817117	SU2936	KTS1412513		
SU1860	KTS0818017	SU2937	KTS1412110		
SU1862	KTS0818117	SU2938	KTS1412113		
SU1945	KTS0828318	SU3000	KTS1876418		
SU1952	KTS0801918	SU3008	KTS1876690		
SU1980	KTS0883112	SU3055	KTS 1875400		
SU1988	KTS0883512	SU3060	KTS1891800		
SU2003	KTS0890111	SU3110	KTS1867013		
SU2005	KTS0890114	SU3110-001	KTS1867113		
SU2005-001	KTS0890518	SU3112	KTS1867016		
SU2007	KTS0890118	SU3112-001	KTS1867116		
SU2009	KTS0890120	SU3114	KTS1867020		
SU2012-001	KTS0895014	SU3114-001	KTS1867120		
SU2014-001	KTS0895018	SU3116	KTS1867025		
SU2015	KTS0891322	SU3116-001	KTS1867125		
SU2086	KTS0816423	SU3145	KTS1863410		
SU2096	KTS0816523	SU3146	KTS1863510		
SU2301	KTS1010213	SU3147	KTS1862107		
SU2302	KTS1010214	SU3190	KTS1869031		
SU2303	KTS1010216	SU3303	KTS1818325		
SU2307	KTS1010225	SU3305	KTS1818750		
SU2332	KTS1012014	SU3340	KTS1817025		
SU2333	KTS1012016	SU3342	KTS1817032		
SU2337	KTS1012025	SU3344	KTS1817040		
SU2390	KTS1015212	SU3346	KTS1817050		
SU2391	KTS1015412	SU3460	KTS1815851		
SU2450	KTS1025615	SU3462	KTS1815626		
SU2452	KTS1025620	SU3464	KTS1815876		
SU2454	KTS1025625	SU3470	KTS1815001		
SU2500	KTS1018915	SU3472	KTS1815002		
SU2504	KTS1018812	SU3474	KTS1815003		
SU2510	KTS1024818	SU3478	KTS1815004		
SU2511	KTS1024618	SU3490-001	KTS1814001		
		SU3490-002	KTS1814002		
		SU3491	KTS1814000		
		SU3550-0	KTS1828102		
		SU3551-01	KTS1828103		
		SU3552-01	KTS1828104		





K.T. SURGICO

K.T.Surgico Instruments Care & Handling Specifications

- Brand new instruments should be cleaned prior to first sterilization.
- Contaminated instruments should be processed as soon as possible.
- Stubborn protein particles can be removed with a scrub brush. Do not use steel wool, abrasives or an acid rinse.
- All joints on instruments prior to preparation.
- Lubrication is vital to a long instrument life. Avoid silicon lubricants because they tend to build up and mix with debris to clog moving parts, which becomes almost impossible to remove. Debris buildup can have a "rust-like" appearance.
- Rinse your cleaned instruments in demineralized water. Be sure to remove all residual-cleaning compounds before sterilization, as they can cause stains.
- Dry the instruments thoroughly after rinsing.

K.T.SURGICO WARRANTY PROGRAM

Corporate Mission

As a worldwide exporter / manufacturers, providing the healthcare professional with the highest level of quality instruments, products and services through our business partners.


We are dedicated to always surpass the expectations of our customer.

Uncompromised quality, each time, every time, guaranteed.

Guarantee

All K. T. Surgico Instruments are guaranteed to be free from defects in workmanship and material. Any K. T. Surgico Instruments which proves defective in workmanship or material, will either be repaired or replaced without charge. Our liability under this guarantee shall be limited to repair or replacement of defective merchandise.

K. T. SURGICO
Sialkot Pakistan.



Certified Under:
cGMP - QSR BY FDA
(U.S.A)
ISO - 9001 : 2000
CE Mark

K.T. SURGICO

Manufacturers & Exporters of Quality Dental Instruments



P.O.Box: 1830 Rangpura
K.T. Road,
Sialkot - Pakistan

Tel : 0092-52-4598939 - 4583501
Fax : 0092-52-4598940 - 4583502
E-mail : kts@skt.comsats.net.pk
info@ktsurgico.com

WWW.KTSURGICO.COM